



MOLECULAR DIAGNOSIS RAW MATERIALS AND SOLUTION

BIORI PRODUCT BROCHURE **2023**

Focus on raw material development and
provide high-quality services

Zhuhai Biori Biotechnology Co.,Ltd

ABOUT US

COMPANY PROFILE


Founded in 2012, Zhuhai Biori Biotechnology Co.,Ltd has been focusing on the core raw materials of nucleic acid diagnostic reagents, and was successively awarded as “The National High-tech Enterprise”, “The Specialized, Refined, Differentiated, and Innovative Little Giant Enterprise in Guangdong Province” and “The Unicorn Seed Enterprise in Zhuhai”. Through the rapid development in recent years, Biori has successfully developed full kinds of raw materials(key enzyme and premix) for nucleic acid diagnosis, which have been widely used by many companies in their development and production of nucleic acid diagnostic reagents. At the same time, Biori is increasing its investment in the development and research of mRNA vaccine raw materials and a variety of constant temperature amplification raw materials, and related products have been launched one after another.

Focus on raw material development and provide professional service. Biori is entering a new era of industry development at a high speed, and we will firmly pore the upstream core technology to facilitate the development of health industry with high-quality products and services!

QUALIFICATION HONOR



National High-Tech Enterprise



Small and Medium-sized the Specialized, Refined, Differentiated, and Innovative Enterprise in zhuhai xiangzhou district--Guangdong Province



Zhuhai "Unicorn" seed Enterprise



ISO9001: Certification of 2015 system



Zhuhai advanced group in fighting against COVID-19




Top 50 High-tech Enterprises in Zhuhai Xiangzhou District



Small and Medium-sized Growth Project Cultivation Enterprise in Zhuhai Xiangzhou District

INTRODUCTION TO CUSTOMIZED SERVICES

• R & D platform



• Customer groups

Biori has been focusing on the research and development of nucleic acid diagnostic raw materials for nearly ten years, and has accumulated rich experience. We can provide customers with a series of solutions in various related fields, like technical demonstration to R&D projects, research optimization, pilot-scale development, and production scale-up.




• Customized product

Biori's products, such as nucleic acid diagnostic enzymes, high-performance premix, direct amplification reagents, rapid amplification reagents, isothermal amplification reagents, extraction reagents, instruments and equipment, can provide customers with OEM and flexible customized development services.



• Service process



CONTENTS

01 LYO-READY RAW MATERIALS 04

Lyo-ready raw materials series results display	05
--	----

02 RAW MATERIALS FOR NUCLEIC ACID DETECTION 12

Novel hot-start Taq DNA Polymerase	13
--	----

NEW
DNA one-tube reagent

MMLV reverse transcriptase	13
----------------------------------	----

Tth DNA Polymerase	14
--------------------------	----

Direct amplification reagents	14
-------------------------------------	----

Fast Amplification Premix For Purified Template	15
---	----

NEW PRODUCTS 20

Universal nucleic acid release reagent Stool Direct Premix RNA one-tube reagent

Isothermal amplification	25
--------------------------------	----

NEW
Bst 2.0 HS

Other products	34
----------------------	----

03 INSTRUMENTS & REAGENTS 35

Centrifugal mixer	36
-------------------------	----

CM-8 / plus

Full automatic nucleic acid extractor	37
---	----

NEW
BTE-8 BTE-32

Sample pretreatment products NEW	38
---	----

Nucleic acid extraction kit (magnetic bead) Release reagent

Sample preservation solution



PART 01

LYO-READY RAW MATERIALS




Lyo-ready raw materials series results display

Lyophilized formulations are available for the all products

Varirty	Diverse	Advantage
Full-line PCR test materials support lyophilized formula	PCRTube (8-Tube Strips), lyophilized microsphere (sphere), Vial	Full stability verification, a variety of products can be customized, there are certified customers.

Since 2014, Biori has been developing lyophilized raw materials for molecular diagnostic reagents, and has taken the lead in developing lyophilized enzymes and premix for DNA and RNA detection in China. At present, Biori has a full set of lyophilized raw materials suitable for fluorescence quantitative detection of DNA and RNA. Including enzymes for lyophilization, premix, lyophilization protectants, rapid and direct amplification reagents for lyophilization, loop-mediated isothermal amplification reagents for lyophilization and many other series of lyophilization raw materials products.

A large number of long-term accelerated stability tests have been carried out on the lyophilized materials, which can provide stronger support for the research and development of lyophilized reagents. At present, the lyophilized reagent products of our customers are the first to obtain NMPA registration approval.



During the cooperation, the client only needs to provide relevant primers or probes to complete the feasibility study of the lyophilization project. Biori will help the client to optimize the lyophilized procedure and guide the client to perform the lyophilized operation.

In order to cooperate with the customer's research and development work, Biori also undertakes the lyophilized service. Available lyophilization forms include PCR Tube(8-Tube Strips), lyophilized microsphere(sphere) and Vial, to help customers quickly obtain the most reliable project results.

■ Stability data of lyophilized product in 8-Strip tubes



01 DNA amplification reagent lyophilized product

Test 1: Accelerated stability (45°C) verification of human genome quadruple detection system


Experimental conditions:

Code	FM5071	System	Human genome quadruple detection system
Template	Human genome template	Template concentration	2.5、25、250 pg/μL
Cycling program	50°C 2min;95°C 5min;50 Cycles (95°C 10s, 55°C 40s)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box, and accelerated at 45°C.

Experimental results:

45°C for 90 days▼



Test 2: Accelerated stability (45°C) verification of human genome quadruple detection system


Experimental conditions:

Code	FM5071	System	Human genome quadruple detection system
Template	Human genome purification template	Template concentration	2.5、25、250 pg/μL
Cycling program	50°C 2min;95°C 5min;50 Cycles (95°C 10s, 55°C 40s)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box, Place at room temperature.

Experimental results:

Room temperature for 1 year▼



02 RNA amplification reagent lyophilized products

Test 1: Accelerated Stability (55°C) Verification of Human Total RNA Purification Template Triple Detection System


Experimental conditions:


Code	FM5134	System	Triple detection system of human total RNA
Template	Human total RNA purification template	Template concentration	0.464、4.64、46.4 pg/μL
Cycling program	50°C 15min; 95°C 1min; 50 Cycles (95°C 15s, 56°C 45s)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box, and accelerated at 55°C.

Experimental results:

55°C for 30 days▼





Test 2: Room temperature Stability Verification of Human Total RNA Purification Template Triple Detection System




Experimental conditions:

Code	FM5134	System	Triple detection system of human total RNA
Template	Human total RNA purification template	Template concentration	0.464, 4.64, 46.4 pg/μL
Cycling program	50°C 15min; 95°C 1min; 50 Cycles (95°C 15s, 56°C 45s)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box, Place at room temperature.

Experimental results:

Room temperature for 1 year▼



03 Loop-mediated Isothermal Amplification (LAMP)

Test 1: Accelerated stability validation of LAMP lyophilized system (accelerated at 55°C for 30 days)


Experimental conditions:

Code	HW201-R01	System	LAMP lyophilized system
Template	λDNA	Template concentration	0.05pg/T, 0.5pg/T, 5pg/T
Cycling program	65°C, 60min(Fluorescence was collected every minute)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box and accelerated at 55°C.

Experimental results:

55°C for 30 days▼



Note: The amount of template addition from left to right is 5pg/T, 0.5pg/T, 0.05pg/T.

Test 2: Accelerated stability validation of RT-LAMP lyophilized system (accelerated at 55°C for 30 days)


Experimental conditions:

Code	HW202-R01	System	RT-LAMP lyophilized system
Template	Human total RNA	Template concentration	1ng/T, 10ng/T, 100ng /T
Cycling program	65°C, 60min(Fluorescence was collected every minute)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box and accelerated at 55°C.

Experimental results:

55°C for 30 days▼



Note: template addition from left to right: 100ng/T, 10ng/T, 1ng/T.

Stability data of lyophilized products in vials



RNA amplification reagent lyophilized product

Test: Accelerated stability (55°C) validation of triple detection system for human total RNA purification template




Experimental conditions:

Code	FM5134	System	Triple detection system for human total RNA
Template	Human total RNA purification template	Template concentration	0.464, 4.64, 46.4 pg/μL
Cycling program	50°C 15min; 95°C 1min; 50 Cycles (95°C 15s, 56°C 45s)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box and accelerated at 55°C.

Experimental results:

55°C for 45 days▼



Stability data of lyophilized microspheres



RNA amplification reagent lyophilized product

Test: Accelerated stability (55°C) validation of triple detection system for human total RNA purification template




Experimental conditions:

Code	FMD5084	System	Triple detection system for human total RNA
Template	Human total RNA purification template	Template concentration	0.464、4.64、46.4 pg/μL
Cycling program	50°C 15min; 95°C 1min; 50 Cycles (95°C 15s, 56°C 45s)		

Note: In the control group, the liquid reagent of the same batch of lyophilized reagent was prepared immediately. The lyophilized reagent was vacuum-packed in aluminum foil bag after leaving the box and accelerated at 55°C.

Experimental results:

55°C for 30 days▼





PART 02

RAW MATERIALS FOR NUCLEIC ACID DETECTION



- Novel hot-start Taq DNA Polymerase
- MMLV reverse transcriptase
- Tth DNA Polymerase
- Direct amplification reagents
- Fast Amplification Premix For Purified Template
- New Products
- Isothermal amplification
- Other products


All products can be customized / have anti-pollution system

The following icons (without text description) will appear below the product name to annotate its features:

- 1 ▶ Suitable for qPCR 2 ▶ Suitable for qRT-PCR 3 ▶ Lyophilizable 4 ▶ Can be used with anti-pollution system
 5 ▶ Customizable concentration 6 ▶ Suitable for multiplex amplification 7 ▶ Stable preserve in one tube

Novel hot-start Taq DNA polymerase and its relevant premix


Biori hot-start Taq, novelty modification, has the characteristics of high specificity, high sensitivity, and complete categories. It is suitable for all kinds of PCR detection methods

Product name	Code	Concentration	Characteristic	Recommended Application
HS Taq ① ④ ⑤ ⑥	E03	5U/μL	①Chemically-modified and incubate at 95°C for 10 min for polymerase activation. ②The activity of polymerase is gradually released, and performs high specificity.	①End-Point PCR and quantitative PCR. ②Multiplex PCR and genotyping test.
2× HS Premix (Probe qPCR) ① ④ ⑤ ⑥	M203	\		
Superstart Taq plus ① ② ③ ④ ⑤ ⑥	E07	5U/μL	①Modified by antibody and hot start at 95°C for 1~5 min. ②High sensitivity and specificity. ③Stable detection of low concentration samples, and performs high fluorescence value.	①End-Point PCR and quantitative PCR. ②Multiplex PCR, genotyping test and high sensitivity detection of viral nucleic acid.
2× Superstart Premix plus (Probe qPCR) ① ② ④ ⑤ ⑥	M207	\		
Hyperstart® Taq ① ② ③ ④ ⑤ ⑥	E09	5U/μL	①Modified by aptamer, and reversible hot start above 60°C. ②Dynamically modified during PCR process, with high sensitivity and specificity. ③ It can be stored stably at 2~8°C.	①End-Point PCR and quantitative PCR. ②Multiplex PCR, genotyping test and high sensitivity detection of viral nucleic acid.
2× Hyperstart® Premix (Probe qPCR) ① ② ④ ⑤ ⑥	M209	\		
Robustart Taq ① ② ③ ④ ⑤ ⑥	E16	5U/μL	①Modified by antibody and hot start at 95°C for 1~5 min. ②High sensitivity and specificity. ③It has strong applicability to various detections.	①End-Point PCR and quantitative PCR. ②Multiplex PCR, genotyping test and high sensitivity detection of viral nucleic acid.
2× Robustart Premix (Probe qPCR) ① ② ④ ⑤ ⑥	M216	\		
2× Robustart Premix OT (Probe qPCR) ① ② ④ ⑤ ⑥ 	M2211	5U/μL	①Modified by antibody and hot start at 95°C for 1~5min. ②High stability under one-tube conditions (enzyme, buffer and primer probe are mixed in one tube).	①Detection of human and animal-related DNA viruses. ②POCT
Robustart Taq QS II ① ② ③ ④ ⑤ ⑥	E20	5U/μL	①Modified by antibody and hot start at 95°C for 1~5 min. ②Amplification rate of rapid DNA polymerase is no less than 1kb/10s. ③It exhibits extremely high resistance to inhibitors found in blood, swabs, etc.	①End-Point PCR and quantitative PCR. ②Rapid amplification, direct amplification and POCT nucleic acid detection.

MMLV Reverse Transcriptase and its relevant premix

Biori MMLV Reverse Transcriptase has the characteristics of high efficiency reverse transcription, good stability, excellent performance and strong applicability. It can match all kinds of amplification enzymes and is suitable for all kinds of RT-PCR detection methods.

Product name	Code	Concentration	Characteristic	Recommended Application
Celascript® RTase ② ③ ④ ⑤ ⑥	E12	200U/μL	①Reverse transcription temperature range from 42°C to 50°C. ②It has high efficiency in cDNA synthesis. ③Strong substrate affinity.	①cDNA synthesis(E12), quantitative RT-PCR. ②Detection of human and animal RNA viruses. High sensitivity detection of novel coronavirus, influenza virus and respiratory tract virus in multiplex detection.
5× Celascript® RT Premix (Probe qRT-PCR) ② ③ ④ ⑤ ⑥	M512	\		
Neoscript® RTase ② ③ ④ ⑤ ⑥	E13	200U/μL	①Reverse transcription temperature range from 42°C to 55°C. ②Excellent specificity and sensitivity. ③High applicability, and suitable for all kinds of high sensitivity detection.	①cDNA synthesis(E13), quantitative RT-PCR. ②Detection of human and animal RNA viruses. High sensitivity detection of novel coronavirus, influenza virus and respiratory tract virus in multiplex detection.
5× Neoscript® RT Premix (Probe qRT-PCR) ② ③ ④ ⑤ ⑥	M513	\		

Product name	Code	Concentration	Characteristic	Recommended Application
Neoscript® RT Premix OT (Probe qRT-PCR) ② ③ ④ ⑤ ⑥ ⑦ 	M527	\	① Reverse transcription temperature range from 42°C to 55°C. ② High stability under one-tube conditions (enzyme, buffer and primer probe are mixed in one tube). ③ High applicability, suitable for all kinds of high sensitivity detection.	① Detection of human and animal RNA viruses. High sensitivity detection of novel coronavirus, influenza virus and respiratory tract virus in multiplex detection. ② POCT.

Tth Polymerase and its relevant premix



Biori Tth polymerase with the functions of RNA reverse transcription and DNA amplification, the inhibitor is well tolerated. We provide a variety of modification schemes for it and it is suitable for various PCR and RT-PCR detection methods.

Product name	Code	Concentration	Characteristic	Recommended Application
Hyperstart® Tth ① ② ③ ④ ⑤ ⑥	E10	5U/μL	① Modified by aptamer, and reversible hot start above 60°C. ② It has both polymerase and reverse transcriptase activities. ③ In the presence of Mn ²⁺ , the reverse transcription can be carried out at a high temperature from 60 to 70°C. The reverse transcription efficiency of complex RNA templates is higher and non-specificity is reduced. ④ It preforms high inhibitor-tolerance.	① Quantitative PCR. ② Multiple PCR, Genotyping, Highly sensitive virus detection. ③ Preserve at 2-8°C.
5× Hyperstart® Tth Premix (Probe qRT-PCR) ① ② ③ ④ ⑤ ⑥	M510	\		
Robustart Tth ① ② ③ ④ ⑤ ⑥	E17	5U/μL	① Modified by antibody and hot start at 95°C for 1~5 min. ② It has both polymerase and reverse transcriptase activities. ③ In the presence of Mn ²⁺ , the reverse transcription can be carried out at a high temperature from 60 to 70°C. The reverse transcription efficiency of complex RNA templates is higher and non-specificity is reduced. ④ It preforms high inhibitor-tolerance.	① Quantitative PCR. ② Multiple PCR, genotyping, and highly sensitive virus detection.

Direct Amplification Premix

Biori Direct Amplification Premix use innovative and unique formulations containing antibody-modified polymerase, have the characteristics of high sensitivity and inhibitor-tolerance. It can be lyophilized and used for in-situ lysis and amplification of biological samples, and ensure the accuracy of results. It is suitable for extraction-free nucleic acid detection of blood, swabs and other samples.

Product name	Code	Characteristic	Recommended Application
2× Superstart™ Direct Premix (Probe qPCR) ① ③ ④ ⑤ ⑥	MD201	① It exhibits extremely high resistance to inhibitors found in blood, swabs, tissue homogenate, etc. The recommended volume of different samples is whole blood ≤5%, plasma ≤30% and serum ≤30%. ② It is suitable for direct fluorescence quantitative PCR detection of DNA without extraction.	① Human and animal related DNA virus detection. ② Genomic DNA amplification and genotyping.
2× Hyperstart® Direct Premix (Probe qPCR) ① ③ ④ ⑤ ⑥	MD202	① It exhibits extremely high resistance to inhibitors found in blood, swabs, tissue homogenate, etc. The recommended volume of different samples is whole blood ≤5%, plasma ≤30% and serum ≤30%. ② It is suitable for direct fluorescence quantitative PCR detection of DNA without extraction.	① Human and animal related DNA virus detection. ② Genomic DNA amplification and genotyping.
2× SensiDirect® Premix (Probe qPCR) ① ③ ④ ⑤ ⑥	MD209	① It exhibits extremely high resistance to inhibitors found in blood, swabs, tissue homogenate, etc. The recommended volume of different samples is whole blood ≤5%, plasma ≤30% and serum ≤30%. ② It is suitable for direct fluorescence quantitative PCR detection of DNA without extraction.	① Human and animal related DNA virus detection. ② Genomic DNA amplification and genotyping.
2× Fast Direct Premix (Probe qPCR) ① ③ ④ ⑤ ⑥	MD207	① The amplification rate is not less than 1kb/5s, which is suitable for rapid extraction-free fluorescence quantitative PCR, and the detection can be completed within 30min. ② It exhibits extremely high resistance to inhibitors found in blood, swabs, tissue homogenate, etc. The recommended volume of different samples is whole blood ≤5%, plasma ≤30% and serum ≤30%.	① Human and animal related DNA virus detection. ② Genomic DNA amplification and genotyping.

Product name	Code	Characteristic	Recommended Application
2×SensiDirect RT Premix (Probe qRT-PCR) ② ③ ④ ⑤ ⑥	MD210	①It exhibits extremely high resistance to inhibitors found in blood, swabs, tissue homogenate, etc. The recommended volume of different samples is whole blood ≤5%, plasma ≤30% and serum ≤30%. ②It is suitable for direct fluorescence quantitative PCR detection of DNA without extraction.	①Human and animal related RNA virus detection. ②Direct amplification detection without extraction.
2×Fast Direct RT Premix (Probe qRT-PCR) ② ③ ④ ⑤ ⑥	MD208	①It is suitable for rapid extraction-free fluorescent quantitative RT-PCR, and the detection can be completed within 40min. ②It exhibits extremely high resistance to inhibitors found in blood, swabs, tissue homogenate, etc. The recommended volume of different samples is whole blood ≤5%, plasma ≤30% and serum ≤30%.	①Human and animal related RNA virus detection. ②Direct amplification detection without extraction.
Stool Direct Premix (Probe qPCR) ② ③ ④ ⑤ ⑥ 	MD211	①It has a high tolerance for samples derived from stool and anal swabs, and the proportion of stool directly added to the sample is ≤0.5%. ②Suitable for DNA extraction-free direct fluorescence qPCR detection of stool samples.	① Detection of DNA viruses and pathogens related to humans and animals. ②Amplification and typing of genomic DNA.
通用型核酸释放剂 ① ② ⑤ 	AS17	① Pre-treatment of biological samples, such as swabs, anticoagulant whole blood, stool, etc., swab sample ≤ 20%, anticoagulation whole blood sample ≤ 20%, stool fluid ≤20%. ②Efficient and time-saving, one-step cracking and binding, instant release, no heating and no centrifugation during processing.	① Detection of DNA/RNA viruses and pathogens related to humans and animals. ②Fast amplification, POCT.

Fast Amplification Premix For Purified Template

Product name	Code	Characteristic	Recommended Application
2×Superstart Fast Premix (Probe qPCR) ① ③ ④ ⑤ ⑥	M222	①The amplification rate is not less than 1kb /10s, which is suitable for rapid fluorescent quantitative PCR of purified DNA template, and the reaction can be completed within 30min. ② It has excellent specificity and sensitivity.	①Human and animal related DNA virus detection. ②Genomic DNA amplification and genotyping. ③Rapid amplification detection.
5×Neoscript® Fast RT Premix (Probe qRT-PCR) ② ③ ④ ⑤ ⑥	M524	①The amplification rate is not less than 1kb /10s, which is suitable for rapid fluorescent quantitative PCR of purified RNA template, and the reaction can be completed within 30min. ② It has excellent specificity and sensitivity.	①Human and animal related RNA virus detection. ②Genomic DNA amplification and genotyping. ③Rapid amplification detection.



01

Testing data of direct DNA amplification reagent

Experimental conditions:

Code	MD2011	System	HBV, human genomic DNA detection system
Template	HBV-positive plasma, whole blood, Oral swabs	Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)

● Direct DNA amplification of plasma

HBV-positive plasma: 48, 96, 1.6×10^3 , 1.6×10^4 IU/mL. All concentrations of HBV positive plasma were directly added into the reaction system for amplification without pretreatment. And plasma is added 25% of the total PCR volume into the reaction system before amplification, 6.25μL HBV-positive plasma was added into 25μL reaction system.




Fig.A Amplification curve of HBV





Fig.B Amplification curve of Internal standard


● Direct DNA amplification of human whole blood

Human whole blood: Human whole blood was directly added into the reaction system for amplification without pretreatment. And human whole blood is added 5% of the total PCR volume into the reaction system before amplification, 1.25μL Human whole blood was added into 25uL reaction system.



● Direct DNA amplification of Oral swabs

Oral swabs: After brushing in the oral cavity for about 10 times, we wash the swab with 1000μL H₂O, centrifuge at 2500rpm for 10min and discard the supernatant, then there is about 100μL remaining liquid, blow and mix well the remaining liquid, take 10μL and add into reaction system for amplification.




02 Testing data of amplification efficiency

Experimental conditions:

Code	MD2091	System	HBV detection system
Template	HBV-positive plasma, which added 25% of the total PCR volume into the reaction system	Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)

● Direct DNA amplification of plasma

HBV-positive plasma: 1.6×10^3 、 1.6×10^4 、 1.6×10^5 、 1.6×10^6 、 1.6×10^7 IU/mL. All concentrations of HBV positive plasma were directly added into the reaction system for amplification without pretreatment. And plasma is added 25% of the total PCR volume into the reaction system before amplification



03 Testing data of amplification comparison between biological samples without extraction and purified biological samples after extraction


Experimental conditions:

Code	MD2091	System	HBV detection system
Template	HBV-positive plasma, HBV-positive serum, Magnetic bead extraction HBV positive plasma template	Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)

Magnetic bead extraction HBV positive plasma template: 200μL HBV positive plasma was extracted and eluted with 85μL 1× TE Buffer. And template is added 40% of the total PCR volume into the reaction system before amplification.

HBV-positive plasma: 1.6×60 IU/mL, 4T; 1.6×10^3 IU/mL, 2T; 1.6×10^4 IU/mL, 2T. All concentrations of HBV positive plasma were directly added into the reaction system for amplification without pretreatment. And plasma is added 25% of the total PCR volume into the reaction system before amplification.

HBV-positive serum: 1.6×60 IU/mL, 4T; 1.6×10^3 IU/mL, 2T; 1.6×10^4 IU/mL, 2T. All concentrations of HBV positive serum were directly added into the reaction system for amplification without pretreatment. And plasma is added 25% of the total PCR volume into the reaction system before amplification




Template	Template concentration (IU/mL)	Equivalent template concentration (IU/T)	Average Ct value	Template	Template concentration (IU/mL)	Equivalent template concentration (IU/T)	Average Ct value
Magnetic bead extraction HBV positive plasma template	1×10^4	235	28.72	HBV-positive plasma	1.6×10^4	100	30.17
	1×10^3	23.50	31.90		1.6×10^3	10	33.45
	50	1.20	35.78		1.6×60	0.6	37.90
HBV-positive serum	1.6×10^4	100	29.87				
	1.6×10^3	10	33.17				
	1.6×60	0.60	38.68				

04

Testing data of DNA rapid direct amplification product

Experimental conditions:

Code	MD2071	Template concentration	0.6、1、10、100IU/T
Template	HBV-positive plasma	Cycling program	50°C 2min; 95°C 1min; 50 Cycles (95°C 1s, 60°C 5s); Complete the amplification within 23min
Sample volume	Added 25% of the total PCR volume into the reaction system before amplification		



05


Testing data of RNA direct amplification product

Experimental conditions:

Code	MD 2104	Cycling program	50°C 15min; 95°C 1min; 50 Cycles (95°C 15s, 56°C 45s)
Template	HCV-positive plasma (Added 20% of the total PCR volume into the reaction system), Magnetic bead extraction HCV positive plasma template		

HCV-positive plasma: 1.3×10^3 U/mL, 6T; 2.61×10^4 IU/mL, 2T. All concentrations of HCV positive serum were directly added into the reaction system for amplification without pretreatment. And plasma is added 20% of the total PCR volume into the reaction system before amplification.


Magnetic bead extraction HCV positive plasma template: 200μL HCV positive plasma (130 IU/mL, 6T; 2.61×10^3 IU/mL, 2T) was extracted and eluted with 85μL $1 \times$ TE Buffer. And template is added 40% of the total PCR volume into the reactionsystem before amplification.




Template	Template concentration (IU/mL)	Equivalent template concentration (IU/T)	Average Ct value	Template	Template concentration (IU/mL)	Equivalent template concentration (IU/T)	Average Ct value
Magnetic bead extraction HCV positive plasma template	2.61×10^3	61.40	31.51	HCV- positive plasma	2.61×10^4	130.50	30.61
	1.30×10^2	3.06	35.97		1.30×10^3	6.50	35.34

06 Comparison with similar products

HCV-positive plasma: 1×10^4 IU/mL, 2T; 1×10^5 IU/mL, 2T. All concentrations of HCV positive serum were directly added into the reaction system for amplification without pretreatment. And plasma is added 20% of the total PCR volume into the reaction system before amplification.




Saliva template(Human total RNA amplification): 10% of the total PCR volume, 2T. 20% of the total PCR volume, 2T. Saliva template was directly added into the reaction system for amplification without pretreatment.



07 Testing data of RNA rapid direct amplification product

Experimental conditions:

Code	MD2084	Template concentration	2.5、5、50、500IU/T
Template	HCV-positive plasma	Cycling program	50°C 7min; 95°C 1min; 50 Cycles (95°C 1s, 56°C 10s) Complete the amplification within 23min
Sample volume	Added 20% of the total PCR volume into the reaction system before amplification		



NEW PRODUCTS

Universal nucleic acid release reagent

Stool Direct Premix

RNA one-tube reagent

NEW Universal nucleic acid release reagent

LYO-READY RAW
MATERIALS

RAW MATERIALS FOR
NUCLEIC ACID DETECTION


INSTRUMENTS
& REAGENTS

01 Comparison test of similar products (anticoagulant human whole blood)

Experimental conditions:

Code	AS17	Detection reagent	MD2091
Template	EDTA Anticoagulated human whole blood	Template amount	5μL/T
Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)		

Whole blood sample processing: Take 20 μL of whole blood and 60 μL of release reagent into an Ep tube, pat it slowly for 15-20 times to avoid the air bubbles, and add the sample directly after mixing for extraction-free nucleic acid amplification detection.




02

Pseudovirus sample processing effect test

Experimental conditions:

Code	AS17	Detection reagent	MD2104
Template	RBCS pseudovirus	Template concentration	1.4×10^5 , 1.4×10^6 , 1.4×10^7 and 1.4×10^8 copies/ μ L
Cycling program	50°C 15min; 95°C 1min; 50cycles (95°C 15s, 56°C 45s)		

Pseudovirus sample processing: Take 30 μ L of RBCS pseudovirus (gradient concentration) and 30 μ L of release reagent in an Ep tube, pat it slowly for 15-20 times to avoid the air bubbles, and add the sample directly after mixing for extraction-free nucleic acid amplification detection. Treated with 1 \times TE Buffer as control.







03

Stool sample treatment effect test

Experimental conditions:

Code	AS17	Detection reagent	MD2084-4
Template	Human genome in stool samples	Template concentration	0.0125 ng/ μ L
Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)		

Stool sample processing: Take 0.2 g of stool sample and add it to 1 mL of TE Buffer, mix and liquefy well; take the supernatant after instant centrifugation to obtain 20% stool sample, which is diluted with sterile water to obtain 2% stool samples; mix 2% stool solution and release reagent in equal volume, and treat with sterile water as a control; dilute the human genome with the above-treated samples to obtain a human genome template with a concentration of 0.0125 ng/ μ L; then 20% sample was directly added to the reaction system for amplification.




04 Accelerated Stability Verification

Experimental conditions:

Code	AS17	Detection reagent	MD2091
Template	Anticoagulated human whole blood	Template amount	5 μ L/T
Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)		

Test method: Mix the nucleic acid release reagent, divide it into two parts, put them at -20°C and 37°C respectively, and take them out on the 30th day. Take 20 μ L of whole blood and 60 μ L of release reagent into an Ep tube, pat it slowly for 15-20 times to avoid the air bubbles, and add the sample directly after mixing for extraction-free nucleic acid amplification detection.

Whole blood sample processing: Take 20 μ L of whole blood and 60 μ L of release reagent into an Ep tube, pat it slowly for 15-20 times to avoid the air bubbles, and add the sample directly after mixing for extraction-free nucleic acid amplification detection.




NEW Stool Direct Premix Comparison of similar products

01 Salmonella detection in stool samples

Experimental conditions:

Code	MD2111	Detection reagent	Salmonella detection system
Template	Salmonella in stool samples	Template concentration	1.496, 14.96, and 149.6 copies/ μ L
Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)		

Stool sample processing: Take 0.2 g of stool samples and add them to 1 mL of TE Buffer, fully mix and liquefy; take the supernatant after instant centrifugation to obtain 20% stool samples, which are diluted with sterile water to obtain 1% stool Samples; serially dilute 10-fold 1% fecal samples of Salmonella to make the Salmonella concentrations at 149.6, 14.96 and 1.496 copies/ μ L respectively; then add 20% of the sample volume directly to the reaction system







02

Stool Background Human Genome Template Quadruple Assay

Experimental conditions:

Code	MD2111	System	Human Genomic DNA Quadruple Detection System
Template	Human genome in stool samples	Template concentration	0.0125 ng/ μ L
Cycling program	50°C 2min; 95°C 5min; 50 Cycles (95°C 10s, 55°C 40s)		

Stool sample processing: Take 0.2 g of stool sample and add it to 1 mL of TE Buffer, mix and liquefy well; take the supernatant after instant centrifugation to obtain 20% stool sample, which is diluted with sterile water to obtain 1% and 2.5% stool samples; The human genome was diluted with the above two concentrations of stool samples to obtain a human genome template at a concentration of 0.0125 ng/ μ L; then 20% sample was directly added to the reaction system for amplification.









NEW RNA one-tube reagent

01 Comparison test of similar products

Experimental conditions:

Code	M5274	System	Respiratory Virus Triple Detection System
Template	Respiratory virus quality control sample purification template	Template concentration	0.25、2.5 copies/μL
Cycling program	50°C 15min; 95°C 1min; 50 Cycles (95°C 15s, 56°C 45s)		

Test method: Mix buffer, enzyme and primer probe to prepare working solution, and place it at 37°C for 3 days. The comparison test was conducted with the instant preparation working solution of the same system as a control.



NEW PRODUCTS

Bst 2.0 HS

Loop-mediated Isothermal Amplification (LAMP)

Enzymes	
Product Name	Code
Bst 2.0 (8U/μL)	E102
Bst 2.0 HS (8U/μL)	E103
Neoscript RTase (200U/μL)	E13
RNase H II	E201
DNA-Lamp NEW	
Product Name	Code
2×Lamp Premix (Bst 2.0 HS)(Eva Green/SYBR Green)	HW205-R01
2×Lamp Premix (Bst 2.0 HS)(Probe)	HW205-P01
Colorimetric Lamp Kit (Bst 2.0 HS)	HW205-M01
Colorimetric Lamp Kit (Bst 2.0 HS)	HW205-M02
Colorimetric Lamp Kit (Bst 2.0 HS)	HW205-M03
RNA-Lamp NEW	
Product Name	Code
2×RT Lamp Premix (Bst 2.0 HS)(Eva Green/SYBR Green)	HW206-R01
2×RT Lamp Premix (Bst 2.0 HS)(Probe)	HW206-P01
Colorimetric RT Lamp Kit (Bst 2.0 HS)	HW206-M01
Colorimetric RT Lamp Kit (Bst 2.0 HS)	HW206-M02
Colorimetric RT Lamp Kit (Bst 2.0 HS)	HW206-M03


M01:Cresol red staining method, The color is purple before reaction. The positive reaction turns to yellow, while the negative reaction will remain the original color.

M02: Calcein staining method, The color is light orange before reaction. The positive reaction turns to green, while the negative reaction will remain the original color.

M03: Hydroxyphenol blue, the color is violet before the reaction, the positive reaction turns to blue while the negative reaction will remain violet.

The above products can be provided in lyophilized version, please inquire for the product code.

Bst 2.0 / Bst 2.0 HS



Biori has obtained Bst2.0 through a high-throughput screening platform. The enzyme has been well received by customers for its excellent sensitivity and good resistance to inhibitors. At the same time, through the research and development of the protein modification technology platform, our hot-start version Bst2.0 was launched firstly in domestic market, further improves the specificity and sensitivity of LAMP reagents. On this basis, Bioi has developed a series of various LAMP and RT-LAMP isothermal amplification reagents such as colorimetric version, fluorescent dye method, and probe method. There are lyophilized version for the above-mentioned raw materials.

Bst2.0HS is a hot-start isothermal polymerase obtained by reversible modification technology on the basis of Bst2.0 DNA polymerase. It can completely block the activity of the enzyme at room temperature and establish a reaction at room temperature with preventing non-specificity and improving reaction efficiency. In addition, Bst2.0HS DNA polymerase does not require a separate activation step.


◆ 37°C accelerated stability data of lyophilized reagent raw materials in Bst2.0HS system

When the lyophilized reagent raw materials Bst2.0HS, Neoscript RTase and RNaseH II were separately accelerated at 37°C for 14 d, the performance of the whole reaction system remained unchanged.

Experimental conditions:

Code	FHW206-P01 (Bst 2.0 HS)	System	RT-LAMP (Probe) 25μL/T
Template	Human total RNA(Detection of β-actin gene)	Template amount	0.01ng/T
Cycling program	65°C 60min; 60 Cycles		

Experimental results:



◆ Performance after modification


1.1 Specificity test

Biori BST2.0HS can perfectly replace imported products (NEB Bst 2.0 WarmStart DNA Polymerase), Biori BST2.0HS can achieve the same specificity as competing products, and performs better than the conventional Bst DNA polymerase.

Experimental conditions:

Code	HW205-R01 (Bst 2.0 HS)	System	Real-time Fluorescent LAMP Kit
Template	Human genome (GAPDH gene)	Template amount	0.1ng/T(Left) NTC(Right)
Cycling program	65°C 60min; 60 Cycles		

Experimental results:



1.2 Room temperature stability test


The Biori Bst 2.0 HS enables room temperature operation.

Experiments show that after adding template to prepare a mixed reaction system, under two conditions of incubation at 25°C for 3 hours (green) and immediate reaction after preparation (blue), Biori Bst 2.0 HS has similar reaction results compared with control Bst 2.0 ws. And the reaction performance of the Bst enzyme without hot-start modification is reduced after being incubated at room temperature for 3 hours, and the hot-start modified Bst enzyme can be equivalent to the performance of instant preparation, effectively maintaining the stability of the reaction performance.

Experimental conditions:

Code	HW205-R01 (Bst2.0 HS)	System	LAMP (Eva green) 25μL/T
Template	Human genome (GAPDH gene)	Template amount	1ng/T
Cycling program	65°C, 60min, 60 Cycles		

Experimental results:




◆ Conventional Lamp performance comparison test

2.1 Performance comparison of DNA purified templates

Experimental conditions:

Code	HW205-R01 (Bst2.0 HS)	System	LAMP (Eva green) 25μL/T
Template	Human genome (GAPDH gene primer1, left; GAPDH gene primer2, right)	Template amount	0.1ng/T
Cycling program	65°C, 60min, 60cycles		

Experimental results:




2.2 Performance comparison of RNA purified templates

Experimental conditions:

Code	HW206-R01 (Bst2.0 HS)	System	RT-LAMP (Eva green) 25μL/T
Template	Human total RNA (β-actin gene primer1, left R18S12 gene primer2, right)	Template amount	0.01ng/T
Cycling program	65°C, 60min 60 cycles		

Experimental results:




◆ Performance comparison of Direct premix

3.1 Performance in whole blood sample

Experimental conditions:

Code	HW206-R01 (Bst2.0 HS)	System	RT-LAMP (Eva green) 25μL/T
Template	Human whole blood (β-actin gene)	Template amount	5%
Cycling program	65°C, 60min, 60 cycles		

Experimental results:




3.2 Performance in oral swab sample

Experimental conditions:

Code	HW206-R01 (Bst2.0 HS)	System	RT-LAMP (Eva green) 25μL/T
Template	Human oral swab (β-actin gene)	Template amount	Swab solution 5μL/T
Cycling program	65°C, 60min, 60 cycles		

Experimental results:




3.3 Performance in saliva sample

Experimental conditions:

Code	HW206-R01 (Bst2.0 HS)	System	RT-LAMP(Eva green)25μL/T
Template	Human saliva (β-actin gene)	Template amount	3μL/T
Cycling program	65°C, 60min, 60cycles		

Experimental results:




◆ Corona virus testing(Probe method) ➦

Experimental conditions:

Code	HW206-P01 (Bst2.0 HS)	System	RT-LAMP (Probe method) 25μL/T
Template	Extracted nucleic acid of corona virus	Template amount	About 31.25 copies/T (red) About 15.63 copies/T (blue) About 7.82 copies/T (purple) NTC (grey)
Cycling program	65°C, 60min, 60 cycles		

Experimental results:




◆ Corona virus testing(Eva green) ❷

Experimental conditions:

Code	HW206-P01 (Bst 2.0 HS)	System	RT-LAMP (Eva green) 25μL/T
Template	Extracted nucleic acid of corona virus	Template amount	About 31.25 Copies/T (red) About 15.63 Copies/T (blue) About 7.82 Copies/T (purple) NTC (grey)
Cycling program	65°C, 60min, 60 cycles		

Experimental results:





◆ Testing result of Colorimetric LAMP Kit ❷

6.1 Colorimetric LAMP Kit (RY)

Experimental conditions:

Code	HW206-M01 (Bst 2.0 HS)	System	Colorimetric RT-LAMP (RY)
Template	Human total RNA	Template amount	0.01ng/T、0.1ng/T、1ng/T
Cycling program	65°C, 60min ,60 cycles		

Experimental results:



	Tube5: Negative; Tube 6: 0.01ng/T; Tube 7: 0.1ng/T; Tube 8: 1ng/T							
Before								
After								

6.2 Colorimetric LAMP Kit(OG):

Experimental conditions:

Code	HW206-M02 (Bst 2.0 HS)	System	Colorimetric RT LAMP Kit(OG)
Template	Human total RNA	Template amount	0.01ng/T、0.1ng/T、1ng/T
Cycling program	65°C, 60min ,60 cycles		

Experimental results:



	Tube5:Negative; Tube 6:0.01ng/T; Tube 7:0.1ng/T; Tube 8:1ng/T							
Before								
After								

6.3 Colorimetric LAMP Kit (HNB):

Experimental conditions:

Code	HW206-M03 (Bst 2.0 HS)	System	Colorimetric RT LAMP Kit (HNB)
Template	Human total RNA	Template amount	0.01ng/T、0.1ng/T、1ng/T
Cycling program	65°C, 60min ,60 cycles		

Experimental results:

	Tube 5: Negative; Tube 2: 0.01ng/T; Tube 3: 0.1ng/T; Tube 4: 1ng/T
Before	
After	

◆ Other Products

01. Reagents using dye-binding method

Product Name	Code	Characteristics	Recommended Application
Accustart Premix (SYBR qPCR) ① ③ ④ ⑤	M203S	① Contained chemical modified enzyme, which activated at 95°C for 5~10 minutes. ② Gradually release enzyme activation; high specificity.	① qPCR. ② High sensitivity virus detection.
Robustart Premix (SYBR qPCR) ① ③ ④ ⑤	M216S	① Contained antibody modified enzyme, which activated at 95°C for 1~5 minutes. ② High sensitivity and specificity. ③ High applicability in various detection.	① qPCR. ② High sensitivity virus detection.
Novel Premix (SYBR qPCR) ① ④ ⑤	PM202S	① Contained antibody modified enzyme, which activated at 95°C for 1~5 minutes. ② High sensitivity and specificity. ③ Higher amplification efficiency to templates contained 40% to 80% GC.	① qPCR. ② High sensitivity virus detection. ③ Adapted to templates with high GC content.
Neoscript RT Premix (SYBR qRT-PCR) ② ③ ④ ⑤	M513S	① Reverse transcription temperature range is between 42°C to 55°C. ② Excellent specificity and sensitivity. ③ High applicability in various high sensitivity detection.	① qRT-PCR. ② High sensitivity detection of human- and animal- derived RNA viruses.

PCR Premix

Product Name	Code	Characteristics	Recommended Application
2×Robustart Multi Premix ④ ⑥	PM201	① Applied to DNA multiplex amplification in End-point PCR. ② Applied to amplification of templates contained 40% to 80% GC.	① Multiplex PCR. ② DNA Microarrays. ③ Reverse Dot Blot.
2×Novel GC Rich Premix ④ ⑥	PM202	Applied to amplification of High GC templates, which contained 40% to 85% GC.	① Amplification of templates with high GC content. ② DNA Microarrays. ③ Reverse Dot Blot.

Product Name	Code	Characteristics	Recommended Application
2× Blood Direct Premix ④	MD205	① Successful tolerance to blood-derived inhibitors, available to extraction-free End-point PCR. ② Good amplification capability of blood with different anticoagulants, even the addition volume up to 45% in reaction system. ③ Applied to amplification of templates contained 40% to 60% GC. When GC content is 60% to 80%, keep blood addition between 4% to 10% in reaction system and the length of amplified fragments no longer than 1.2kb.	① Extraction-free amplification of blood samples. ② DNA Microarrays. ③ Reverse Dot Blot.
2× Blood Direct Multi Premix ④ ⑥	MD206	① Successful tolerance to blood-derived inhibitors, available to End-point Multiplex PCR. ② Good amplification capability when directly add blood with different anticoagulants, the maximum dosage is up to 20% and keep addition no more than 8% in multiplex PCR. ③ Applied to amplification of templates contained 40% to 60% GC. When GC content is 60% to 80%, keep blood addition between 4% to 10% in reaction system and the length of amplified fragment no longer than 1.2kb.	① Extraction-free amplification of blood samples. ② DNA Microarrays. ③ Reverse Dot Blot.
Neoscript RTase 1ST Strand cDNA Synthesis Kit ③ ⑤ ⑥	PM05	cDNA first strand synthesis	① cDNA first chain synthesis. ② Molecular Hybridization. ③ PCR or qPCR.
5× Neoscript RT Premix (RT-PCR) ③ ④ ⑤ ⑥	PM504	① Single or multiplex amplification in End-point One-Step RT-PCR. ② Applied to amplification of templates contained 40% to 60% GC.	① Single or multiplex RT-PCR. ② DNA Microarrays. ③ Reverse Dot Blot.

02. STR Multiplex PCR Premix

Product Name	Code	Characteristics	Recommended Application
2× STR Premix Fast	STR201	① Suitable for STR multiplex PCR in medicolegal investigation. ② Compatible with fast amplification procedure, equilibrium of performance in multiplex amplification. It performs high sensitivity and specificity. ③ It exhibits extremely high resistance to inhibitors found in blood, and available to Direct PCR of extraction-free blood and blood-collected card.	Paternity test, medicolegal investigation and so on.

03. Other PCR-related products

Product Name	Code	Characteristics	Recommended Application
UNG (Uracil-DNA Glycosylase)	E01	Digest templates contained dUTP at 50°C and lose activity when heated to 95°C.	Anti-pollution amplification.
TS-UNG (Temperature Sensitive UNG)	E04	① Digest templates contained dUTP between 20°C to 37°C and lose activity when heated to 50°C. ② Suitable for anti-pollution reverse transcription.	Anti-pollution amplification.
T4 Gene 32 Protein	AS02	A single strand DNA-binding protein (ssDNA) derived from phage T4.	① Facilitate digestion of restriction endonuclease; ② Enhance DNA synthesis in vitro; ③ Improve output and efficiency of reverse transcription in RT-PCR.

Product Name	Code	Characteristics	Recommended Application
Rnase Inhibitor	AS05	①No inhibition for Taq, AMV, MMLV and Phage RNA Polymerases (SP6, T7, or T3). ②It can keep stable at a low concentration of DTT (less than 1mM) because of stronger antioxidant capacity.	①cDNA first chain synthesis. ②RT-PCR or RT-qPCR. ③In vitro transcription and translation. ④Separation and purification of mRNA. ⑤Related application needed integrated RNA.
dNTP (A/G/C/T/U)	As08 (As08-1/-2/-3/-4/-5)	/	PCR Reaction





PART 03

INSTRUMENTS & REAGENTS



CM-8 / CM-8 Plus Centrifugal mixer

BTE-32 BTE-8 Automatic nucleic acid extraction instrument

Sample pretreatment products

Nucleic acid extraction kit (magnetic bead) Release reagent

Sample preservation solution



Zhuhai BTE Biotechnology Co.Ltd., founded in March 2021, is an automatic life science instrument and equipment supplier. The company has an experienced R & D team, and has declared 7 related product patents since its establishment.

With enterprising and innovative, integrity and responsible business philosophy, BTE will layout molecular diagnosis, immunology, scientific instruments and other fields. We commit to providing customers with high-quality product solutions and professional service.

LYO-READY RAW
MATERIALS

RAW MATERIALS FOR
NUCLEIC ACID DETECTION

INSTRUMENTS
& REAGENTS

CM-8 / CM-8 Plus

Centrifugal mixer

Ingenious design

Exquisite appearance
drill-cut design [Patent No.:202130389474.0](#)
high quality

Graded locking
Smooth operation [Patent No.:202110698687.0](#)
Small size and portable

Integrated pallet design
Good uniformity between tubes
Stable baseline in PCR



Intelligent actions

Compound knob [Patent No.:202121409450.8](#)
Smart Halo
Easy-to-use with start/stop combo button

Low-noise operation
Polarization shutdown
Automatically operation and reliable

Inductive control
Open cover protection [Patent No.:202121406810.9](#)
Smart and versatile

A new option to improve PCR hands-on efficiency and uniformity between tubes

Technical Specifications

CM-8 Technical Specifications	
Specifications	Parameter
Centrifugation RPM	1500-4000rpm
Minimum Speed Regulation	100rpm
Operating mode	Centrifugal and mix
Mixing Intensity Level	Strong, medium and soft
Rotor Capacity	8-strip tubes
Lock Mode of Rotor	Quick-release knob
Lock Force	Multilevel adjustment
Procedural memory	Yes
Program Cycle	1-999
Open cover protection	Induction Type
Abnormal vibration Halting	Yes
Input	AC 100~240V, 1.4A, 50/60HZ
Maximum input power	60W
Size(mm)	218mm×234mm×136mm
Weight	≤1.5kg

CM-8 Plus Technical Specifications	
Specifications	Parameter
Centrifugation RPM	1500-4000rpm
Minimum Speed Regulation	100rpm
Operating mode	Centrifugal, mix, Centrifugal and mix
Mixing Intensity Level	Strong, medium and soft
Rotor Capacity	8-strip tubes, 1.5/2mL Centrifuge tube
Lock Mode of Rotor	Quick-release knob
Lock Force	Multilevel adjustment
Procedural memory	Yes
Program Cycle	1-999
Open cover protection	Induction Type
Abnormal vibration Halting	Yes
Input	AC 100~240V, 1.4A, 50/60HZ
Maximum input power	60W
Size(mm)	218mm×234mm×136mm
Weight	≤1.5kg

BTE-8 / BTE-32

Full automatic nucleic acid extractor

NEW

BTE-8

High-efficiency operation

Operating 1-8 samples flexibly
Rotary mix, efficient and mute

Convenient operation

The consumables compartment
ejects automatically



Hardcore system

High precision & high efficiency
temperature control
50μl-3200μl huge operating volume

Intelligent design

Operation by one smart knob
Novel vertical design to save space

Technical Specifications

BTE-32 Technical Specifications	
Specifications	Parameter
Sample throughput	1~32
Process volume	20~1000μL
Bead recovery rate	≥98%
Inter-pore uniformity	CV<3%
Rate of recovery	70%-130%
Extraction capacity	Class fg
Compatibility	Suitable for a variety of magnetic bead method nucleic acid extraction and purification reagents
Program management	Create, edit, save, delete, and more
Pollution control	Exhaust ventilate, ultraviolet disinfection function, self-cleaning
Consumables	Standard 96 deep-well plates and 8-link magnetic rod sleeves
Maximum output power	350W
Data storage	Built-in SD card; 255 groups of programs can be stored
Size(mm)	360×382×480
Weight	≤20kg

BTE-8 Technical Specifications	
Specifications	Parameter
Sample throughput	1~8
Process volume	50μl~1000ul(standard), 1000~3200μl
Product name	Automatic nucleic acid extraction instrument
Code	BTE-8
Inter-pore uniformity	CV<3%
Rotate speed	Multi-position
Extracting time	8-20mins/time (typical value)
Illuminating system	Internal LED
Sterilization	Ultraviolet ray
Lysis/elution temperature	Room temperature~120°C
Maximum output power	90W
Power input	DC 24V
Size(mm)	170×360×356
Weight	<7kg

BTE-32

Intelligent recyclable magnetic adsorption, fast and efficient operation

The instrument can complete the nucleic acid extraction of 1-32 samples within 12-20 minutes, and the magnetic suction can be repeatedly cycled, and the recovery rate of magnetic beads is ≥98%.

Simple structure, unique appearance and powerful

10-inch color low-position touch screen, ergonomic design, silent electronic control cabin design, and card-type sliding interface, improving aesthetics and operating comfort.



Using BTE nucleic acid extraction kit (Magnetic Beads-based Method), the procedure of nucleic acid extraction can be completed within 8 minutes!

Convenient operation, and one-click disassembly

Simplified operation steps, user-friendly function settings, removable magnetic rods for easy maintenance and replacement, effectively extend the service life of the instrument.

Scientific and technological wisdom anti-fouling, stable and reliable

Automatic lighting control, intelligent matching of "disinfection" and "extraction" status, and intelligent negative pressure system, reducing cross-contamination risk.

Sample pretreatment products



LYO-READY RAW MATERIALS

RAW MATERIALS FOR NUCLEIC ACID DETECTION

INSTRUMENTS & REAGENTS

Nucleic acid extraction reagent (magnetic bead)

PRODUCT DETAILS			
Name	Features	Specification	Code
MagBio Viral DNA/RNA Kit I	①High sensitivity: DNA virus nucleic acid extraction concentration can be as low as 10IU/mL; RNA virus nucleic acid extraction concentration can be as low as 50IU/mL. ②Wide application: suitable for the extraction of various DNA and RNA; it can easily extract oral/pharyngeal swab washes, whole blood, plasma,serum, saliva, bronchoalveolar lavage fluid and other liquid samples nucleic acid. ③Simple and safe: avoid centrifugation, easy to operate; no toxic reagents such as chloroform and phenol, and high safety.	IME-DR0601-I	
		Bottled package (32/48/64/96T)	IME-DR0601-I-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0601-I-W32T/48T/64T/96T
MagBio Fast Viral DNA/RNA Kit	①Extensive range of samples: It can easily extract nucleic acids from liquid samples such as oral and nose swabs, whole blood, plasma, serum, saliva, and bronchoalveolar lavage fluid. ②Efficient and fast: It only takes 10 minutes for automated equipment to extract 96 sample programs at a time.	IME-DR0601-III	
		Bottled package (32/48/64/96T)	IME-DR0601-III-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0601-III-W32T/48T/64T/96T
MagBio Swab DNA/RNA Kit II	①Efficient and time-saving: The automated equipment takes only 8.5 minutes to extract swab samples. ②High sensitivity: RNA virus extraction concentration can be as low as 100copies/mL. ③ Simple operation: one-step cleavage and binding, one-step cleaning, one-step elution, and three-step method to complete nucleic acid purification. ④Stable and reliable: The reagents are equipped with automatic extraction equipment, and the test results are more stable and reliable.	IME-DR0101-II	
		Bottled package (32/48/64/96T)	IME-DR0101-II-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0101-II-W32T/48T/64T/96T

Name	Features	Specification	Code
MagBio Swab DNA/RNA Kit I	①Proteinase K-Free: No need to add proteinase K, the operation is easier. ②High sensitivity: RNA virus nucleic acid extraction concentration can be as low as 200copies/mL. ③ Efficient and time-saving: the automated equipment takes only 10.5 minutes to extract swab samples; ④Stable and reliable: The reagents are equipped with automatic extraction equipment, and the test results are more stable and reliable. ⑤High safety: It does not contain toxic reagents such as chloroform and phenol, and the operation is safer. ⑥ Simple operation: one-step cleavage and binding, one-step cleaning, one-step elution, and three-step method to complete nucleic acid purification.	IME-DR0101-I	
		Bottled package (32/48/64/96T)	IME-DR0101-I-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0101-I-W32T/48T/64T/96T
MagBio Swab DNA/RNA Kit III	①Efficient and time-saving: The automated equipment takes only 10.5 minutes to extract swab samples. ②High sensitivity: RNA virus extraction concentration can be as low as 200copies/mL. ③ Simple operation: one-step cleavage and binding, one-step cleaning, one-step elution, and three-step method to complete nucleic acid purification. ④Full-component alcohol-free: The whole reagent is alcohol-free, free from drying and safer.	IME-DR0101-III	
		Bottled package (32/48/64/96T)	IME-DR0101-III-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0101-III-W32T/48T/64T/96T
MagBio Plasma/Serum DNA/RNA Kit	①Efficient and time-saving: The automated equipment takes only 8.5 minutes to extract plasma samples. ②High sensitivity: DNA virus extraction concentration can be as low as 50IU/mL. ③ Simple operation: one-step cleavage and binding, one-step cleaning, one-step elution, and three-step method to complete nucleic acid purification. ④ Alcohol-free cleaning solution: Cleaning is alcohol-free, free from drying, and faster.	IME-DR1001-I	
		Bottled package (32/48/64/96T)	IME-DR1001-I-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0101-III-W32T/48T/64T/96T
MagBio Blood DNA/RNA Kit II	①Convenient and safe: avoid centrifugation, easy to operate, and high safety; no toxic reagents such as chloroform and phenol. ②High sensitivity: The DNA virus extraction concentration in the whole blood sample can be as low as 50IU/mL. ③ Alcohol-free cleaning solution: Alcohol-free cleaning, free from drying, and more stable.	IME-DR1001-II	
		Bottled package (32/48/64/96T)	IME-DR1001-II-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR1001-II-W32T/48T/64T/96T

Name	Features	Specification	Code
MagBio Viral DNA/RNA Kit II	①Universal application: suitable for the extraction of various DNA and RNA; it can easily extract nucleic acid from oral/pharyngeal swab washing fluid, whole blood, plasma,serum, saliva, bronchoalveolar lavage fluid and other liquid samples . ②High sensitivity: RNA virus extraction concentration can be as low as 200copies/mL. ③ Enzyme-free treatment: no need to add proteinase K. ④Stable and reliable: The test results are more stable and reliable when using with automatic extraction equipment .	IME-DR0601-II	
		Bottled package (32/48/64/96T)	IME-DR0601-II-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0601-II-W32T/48T/64T/96T
MagBio Tissue DNA/RNA Kit	①It can process various types of animal tissue samples such as genitals, viscera, and skin. ②Flexible sample pre-processing methods can be provided. ③Easy to use and can be adapted to automated instruments for nucleic acid extraction.	IME-DR0401	
		Bottled package (32/48/64/96T)	IME-DR0401-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0401-W32T/48T/64T/96T
MagBio Plant DNA/RNA Kit	①It can process various types of plant samples such as leaves, seeds, rhizomes, etc. ②It is easy to use and can be adapted to automated instruments for nucleic acid extraction.	IME-DR0501	
		Bottled package (32/48/64/96T)	IME-DR0501-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-DR0501-W32T/48T/64T/96T
MagBio Starch DNA Kit	①Suitable for DNA extraction from sweet potato, cassava, corn and other edible starches. ②It is easy to use and can be adapted to automated instruments for nucleic acid extraction. ③ Safe and non-toxic.	IME-D1201	
		Bottled package (32/48/64/96T)	IME-D1201-B32T/48T/64T/96T
		Deep-well Multiwell Plate package (32/48/64/96T)	IME-D1201-W32T/48T/64T/96T
MagBio Free DNA Kit	①It is suitable for extracting free nucleic acid from plasma/serum samples. ②The sample processing volume is flexible and optional, maximum 10ml. ③Easy to use, safe and non-toxic.	IME-D1701	
		Bottled package (50/100)	IME-D1701-B50T/100T
		Deep-well Multiwell Plate package (50/100)	IME-D1701-W50T/100T
		* Welcome to consult the detailed sample size processing range	

Disposable virus sampling tube

NEW

Disposable virus sampling tube products are mainly used for sample collection, transportation and storage. By fully mixing the collected samples with the virus nucleic acid preservation solution, the virus in the samples can be inactivated, and the virus in the samples can be effectively guaranteed. The integrity of the viral nucleic acid is maintained, the collected samples can be transported at room temperature and stored for a long time, and the preserved viral nucleic acid samples can be widely used in PCR detection, etc.

Product Features

Safe and stable

Inactivated samples to protect operator safety; sample nucleic acids can be stored stably for 14 days at room temperature

High compatibility

Compatible with a variety of nucleic acid extraction reagents.

Simple operation

Easy to use and operate

Wide applicability

Suitable for mouth, nose and throat swab samples, serum/plasma samples, virus culture fluid, saliva and other samples.

Product model number / Specification

Product composition/ Model number	SP-1	SP-2	SP-3	AP-1	AP-2
Code	CYG01-SP-1	CYG01-SP-2	CYG01-SP-3	CYG01-AP-1	CYG01-AP-2
Viral nucleic acid preservation solution	2mL/Tube×50	6mL/Tube×50	11mL/Tube×25	3mL/Tube×50	5mL/Tube×50
sampling tube	10mL Tube×50	10mL Tube×50	30mL Tube×25	10mL Tube×50	10mL Tube×50
Sampling swab	1 swab×50	1 swab×500	1 swab×500	1 swab×50	1 swab×50
Product instruction	1	1	1	1	1
Remark	One human sample	10 in 1	20 in 1	Animal sample	Environmental sample
Special Note	Sampling swabs are optional and packaged separately, which is not included in the kit.				





qRT-PCR&RT-PCR raw materials

- MMLV reverse transcriptase
 - Tth
 - One-step premix
- can match with UNG+dUTP

qPCR&PCR raw materials

- Hot start Taq polymerase
 - Premix
- can match with UNG+dUTP

PCR ancillary product

- UNG
- TS-UNG
- dNTP

Isothermal amplification

hot start version lyophilized

- LAMP
- RT-LAMP
- Multiple enzyme fast isothermal amplification

Biomedicine

- mRNA vaccine raw materials enzyme
- Quality and safety control of biological products

Mycoplasma CHO E.coli Vero Pichia
Automatic extraction instrument and reagent

Molecular diagnostic POCT raw materials

lyophilized fast resistance to inhibitor
can match with UNG+dUTP

- Lyophilized enzyme
- Premix
- Direct, fast PCR raw materials

(magnetic bead) Nucleic acid extraction kit

- Universal blood, swab
- Quick swab
- Alcohol-free blood, swab


Instrument

CM-8

All-in-one Centrifugal Mixer

BTE-8/BTE-32

Automatic nucleic acid extraction instrument



Enterprise duty

Develop upstream technology and service diagnosis industry.
Create value for customers, teams and society.



Enterprise vision

Pursuit of excellence,
Guardian Life



Enterprise values

Honest and responsible,
enterprising and innovative

CORPORATE CULTURE



Official website



LinkedIn

Zhuhai Biori Biotechnology Co.,Ltd

Tel: +86-0756-8699969

Email: marketing@biori.com.cn

Web: www.biori-en.com

Address: No. 333, Pingbei 1st Road, Nanping Science Park, Xiangzhou District, Zhuhai