Ultrasonic Skin Scrubber

FEATURES

#Perfect cleaning effect
High frequency vibration, Minimize the irritation to the skin and helps the deep cleansing inside the pores.
#Convenient LCD screen
It will be very conveniently & clearly to know the operation mode & method of the scrubber by the LCD screen, even without reading the manual.
#Effective EMS function
The EMS function can wake up the sleeping collagen and help pore shrinking at the same time.
#Nutrients penetrate the skin
By generating galvanic ion, it helps the skin absorb nutrients better.

01 This manual including the composition, usage, precautions of the product.
02 Please fully using this instruction manual and precautions before using this product and familiarize yourself with it.
03 The product's waterproof grade is IP X-5
04 Please charging before the first use of the product.
05 This product is not a medical device, but a beauty device.
06 This product is a beauty device product that has been tested by the Korea Institute of Dermatological Sciences for human use.
07 The effect of this product varies depending on the situation of each person and the use of skincare product.
08 Be careful to keep and use.This product is not a baby product. Please keep it away from the children .
09 It is not an edible product. Please do not eat or bite.
10 Please put the product into the product box during storage. Do not expose in sunlight directly, keep it in a cool place, free of moisture.

Page 1

CONTENTS

1.Product Information
1-1 Components ----- Page 3
1-2 Name of each part ----- Page 3
1-3 Product Specifications ----- Page 4
2.Using Method
2-1 Charging method ----- Page 5
2-2 How to use ----- Page 6
3.Precautions
3-1. Notes ----- Page 8

Page 2

1.Product Information

1-1 Components
This product includes the following components

Product Box

One instruction manual

USB magnetic charging cable

An Ultrasonic Scrubber Device

1-2 Name of each part

Outproof cap

Double curve shovel

LCD Screen

Electric pole

ON/OFF switch

USB magnetic charging cable

Charging connector

Page 3

1.Product Information

1-3. Product Specifications
The specifications of this product are as follows:

Name:	Ultrasonic Skin Scrubber
Battery:	330mAh@3.7V
Size:	24*12*118mm (Width X Length X Height)
Weight:	120g
Material:	ABS + 304 stainless steel
Year of manufacture:	
Country of manufacture:	People's Republic of China.
Main function:	Skin Cleaning, Lifting, Help nutrient absorb.

Page 4

2. Using Method

2-1 Charging Method
Please charging before the first use of the product:
(The battery charging status is different, the charging time may vary.Maximum 2 hours)

① Connect the micro USB magnetic cable with the charging connector on the product.
② Connect the other side of the cable to a 5V 1A plug, or connect it to your computer for charging.
③ The battery pattern on LCD screen is flashing while charging, but if you want to use the product during charging, please disconnect the USB cable.

Please refer to!

If the battery icon is blinking during use, it means the battery power is less, so you need to charge it.
I The product cannot be used during charging.
II The battery pattern on LCD screen is flashing while charging.
III After the battery is fully charged, the maximum usage time is about 1 hour. Depending on the usage mode, the battery status, charging status and environment will be different.

Page 5

2. Using Method

2-2. How to use

Please familiarize yourself with the function of each mode before using it.

This product has the function of automatically shutting down after 5 minutes working.
If you need to use it extra, please press the power button again and select the mode you need.
When using the EMS mode, it is recommended to start with EMS/L and change to EMS/H or keep it depending on the feeling of use.

1.Sonic Cleansing
Function Intro
High frequency vibration to remove the dead skin

How To Use
1. Please wash your face first and ensure that your face is fully wetted.
2. Turn on the switch and choose the Sonic cleansing mode.Move the shovel vertically and gently on the skin.

2.Sonic Cleansing(+ION)
Function Intro
The positive ions will help to remove cosmetics residual from the skin.

How To Use
1. Please wash your face first and ensure that your face is fully wetted.
2. Turn on the switch and choose the Sonic cleansing (+ION) mode. Move the shovel vertically and gently on the skin.

3.Galvanic (-ION)
Function Intro
The skincare products themselves have negative ions, we use the characteristic of mutual exclusion of negative ions to help the skincare product be absorbed deeper.

How To Use
1. Apply skincare product like cream on the face.
2. Choose the Galvanic mode and hold the electric pole position on the product, then use the back side of the shovel and push it gently on the skin.

4.EMS/L(EMS LOW)
Function Intro
The electrical stimulation can shrink the skin immediately, the shrink of the muscle can help the tightening of the skin's face(cheek).

How To Use
1. Apply skincare product like cream on the face.
2. Choose the EMS/L mode and hold the electric pole position on the product, then use the back side of the shovel and push it gently on the skin.

5.EMS/H(EMS HIGH)
Function Intro
The electrical stimulation can shrink the skin immediately, the shrink of the muscle can help the tightening of the skin's face(cheek).

How To Use
1. Apply skincare product like cream on the face.
2. Choose the EMS/H mode and hold the electric pole position on the product, then use the back side of the shovel and push it gently on the skin.

Page 6

3. Precautions

3-1 Notes

1. Before using and storing, you must fully understand the precautions using and store product according to its intended use.
2. The product's waterproof grade is IP X-5. Please be careful not to let water or liquid enter inside of product .
3. This product is not a medical device.
4. This product is a beauty device product that has been tested by the Korea Institute of Dermatological Sciences for human use.
5. The effect of this product varies depending on the situation of each person and the use of skincare products.
6. If the skin is abnormal due to using of product, stop using it immediately and consult a professional doctor.
7. If there is any problem with the product, please stop using it immediately and consult the manufacturer or seller or the place of purchase immediately.
8. Don't throw or hit products.
9. Do not store in a place where the temperature is too low or too high and in a damp place.
10. Do not put the product in quilts, electric blankets and cars.
14. Please pay attention do not let other metal products or liquids come into or contact with the charger of the product.
15. The user's behavior of disassembling, assembling, repairing and modifying the product may cause of the malfunction, and the resulting product failure cannot be repaired.

Page 7

3. Precautions

3-1 Notes

1. Before using and storing, you must fully understand the precautions using and store product according to its intended use.
2. The product's waterproof grade is IP X-5. Please be careful not to let water or liquid enter inside of product .
3. This product is not a medical device.
4. This product is a beauty device product that has been tested by the Korea Institute of Dermatological Sciences for human use.
5. The effect of this product varies depending on the situation of each person and the use of skincare products.
6. If the skin is abnormal due to using of product, stop using it immediately and consult a professional doctor.
7. If there is any problem with the product, please stop using it immediately and consult the manufacturer or seller or the place of purchase immediately.
8. Don't throw or hit products.
9. Do not store in a place where the temperature is too low or too high and in a damp place.
10. Do not put the product in quilts, electric blankets and cars.
12. This product is not a baby product. Be careful to keep and use. Please keep it away from the children .
13. This product is not an edible product. Do not bite.
14. Please pay attention do not let other metal products or liquids come into or contact with the charger of the product.
15. The user's behavior of disassembling, assembling, repairing and modifying the product may cause of the malfunction, and the resulting product failure cannot be repaired.

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