

地址/Add:中国浙江省杭州市临安区玲珑工业园锦溪南路1288号 No.1288 South Jinxi Road, Linglong Industry Park, Lin'an District, Hangzhou City, Zhejiang Province, 311301, China. 电话/Tel:+86-571-23652388 传真/Fax:+86-571-23651000 网址/Web:http://www.fushanhz.com http://kingmedical.en.alibaba.com

邮箱/Email:export@fushanhz.com

杭州富善医疗器械有限公司

Hangzhou Fushan Medical Appliances Co., Ltd.



杭州富善医疗器械有限公司创立于2005年7月,坐落在杭州市西郊风景优美的临安区玲珑工业园内。

公司占地面积13000平方米,按照GMP要求设计建造了多个万级和十万级净化车间。车间配备 了多套精密的生产和检验设备,可以生产硅橡胶,聚氨酯,聚氯乙烯等医用高分子材料制作的医用 耗材。

同时,公司设计建造了专业的实验室,方便进行无菌、理化和电气试验等。公司配备了符合 CE标准的环氧乙烷灭菌柜,可以在公司内实现产品的灭菌操作并进行灭菌确认,为客户提供精确 的灭菌验证报告。

公司生产销售产品主要包括:呼吸麻醉类、肠胃营养类、创伤引流类、泌尿及妇科类相关耗材。所有产品都已获得欧盟CE认证,部分产品在美国FDA备案。公司还为广大客户提供专业的 OEM和ODM服务。

"规范、严谨、专注、创新"一直是公司发展的指导方针。富善人始终致力于关怀人类健康生活,提供优质医疗用品。

Hangzhou Fushan Medical Appliances Co., Ltd. was established in July 2005. It locates at Lin'an Linglong Industrial Park in the west of Hangzhou.

The factory of Fushan Medical covers an area of 13,000 sqms. Both class 10,000 and class 100,000 clean workshops are designed and constructed according to GMP standards. With various precise manufacture and inspection facilities, Fushan provides different kinds of medical devices by silicone, polyurethane (PU) and polyvinyl chloride (PVC).

Meanwhile, Fushan designed and built professional labs for sterile, physicochemical and electrical tests. It also equipped with CE verified ETO sterilization chamber. Sterilization verification report can be issued as request.

Fushan manufactures and sells various kinds of products, which refers to anesthesia series, gastroenterology series, surgical drainage series, urology and gynecology series etc. All products are granted with CE certificate, and most of them are registered in FDA. Fushan also provides OEM/ODM service to the market.

With the philosophy of "standard, preciseness, concentration and innovation", Fushan is always keeping providing high grade medical devices to the world.



麻醉呼吸类 ANESTHESIA SERIES

麻醉呼吸类



一次性硅胶喉罩 Disposable Silicone Laryngeal Mask Airway

100%采用医用硅胶生产,具有良好的生物相容 性,无毒无刺激。

Made of 100% medical grade silicone, with good bio-compatibility, non-toxic.

- 重复性硅胶喉罩 Reusable Silicone Laryngeal Mask Airway
- 经高压蒸汽灭菌可重复使用40次,每个喉罩都配有一个独特的序列号和对应的灭菌记录卡。 Suitable for autoclave sterilization. Reusable up to 40 times. With trackable serial number.

● 提供带格栅和无格栅罩囊供客户选择;提供透明罩囊和肉色罩囊供客户选择。

Cuff bar: with/without bar for choice.

Cuff color: transparent/pink for choice.

产品编号	Product No.	规格	罩囊最大充气量(ML)	放气板颜色	建议适用患者体重
Disposable	Reusable	Size	Maximum Cuff Inflation Volumes(ML)	Deflating Cap Color	Suggested Patient Weight (KG)
LM1004S1	LM1004S2	1	4	粉色/Pink	0-5
LM1507S1	LM1507S2	1.5	7	蓝色/Blue	5-10
LM2010S1	LM2010S2	2	10	绿色/Green	10-20
LM2514S1	LM2514S2	2.5	14	橙色/Orange	20-30
LM3020S1	LM3020S2	3	20	红色/Red	30-50
LM4030S1	LM4030S2	4	30	黄色/Yellow	50-70
LM5040S1	LM5040S2	5	40	紫色/Purple	70-100

一次性加强型硅胶喉罩

Disposable Reinforced Silicone Laryngeal Mask Airway

◆ 采用进口的医用级硅胶生产, 气管经过钢丝加强, 可 以承受头部和颈部的压力,防止气道扭折,压扁或堵 塞。

Made of imported medical grade silicone, with stainless steel spring. The airway is kink-resistant with reinforced tube.



重复性加强型硅胶喉罩

Reusable Reinforced Silicone Laryngeal Mask Airway

- ◆ 经高压蒸汽灭菌可重复使用40次,每个喉罩都配有一个独特的序列号和对应的灭菌记录卡。 Suitable for autoclave sterilization. Reusable up to 40 times. With trackable serial number.
- 是供带格栅和无格栅罩囊供客户选择:提供透明罩囊和肉色罩囊供客户选择。 Cuff bar: with/without bar for choice. Cuff color: transparent/pink for choice.

产品编号	Product No.	Product No. 规格 罩囊最大3		放气板颜色	建议适用患者体重
Disposable	Reusable	Size	Maximum Cuff Inflation Volumes(ML)	Deflating Cap Color	Suggested Patient Weight (KG)
LM1004S1R	LM1004S2R	1	4	粉色/Pink	0-5
LM1507S1R	LM1507S2R	1.5	7	蓝色/Blue	5-10
LM2010S1R	LM2010S2R	2	10	绿色/Green	10-20
LM2514S1R	LM2514S2R	2.5	14	橙色/Orange	20-30
LM3020S1R	LM3020S2R	3	20	红色/Red	30-50
LM4030S1R	LM4030S2R	4	30	黄色/Yellow	50-70
LM5040S1R	LM5040S2R	5	40	紫色/Purple	70-100

麻醉呼吸类 ANESTHESIA SERIES

一次性PVC喉罩 Disposable PVC Laryngeal Mask Airway

♥ 采用医用级PVC生产。 Made of medical grade PVC.

● 罩囊外形符合人体咽喉部解剖结构,最大限度 减少对患者咽喉部的刺激,气密性好。

The cuff is designed according to anatomical structure of human pharynx and larynx, which reduces stimulation to patient and leakage risk.



一次性PVC硅胶喉罩 Disposable PVC-Silicone Combined Laryngeal Mask Airway

◆ 喉罩气管采用医用级PVC生产,罩囊采用医用 级硅胶生产。

Its tube is made of medical grade PVC, and its cuff is made of medical grade silicone.

◆ 结合硅胶喉罩和PVC喉罩的特点,兼顾舒适性 及经济性。

It combines the advantage of PVC LMA and silicone LMA, with good comfortability and economy.

产品编号 Product No.	<mark>规格</mark> Size	罩囊最大充气量(ML) Maximum Cuff Inflation Volumes(ML)	放气板颜色 Deflating Cap Color	建议适用患者体重 Suggested Patient Weight (KG)
LM1004P	1	4	粉色/Pink	0-5
LM1507P	1.5	7	蓝色/Blue	5-10
LM2010P	2	10	绿色/Green	10-20
LM2514P	2.5	14	橙色/Orange	20-30
LM3020P	3	20	红色/Red	30-50
LM4030P	4	30	黄色/Yellow	50-70
LM5040P	5	40	紫色/Purple	70-100

产品编号 Product No.	规格 Size	罩囊最大充气量(ML) Maximum Cuff Inflation Volumes(ML)	放气板颜色 Deflating Cap Color	建议适用患者体重 Suggested Patient Weight (KG)
LM1004P1	1	4	粉色/Pink	0-5
LM1507P1	1.5	7	蓝色/Blue	5-10
LM2010P1	2	10	绿色/Green	10-20
LM2514P1	2.5	14	橙色/Orange	20-30
LM3020P1	3	20	红色/Red	30-50
LM4030P1	4	30	黄色/Yellow	50-70
LM5040P1	5	40	紫色/Purple	70-100

麻醉呼吸类

一次性加强型硅胶气管插管 Disposable Reinforced Silicone Endotracheal Tube

- ◆ 100%采用医用级硅胶生产;管体柔软透明。 Made of 100% medical grade silicone, the tube is soft and clear.
- 📀 带高容量低压力气囊。 With high volume and low pressure silicone cuff.
- 通气管经过弹簧加强,可防止管体压扁、扭曲。 The main tube is reinforced with stainless steel spring to
- prevent crushing and kinking. ◇ 带墨菲孔。

With Murphy eye design. 立口炉豆 切枝 立口炉豆



PVC气管插管 PVC Endotracheal Tube

- ◎ 采用医用级PVC生产,管体透明度高。 Made of medical grade PVC, the tube is transparent and clear.
- ♦ X光造影线贯穿管身。
- X-ray opaque line through the whole tube. ● 高容量低气压气囊可供选择。
- High volume and low pressure cuff is available for choice.
- ◆ 带墨菲孔。 With Murphy eye design.

Product No.	Size	Froduct No.	Size	Product No.	Size
ET40SR	4.0	ET65SR	6.5	ET90SR	9.0
ET45SR	4.5	ET70SR	7.0	ET95SR	9.5
ET50SR	5.0	ET75SR	7.5	ET100SR	10.0
ET55SR	5.5	ET80SR	8.0		
ET60SR	6.0	ET85SR	8.5		

重复性加强型硅胶气管插管 Reusable Reinforced Silicone Endotracheal Tube

通气管经 The main	低压力 volume 过弹簧 tube is	气囊。 and low pres 加强,可防1	sure silic 上管体压			
带墨菲孔 With Murr				更用达10次。 ple up to 10 tir	nes.	
产品编号 Product No.	规格 Size	产品编号 Product No.	规格 Size	产品编号 Product No.	规格 Size	
T40S2R	4.0	ET65S2R	6.5	ET90S2R	9.0	
	4.0 4.5	ET65S2R ET70S2R	6.5 7.0	ET90S2R ET95S2R	9.0 9.5	
T40S2R			0.0			
T40S2R	4.5	ET70S2R	7.0	ET95S2R	9.5	



产品编号	Product No.	规格	产品编号	Product No.	规格
Without Cuff	With Cuff	Size	Without Cuff	With Cuff	Size
ET20P	1	2.0	ET65P	ET65PC	6.5
ET25P	ET25PC	2.5	ET70P	ET70PC	7.0
ET30P	ET30PC	3.0	ET75P	ET75PC	7.5
ET35P	ET35PC	3.5	ET80P	ET80PC	8.0
ET40P	ET40PC	4.0	ET85P	ET85PC	8.5
ET45P	ET45PC	4.5	ET90P	ET90PC	9.0
ET50P	ET50PC	5.0	ET95P	ET95PC	9.5
ET55P	ET55PC	5.5	ET100P	ET100PC	10.0
ET60P	ET60PC	6.0			

麻醉呼吸类

口腔预成型PVC气管插管 Oral Preformed PVC Endotracheal Tube

◆ 采用医用级PVC生产,管体透明度高。

Made of medical grade PVC, the tube is transparent and clear.

♦ X光造影线贯穿管身。

X-ray opaque line through the whole tube.

- ◆ 高容量低气压气囊可供选择。 High volume and low pressure cuff is available for choice.
- 带墨菲孔。
 With Murphy eye design.
- 口腔预成型,经口使用最佳。

 With oral preformed design, it is better to
 use for oral insertion.



	产品编号 Product No.		产品编号 Product No. 规格 产品编号 Product		Product No.	/201H	
able	Without Cuff	With Cuff	Size	Without Cuff	With Cuff	Size	
	OT20P	1	2.0	OT65P	OT65PC	6.5	
	OT25P	OT25PC	2.5	OT70P	OT70PC	7.0	
	OT30P	OT30PC	3.0	OT75P	OT75PC	7.5	
	OT35P	OT35PC	3.5	OT80P	OT80PC	8.0	
	OT40P	OT40PC	4.0	OT85P	OT85PC	8.5	
	OT45P	OT45PC	4.5	OT90P	OT90PC	9.0	
	OT50P	OT50PC	5.0	OT95P	OT95PC	9.5	
	OT55P	OT55PC	5.5	OT100P	OT100PC	10.0	
	OT60P	OT60PC	6.0				

鼻腔预成型PVC气管插管 Nasal Preformed PVC Endotracheal Tube

● 采用医用级PVC生产,管体透明度高。

•

0

Made of medical grade PVC, the tube is transparent and clear.

X光造影线贯穿管身。 X-ray opaque line through the whole tube.	产品编号	Product No.	规格	产品编号	Product No.	规格
高容量低气压气囊可供选择。	Without Cuff	With Cuff	Size	Without Cuff	With Cuff	Size
High volume and low pressure cuff is available	NT20P	1	2.0	NT65P	NT65PC	6.5
for choice.	NT25P	NT25PC	2.5	NT70P	NT70PC	7.0
· 带墨菲孔。 With Murphy eye design.	NT30P	NT30PC	3.0	NT75P	NT75PC	7.5
鼻腔预成型,经鼻使用最佳。	NT35P	NT35PC	3.5	NT80P	NT80PC	8.0
With nasal preformed design, it is better to	NT40P	NT40PC	4.0	NT85P	NT85PC	8.5
use for nasal insertion.	NT45P	NT45PC	4.5	NT90P	NT90PC	9.0
\frown	NT50P	NT50PC	5.0	NT95P	NT95PC	9.5
	NT55P	NT55PC	5.5	NT100P	NT100PC	10.0
	NT60P	NT60PC	6.0			

加强型PVC气管插管 Reinforced PVC Endotracheal Tube

- 采用医用级PVC生产,管体透明度高。
 Made of medical grade PVC, the tube is transparent and clear.
- 通气管经过弹簧加强,可防止管体压扁、扭曲。 The main tube is reinforced with stainless steel spring to prevent crushing and kinking.
- ◆ 高容量低气压气囊可供选择。 High volume and low pressure cuff is available for choice.
- 带墨菲孔。
 With Murphy eye design.



产品编号	item No.	规格	产品编号	Product No.	规格
Without Cuff	With Cuff	Size	Without Cuff	With Cuff	Size
ET20PR	1	2.0	ET65PR	ET65PCR	6.5
ET25PR	ET25PCR	2.5	ET70PR	ET70PCR	7.0
ET30PR	ET30PCR	3.0	ET75PR	ET75PCR	7.5
ET35PR	ET35PCR	3.5	ET80PR	ET80PCR	8.0
ET40PR	ET40PCR	4.0	ET85PR	ET85PCR	8.5
ET45PR	ET45PCR	4.5	ET90PR	ET90PCR	9.0
ET50PR	ET50PCR	5.0	ET95PR	ET95PCR	9.5
ET55PR	ET55PCR	5.5	ET100PR	ET100PCR	10.0
ET60PR	ET60PCR	6.0			

麻醉呼吸类

ANESTHESIA SERIES

氧气面罩 Oxygen Mask

采用医用级PVC生产,面罩柔软且透明度高。
 Made of medical grade PVC, the mask is soft and clear.

配备可调节鼻夹。供氧管长度可定制。
 Nose clip is adjustable. The length of oxygen tube can be customized.

产品编号 Product No.	型号 Type
102MM-E	XL: Adult Elongated
102MM	L: Adult Standard
222MM-E	M: Pediatric Elongated
222MM	S: Pediatric Standard



吸痰管 Suction Tube

◆ 采用医用级PVC生产,管体透明度高。 Made of medical grade PVC, the tube is transparent and clear.

多种类型接头可供选择。
 Different kinds of connectors are available.



麻醉呼吸类

ANESTHESIA SERIES

产品编号 Product No.	规格(Fr) Size(Fr)	色标号 Color Code	长度(mm) Length(mm)	备注 Remark			
SC05	5	灰色/Gray					
SC06	6	浅绿色/Light Green					
SC08	8	浅蓝色/Light Blue					
SC10	10	黑色/Black	450	Cap-cone type Plain type			
SC12	12	白色/White					
SC14	14	绿色/Green		Plain type			
SC16	16	橙色/Orange					
SC18	18	红色/Red					
SC20	20	黄色/Yellow					

雾化面罩 Nebulizer Mask

采用医用级PVC生产,面罩柔软且透明度高。
 Made of medical grade PVC, the mask is soft and clear.

配备6ml雾化罐。
 6ml nebulizer jar is available.

产品编号 Product No.	型号 Type
AK06AE	XL: Adult Elongated
AK06AS	L: Adult Standard
AK06PE	M: Pediatric Elongated
AK06PS	S: Pediatric Standard

吸引连接管 Suction Connecting Tube

- ◆ 采用医用级PVC生产,管体透明度高。 Made of medical grade PVC, the tube is transparent and clear.
- ◆ 管体长度可定制。常规长度有: 1.8米, 2米, 2.5米, 3米等。 The length of tube can be customized. 1.8m, 2.0m, 2.5m and 3.0m are regular lengths for choice.

产品编号 Product No.	尺寸(m) _{Size(m)}	产品编号 Product No.	尺寸(m) _{Size(m)}	备注 Remark
ST18T	3/16" 1.8m	ST18H	1/4" 1.8m	
ST20T	3/16" 2.0m	ST20H	1/4" 2.0m	Customized
ST25T	3/16" 2.5m	ST25H	1/4" 2.5m	Length
ST30T	3/16" 3.0m	ST30H	1/4" 3.0m	



麻醉呼吸类 ANESTHESIA SERIES



口咽通气道 Oropharyngeal Airway / Guedel Airway

采用无毒聚乙烯材料生产。

Made of non-toxic PE.

● 配备不同颜色的色标以区分型号。

Equipped with different color code to indentify different sizes.

产品编号 Product No.	规格 Size	色标号 Color Code	长度(mm) Length(mm)
OA40A	000	粉红/Pink	40
OA50A	00	蓝色/Blue	50
OA60A	0	黑色/Black	60
OA70A	1	白色/White	70
OA80A	2	绿色/Green	80
OA90A	3	黄色Yellow	90
OA100A	4	红色/Red	100
OA110A	5	浅蓝色/Light Blue	110
OA120A	6	橙色/Orange	120



无针接头 Needle Free Connector

● 采用进口PC材料和医用级硅胶生产,抗感染能力强。 Made of imported PC and medical grade silicone, with high anti-infectivity.

◎ 密封性好,普通款耐受150KPa压力不漏液,正压款耐受 300KPa压力不漏液。

The connector is with high sealability. Standard type connector will not leak under pressure less than 150KPa; positive pressure type connector will not leak under pressure less than 300KPa.

● 6%鲁尔锥度接头设计,方便与市面上的留置针,导管, 注射器,输液器等配套使用。

With 6% Luer conical design, it can be used with regular IV cannula, CVC, syringe, infusion set and other kits conveniently.

◆ 给药和采样无针化,全面保护医务人员,避免职业暴露。

With needle free design for both drug delivery and sampling, it prevent the medical staffs properly from infection damage.



鼻氧管 Oxygen Cannula

● 采用医用级PVC生产,管体柔软且透明度高。 Made of medical grade PVC, the tube is soft and clear.

◆ 有透明管体和绿色管体可供选择。 Transparent tube and green tube are both available.

产品编号 Product No.	型号 Type
NOC-L	成人/Adult
NOC-M	儿童/Pediatric
NOC-S	婴儿/Infant
NOC-XS	新生儿/Neonate



规格 Product No.	5302	5303
分类 Type	普通型 Standard Type	正压封管型 Positive Pressure Type
多路设计 Multi−lumens Design	单路、双路、三路 Single lumen Double lumens Triple lumens	单路、双路、三路 Single lumen Double lumens Triple lumens
延长管长度(cm) Length of Extension Tube	10-30	10-30
尺 寸 Measure	12*24mm	10*36mm
流 速 Flow Velocity	220-260ml/min	160-220ml/min

血管通路类 VASCULAR ACCESS SERIES

血管通路类 VASCULAR ACCESS SERIES

一次性使用中心静脉导管包 Disposable Central Venous Catheter Kit

采用符合人体留置医用级专用PU材质,具备优良的生物相容性和化学稳定性优越的弹性,人体温度下自动软化更利于保护体血管组织,免于受损。

Made of medical grade PU material which is specifically used for human body indwelling. It is with excellent biocompatibility and chemical stability, as well as superior elasticity. The material will soften itself automatically to protect the vascular tissue under body temperature.



曾 腔 Lumen Design	导管规格x长度(cm) Catheter Specification* Length	扩张器外径(mm) Dilator OD (mm)	导丝(外径x长度) Guide wire (CD*Length)	Y针配置 Y-Shaped Syringe Needle	监空针配置 RaulersonSpring-Wire Introduction Syringe	备注 Remarks
	14G×16	2.70	0.032×60cm	availat	ble	
	146×20	2.70	0.032×60cm	availat	le	
0	14G×30	2.70	0.032x70cm	availat	ble	
()	16Gx16	2.00	0.032×60cm	availat	ole	
()	16G×20	2.00	0.032×60cm	availat	ole	
9	16Gx30	2.00	0.032×70cm	availat	ole	
单腔/Single	18G×13	1.70	0.025x50cm	availat	ole	少儿 Pediatric
	18G×16	1.70	0.025×50cm	availat	ole	少儿 Pediatric
	18G×20	1.70	0.025x50cm	availat	ble	少儿 Pediatric
	4 F ×8	1.70	0.018x50cm			婴儿 Infant
	4Fx13	1.70	0.018x50cm			婴儿 Infant
0	5Fx13	2.00	0.021x50cm			少儿 Pediatric
(\mathcal{L})	5Fx16	2.00	0.021x50cm		~~	少儿 Pediatric
	7F×16	2.70	0.035×60cm	availat	ble	
双腔/Double	7Fx20	2.70	0.035×60cm	availat	le	
	7F×30	2.70	0.035x70cm	availat	ble	
	4.5Fx10	1.70	0.021x50cm			少儿 Pediatric
	4.5Fx13	1.70	0.021x50cm			少儿 Pediatric
	5Fx10	2.00	0.021x50cm			少儿 Pediatric
60	5Fx13	2.00	0.021×50cm			少儿 Pediatric
(2)	5.5Fx13	2.00	0.021x50cm		-	
S	5.5F×16	2.00	0.021x50cm		-	
三腔/Triple	7Fx16	2.70	0.035x60cm	availat	ole	
	7Fx20	2, 70	0.035x60cm	availat	ole	
	7Fx30	2.70	0.035x70cm	availat	ble	



经皮内窥镜胃造瘘管/套装

Percutaneous Endoscopic Gastrostomy (PEG) Tube / Kit

全新的设计

Brand new design

根据市场反馈,富善从配件,托盘,到吸塑盒都进行了重新设计。 新套装使用更加便捷。

According to the feedback from market, Fushan redesigned the PEG kit totally from components to package. The new kit is more convenient to handle.

♥ 富善提供两种不同配置的PEG套装:标准装和豪华装。

Both economic pack and luxury pack PEG kit are available for choice.



富善的承诺 Fushan's commitment

◆ 富善医疗致力于开发优质的产品,提供可靠的服务,满足市场的需求。

Fushan provides products with excellent quality and reliable service to meet the requirement of market.

优异的品质

Excellent Quality

● 富善始终坚持着执行最高标准的生产工艺。产品特点包括:

Fushan insists in manufacturing products with the highest standards. The Features of PEG kit include:

· Medical Grade Silicone Construction	·医用级硅胶制造
· External Retention Ring	·固定圈
· Universal Feeding Port Adapters	通用喂食端连接头
· Collapsible Internal Retention Bumper	可折叠内部缓冲垫
· Radiopaque Stripe and Bumper	·X光造影线和缓冲垫
· Tubing Clamp	- 导管夹
· Pull Placement Methods	• 牵引替换方法

产品编号 Product No.		规格	外径	长度	备注
Economic kit	Luxury kit	Size	O.D.(mm)	Length(mm)	Remark
PEG-16SE	PEG-16SL	16FR	5.3		
PEG-18SE	PEG-18SL	18FR	6.0		Length could be customized
PEG-20SE	PEG-20SL	20FR	6.7	600	
PEG-22SE	PEG-22SL	22FR	7.3		customized
PEG-24SE	PEG-24SL	24FR	8.0		

PEG套装 (豪华装)

PEG Kit (Luxury Pack)

Traction Removable PULL PEG Tube ×1 PEG导管 圈套器 5 Fr Retrieval Snare (1.8mm x 160cm) × 1 Feeding Adapter ×1 Y接头 Feeding Adapter with Luer lock ×1 带鲁尔Y接头 Retention Ring ×1 固定圈 25Ga Needle ×1 25G针头 19/21Ga Needle ×1 19/21G针头 Stainless Steel Scissors ×1 不锈钢手术剪 Stainless Steel Hemostat ×1 不锈钢止血钳 Tubing Clamp ×1 导管夹 Water Soluble Lubricant ×1 外科润滑剂 纱布 10*10 Gauze Sponges ×4 带孔纱布 5*5 Slit Gauze Sponges ×2 Placement Wire (0.36/0.6mm × 260cm) × 1 导丝 10ml Syringe ×1 注射器 Disposable Scalpel with #11 Blade ×1 长柄手术刀 Introducer Needle ×1 套管针 Fenestrated Drape ×1 洞巾 Directions for Use ×1 使用说明书

PEG套装(标准装)

PEG Kit (Economic Pack)

Traction Removable PULL PEG Tube × 1 Feeding Adapter × 1 Feeding Adapter with Luer lock × 1 Retention Ring × 1 19/21Ga Needle × 1 25Ga Needle × 1 Placement Wire (0.36/0.6mm × 260cm) × 1 10ml Syringe × 1 Disposable Scalpel with #11 Blade × 1 Introducer Needle × 1 Tubing Clamp × 1 10*10 Gauze Sponges × 4 Fenestrated Drape × 1 Directions for Use × 1



肠胃营养类

GASTROENTEROLOGY SERIES

肠胃营养类 GASTROENTEROLOGY SERIES

硅胶鼻饲管 Silicone Nasogastric Tube / Feeding Tube



● 100%采用医用级硅胶生产,管体柔软透明。 Made of 100% medical grade silicone, the tube is soft and clear. ◆ X光造影线贯穿管身。 X-ray opaque line through the whole tube.

● 管体刻度清晰。 Clear depth marks on the tube.

◆ 尺寸: 5Fr、6.5Fr、8Fr。 Size: 5Fr, 6.5Fr, 8Fr.

硅胶胃造瘘管 Silicone Gastrostomy Tube

◆ 100%采用医用级硅胶生产,管体柔软透明。 Made of 100% medical grade silicone, the tube is soft and clear.

♦ X光造影线贯穿管身。

X-ray opaque line through the whole tube.

内外粘接球囊,柔软且坚韧。

The balloon is glued to the main tube both inside and outside, it is elastic and flexile.

◆ 适用于胃造瘘术的患者。

It is used for gastrostomy operation.

产品编号 Product No.	规格(Fr) Size(Fr)	色标 Color Code	长度(mm) Length(mm)	备注 Remark
GT1205	12Fr/5ml	白色/White		
GT1405	14Fr/5ml	绿色/Green		
GT1620 GT1820 GT2020	16Fr/20ml	橙色/Orange		一次性使用
	18Fr/20ml	红色/Red	240	
	20Fr/20ml	黄色/Yellow		Disposable
GT2220	22Fr/20ml	紫色/Purple		
GT2420	24Fr/20ml	深蓝/Dark Blue		
	Product No. GT1205 GT1405 GT1620 GT1820 GT2020 GT2220	Product No. Size(Fr) GT1205 12Fr/5ml GT1405 14Fr/5ml GT1620 16Fr/20ml GT1820 18Fr/20ml GT2020 20Fr/20ml GT2200 22Fr/20ml	Product No. Size(Fr) Color Code GT1205 12Fr/5ml 白色/White GT1405 14Fr/5ml 绿色/Green GT1620 16Fr/20ml 橙色/Orange GT1820 18Fr/20ml 橙色/Red GT2020 20Fr/20ml 黄色/Yellow GT2220 22Fr/20ml 紫色/Purple	Product No. Size(Fr) Color Code Length(mm) GT1205 12Fr/5ml 白色/White GT1405 14Fr/5ml 绿色/Green GT1620 16Fr/20ml 橙色/Orange GT1820 18Fr/20ml 红色/Red GT2020 20Fr/20ml 黄色/Yellow GT2220 22Fr/20ml 紫色/Purple



硅胶胃管 Silicone Stomach Tube

- 🥥 100%采用医用级硅胶生产,管体柔软透明。 Made of 100% medical grade silicone, the tube is soft and clear.
- ♦ X光造影线贯穿管身。 X-ray opaque line through the whole tube.
- 管体刻度清晰。 Clear depth marks on the tube.

硅胶重头胃管

Silicone Ryle's Stomach Tube

◆ 导管底端带不锈钢珠设计,便于导管在指定位置进行 有效的留置。

Stainless steel balls are sealed in the tip end; it is easier to indwell the tube at a certain position.

ct No	ct No.	规格 (Fr)	长度 (mm) Length(mm)	备注				
Ryle	Ryle's	Size(Fr)	Standard	Ryle's	Remark				
т08	T08SR	8	700 700		700	700	700	700	
T10	T10SR	10	700	700	一次性使用 Disposabe				
T12	T12SR	12							
T14	T14SR	14	1200	1200					
T16	T16SR	16	1200	1200 1200					
T18	T18SR	18							



肠胃营养类 GASTROENTEROLOGY SERIES

PU胃管/鼻饲管 PU Stomach Tube / Nasogastric Tube

♥ 采用医用级PU生产,管体柔软透明,耐腐蚀性好。 Made of medical grade PU, the tube is soft and clear, with high corrosion resistance.

♦ X光造影线贯穿管身。

X-ray opaque line through the whole tube.

● 管体刻度清晰。 Clear depth marks on the tube.

● 配备不同颜色的接头以区分不同 型号的导管。

Nurse can easily indentify different size catheter with different color connectors.



PVC胃管 PVC Stomach Tube

采用医用级PVC生产,管体透明度高。 Made of medical grade PVC, the tube is transparent and clear.

♦ X光造影线贯穿管身。 X-ray opaque line through the whole tube.

< 配备不同颜色的接头以区分不同型号的导管。</p> Nurse can easily indentify different size catheter with different color connectors.



PVC重头胃管 PVC Ryle's Stomach Tube

● 导管底端带不锈钢珠设计,便于导管在指定位置进行有效的留置。 Stainless steel balls are sealed in the tip end; it is easier to indwell the tube at a certain position.

产品编号 Product No.				长度(mm) Length(mm)		刻度位置(cm)	
Standard	Ryle's	Size(Fr)	Color Code	Standard	Ryle's	Marked Position(cm)	
ST10P	ST10PR	10	黑色/Black				
ST12P	ST12PR	12	12 白色/White				
ST14P	ST14PR	14	绿色/Green	1200	1200 1250	45/55 65/75	
ST16P	ST16PR	16	橙色/Orange				
ST18P	ST18PR	18	红色/Red				
ST20P	ST20PR	20	黄色/Yellow				

胃管产品编号 Stomach Tube Product No.	规格(Fr) _{Size(Fr)}	长度(mm) Length(mm)	鼻饲管产品编号 Nasogastric Tube Product No.	规格(Fr) _{Size(Fr)}	长度(mm) Length(mm)
ST08T	8	900	NA40TS	4.0	600
ST10T	10	900	NA50TS	5.0	600
ST12T	12	1250	NA50TL	5.0	1400
ST14T	14	1250	NA65TS	6.5	600
ST16T	16	1250	NA65TL	6.5	1400
ST18T	18	1250	NA80TS	8.0	600
ST20T	20	1250	NA80TL	8.0	1400



PVC Nasogastric Tube

采用医用级PVC生产,管体透明度高。
 Made of medical grade PVC, the tube is transparent and clear.

♦ X光造影线可供选择。

X-ray opaque line is available for choice.

🔷 不同规格采用不同色标区分。

Different sizes tubes are identified with color-coded connector.



硅胶扁平开槽引流管 Silicone Flat Fluted Drain

● 100%采用医用级硅胶生产。

Made of 100% imported medical grade silicone.

扁平引流管由一段透明的硅胶圆管和一段白色扁管组成。内 壁有槽状结构,防止因挤压或纠结引起的管体阻塞。

This Drain combines a transparent silicone round tube and a white flat drain. The inner wall of the implanted flat segment features multiple narrow ducts to prevent possible occlusion when the tube is pressed or kinked.

◆ 全开槽适应最佳引流选择。

Full fluted design is available for optimal drainage with X-ray opaque effectiveness throughout its entire implanted segment.

可使用不锈钢引导针进行手术穿刺放置,并用负压球 收集引流液体,成套使用。

Using Trocar for insertion and combining with silicone Reservoir for fluids suction and collection.

产品编号 Product No.	扁平管宽度 Width of flat part (mm)	扁平管长度 Length of flat part (mm)				
PFD083008	8	300				
PFD112008	11	300				
Remark	Length could	be customized.				

硅胶圆形开槽引流管 Silicone Round Fluted Drain

- 100%采用医用级硅胶生产。
 Made of 100% imported medical grade silicone.
- ◆ 带有X光造影线,用于闭合创伤引流手术。 With X ray opaque line, the round drain is used for closed wound drainage surgical.
- ◆ 管体内部有四个腔道,有利于更快的引流分泌物。 Four channel design is better for a quicker drainage of blood or secretion.
- ◆ 管体采用一体挤出工艺,无多余接口。术后移除引流管时,可减少 患者痛苦。

The tube is extruded in one shot. With hubless design, it reduces patient's pain while removing the drain after surgical.

	品编号 luct No.	规格 Size(Fr)	外径 O.D.(mm)	长度 Length(mm)	备注 Remark
RFD	01030	10Fr	3.3	1200	
RFD	01530	15Fr	5.0	1200	Length could
RFE	01930	19Fr	6.3	1200	be customized
RFD	02430	24Fr	8.0	1200	

产品编号 Product No.		规格(Fr) 色标	长度(mm)	刻度位置(cm)	
Without X-ray Line	With X-ray Line	Size(Fr)	Color Code	Length(mm)	Marked Position(cm)
FT04PS	FT04PSX	4	紫色/Purple	400	14/15/16/20
FT05PS	FT05PSX	F	灰色/Gray	400	14/15/16/20
FT05PL	FT05PLX	5	5 AE/Gray	1000	45/55/65/75
FT06PS	FT06PSX	6	淡绿色	400	14/15/16/20
FT06PL	FT06PLX	D	Light Green	1000	45/55/65/75
FT08PS	FT08PSX		淡蓝色	400	14/15/16/20
FT08PL	FT08PLX	8	Light Blue	1000	45/55/65/75
FT10PL	FT10PLX	10	黑色/Black	1000	45/55/65/75

创伤引流类 WOUND DRAINAGE SERIES

硅胶圆形打孔引流管 Silicone Round Perforated Drain

● 100%采用医用级硅胶生产。

Made of 100% imported medical grade silicone.

♦ X光造影线贯穿管身。

X-ray opaque line through the whole tube.

● 与负压引流球及引导针配合使用。 It is used together with reservvir and trocar.

管体透明,便于观察引流情况。

Transparent tube design makes drainage situation visible.

产品编号 Product No.	规格(Fr) Size(Fr)	外径 O.D.(mm)	
RPD10	10	3.3	
RPD12	12	4.0	
RPD14	14	4.7	
RPD16	16	5.3	
RPD18	18	6.0	
RPD20	20	6.7	
Remark	Length could be customize		



PVC圆形打孔引流管 PVC Round Perforated Tube

- ◆ 采用医用级PVC生产,管体透明度高。 Made of medical grade PVC, the tube is transparent and clear.
- X造影线贯穿管身。 X-ray opaque line through the whole tube.
- 端部32个侧孔设计,引流效率更高。 With 32 lateral eyes design on distal end, the drainage efficiency is much higher.

● 型号: 10Fr-18Fr。 Size: 10Fr-18Fr.



硅胶扁平打孔引流管 Silicone Flat Perforated Drain

- ◆ 100%采用医用级硅胶生产。 Made of 100% imported medical grade silicone.
- 扁平引流管由一段透明的硅胶圆管和一段白色扁管组成。内壁有齿状结构,防止因挤压或纠结引 起的管体阻塞。

This Drain combines a transparent silicone round tube and a white flat drain. The inner wall of the implanted flat segment features tooth-like to prevent possible occlusion when the tube is pressed or kinked.

● 扁平管3/4打孔,或全尺寸打孔,适应最佳引流选择。

3/4 or full perforation are available for optimal drainage with X-ray opaque effectiveness throughout its

● 可使用不锈钢引导针进行手术穿刺放置,并用负压球收集引流液体,成套使用。 Using Trocar for insertion and combining with silicone Reservoir for fluids suction and collection.





硅胶胸腔引流管

Silicone Chest Drain

- 100%采用医用级硅胶生产。 Made of 100% imported medical grade silicone.
- ♦ X光造影线贯穿管身。 X-ray opaque line through the whole tube.

●型号:8Fr-36Fr。 Size: 8Fr-36Fr.

创伤引流类 WOUND DRAINAGE SERIES

创伤引流类 WOUND DRAINAGE SERIES

硅胶T型引流管(胆道引流管) Silicone T-shaped Drain (Kher T Tube)

- 100%采用医用级硅胶生产。
 Made of 100% imported medical-grade silicone.
- 带X光显影线,具有良好的生物相容性。
 With X ray opaque line, and good bio-compatibility.

产品编号 Product No.	规格(Fr) Size(Fr)	外 径 O.D.(mm)
PTD101850	10	3.3
PTD121850	12	4.0
PTD141850	14	4.7



硅胶负压引流球 Silicone Reservoir / Bulb

◆ 采用进口医用级硅胶生产,带止回阀。

Made of imported medical grade silicone, with anti-reflux valve.

● 负压球可以和各种类型的引流管配套使用。

Silicone Reservoir is used with different kind of drains.

加压球体使之产生负压,将球体和引流管连接,在密闭环境下,球体内的负压可以稳定地抽吸血液或分泌物。

Squeeze the reservoir to generate negative pressure inside. Blood or secretion will be sucked inside the bulb due to negative pressure.



产品编号 Product No.	规格(ml) _{Size(ml)}	刻度位置(ml) Marked Position(ml)	备注 Remark	
SR100S	100	25/50/75/100		
SR150S	150	25/50/75/100/125/150	Clamp	
SR200S	200	50/100/150/200	available	
SR400S	400	100/200/300/400		
SR100ST	100	25/50/75/100	With extension	
SR200ST	OST 200 50/100/150/200		tube	

硅胶引流套装 Closed Wound Drainage Kit

- 100%采用医用级硅胶生产。
 Made of 100% imported medical-grade silicone.
- 该套装包含各型号负压球,引流管和不锈钢引导针。
 The kit contains different sizes reservoir, drainage tube and stainless steel trocar.
- 创伤引流套装用于引流伤口处的血液和多余液体,从而减少血肿及消除引流死角。

Used to remove blood and excess fluid from wound, so as to prevent or minimize hematomas and eliminate dead space.

📀 如有必要,可与转换接头配用。

Connectors can be equipped with drains if required.



产品编号 Product	配件 Components	备注 Remark
Closed	100/150/200/400ml silicone reservoir	Connectors
Wound Drainage	flat/round fluted/perforated drain	available
Kit	stainless steel trocar	avaliable

不锈钢引导针

Stainless Steel Trocar

- 不锈钢材质,尖端锋利,针体光滑。
 Made of stainless steel, with sharp pinpoint and smooth body.
- ◆ 适当的弯曲角度,便于创伤穿刺。 With a proper angle, which is more convenience of puncture.
- 尾部光滑设计或阶梯设计可供选择。
 With a smooth end or terraced end for choice.
- 可定做双端弯曲引导针。
 Double bendable trocar could be customized.
- ◆ 尺寸: 10Fr / 15Fr / 19Fr, 尺寸可定做。 Size: 10Fr / 15Fr / 19Fr, size could be customized.



<mark>泌尿类</mark> UROLOGY SERIES

泌尿类 UROLOGY SERIES

双腔硅胶导尿管

2-Way Silicone Foley Catheter

- 100%采用医用级硅胶生产,管体柔软透明。 Made of 100% medical grade silicone, the catheter is soft and clear.
- ◆ 球囊具有良好的均衡性和拉伸性,使用过程中更加安全、可靠。 Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.

● 配备不同颜色的色标环以区分型号。

长度: 400mm ± 10mm。

Equipped with different color rings to indentify different sizes. Length: $400mm \pm 10mm$.

产品编号 Product No.	规格(Fr) _{Size(Fr)}	色标号 Color Code	球囊容量(ml) Balloon Volume(ml)	备注 Remark
FC1210S2	12	橙色/白色/White		
FC1410S2	14	绿色/Green		
FC1610S2	16	橙色/Orange		
FC1810S2	18	红色/Red	5–10ml	
FC2010S2	20	黄色/Yellow		
FC2210S2	22	紫色/Purple		
FC2410S2	24	蓝色/Blue		不带导丝 Without stylet
FC1215S2	12	白色/White		
FC1415S2	14	绿色/Green		
FC1615S2	16	橙色/Orange		
FC1815S2	18	红色/Red	5-15ml	
FC2015S2	20	黄色/Yellow		
FC2215S2	22	紫色/Purple		
FC2415S2	24	蓝色/Blue		
FC1630S2	16	橙色/Orange		
FC1830S2	18	红色/Red		
FC2030S2	20	黄色/Yellow	30ml	
FC2230S2	22	紫色/Purple		
FC2430S2	24	蓝色/Blue		



双腔硅胶导尿管(儿童型)

2–Way Silicone Foley Cathter (Pediatric)

- 100%采用医用级硅胶生产,管体柔软透明。
 Made of 100% medical grade silicone, the catheter is soft and clear
- 球囊具有良好的均衡性和拉伸性,使用过程中更加安全、可靠。 Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.
- 配备不同颜色的色标环以区分型号。长度: 310mm ± 10mm。
 Equipped with different color rings to indentify different sizes. Length: 300mm±10mm.

产品编号 Product No.	规格(Fr) _{Size(Fr)}	色标号 Color Code	球囊容量(ml) Balloon Volume(ml)	备注 Remark
FC0603S2	6	粉色/Pink	3ml	
FC0803S2	8	浅蓝色/Light Blue	. .	带导丝
FC1003S2	10	黑色/Black	3ml	
FC0805S2	8	浅蓝色/Light Blue	2 E	With stylet
FC1005S2	10	黑色/Black	3–5ml	

双腔硅胶导尿管(弯头型/开口型)

2-Way Silicone Foley Catheter (Tiemann Tip/Open Tip)

🔷 100%采用医用级硅胶生产,管体柔软透明。

Made of 100% medical grade silicone, the catheter is soft and clear.

- 球囊具有良好的均衡性和拉伸性,使用过程中更加安全、可靠。 Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.
- ◆ 导头可提供弯头型和开口型设计。
 ◆ 长度: 400mm ± 10mm
 Both Tiemann tip and open tip desin are available.
 ◆ 长度: 400mm ± 10mm



产品编号 Product No.	规格(Fr) _{Size(Fr)}	色标号 Color Code	球囊容量(ml) Balloon Volume(ml)	备注 Remark
FC1215S2T	12	白色/White		
FC1415S2T	14	绿色/Green		
FC1615S2T	16	橙色/Orange		
FC1815S2T	18	红色/Red	5-10ml / 5-15ml	不带导丝
FC2015S2T	20	黄色/Yellow	5-1011175-15111	Without stylet
FC2215S2T	22	紫色/Purple		
FC2415S2T	24	蓝色/Blue		
FC2615S2T	26	粉红/Pink		

泌尿类 UROLOGY SERIES

泌尿类 UROLOGY SERIES

三腔硅胶导尿管 3-Way Silicone Foley Catheter

- 100%采用医用级硅胶生产,管体柔软透明。 Made of 100% medical grade silicone, the catheter is soft and clear.
- 球囊具有良好的均衡性和拉伸性,使用过程中更加安全、可靠。
 Symmetrical balloon expands equally in all directions so as to perform its function safely and efficiently.

配备不同颜色的色标环以区分型号。
 长度: 400mm ± 10mm。
 Equipped with different color rings to indentify different sizes.
 Length: 400mm ± 10mm.



单腔乳胶导尿管

Latex Urethral Catheter

- 采用天然乳胶生产,表面做硅胶涂层处理。 Made of natural latex, with silicone coated.
- 管体有红色和黄色可供选择。
 Both red color and yellow color are available.
- 长度: 400mm ± 10mm。 Length: 400mm ± 10mm.
- 尺寸: 6Fr-30Fr。 Size: 6Fr-30Fr.



双腔乳胶导尿管(弯头型)

2-Way Latex Foley Catheter (Tiemann Tip)

- 采用天然乳胶生产,表面做硅胶涂层处理。
 Made of natural latex, with silicone coated.
- · 弯头设计,更容易置入。
 It is easier to be inserted with a Tiemann tip design.
- ◆ 长度: 400mm±10mm。 Length: 400mm±10mm.



产品编号 Product No.	规格(Fr) _{Size(Fr)}	色标号 Color Code	球囊容量(ml) Balloon Volume(ml)	备注 Remark
FC1210L2T	12	白色/White	5-10ml	
FC1410L2T	14	绿色/Green		
FC1615L2T	16	橙色/Orange		7#8 #
FC1815L2T	18	红色/Red	5–10ml 5–15ml 30ml	不带导丝 Without stylet
FC2015L2T	20	黄色/Yellow		
FC2215L2T	22	紫色/Purple		
FC2415L2T	24	蓝色/Blue		
FC2615L2T	26	粉红/Pink		

<mark>泌尿类</mark> UROLOGY SERIES



双腔乳胶导尿管

2-Way Latex Foley Catheter

- 采用天然乳胶生产,表面做硅胶涂层处理。
 Made of natural latex, with silicone coated.
- 有硬阀和软阀可供选择。
 Both plastic valve and rubber valve are available.
- ◆ 长度:儿童270mm±10mm, 成人400mm±10mm。
- Length: Pediatric 270mm \pm 10mm, Adult 400mm \pm 10mm.

三腔乳胶导尿管

3-Way Latex Foley Catheter

- 采用天然乳胶生产,表面做硅胶涂层处理。 Made of natural latex, with silicone coated.
- 有硬阀和软阀可供选择。
 Both plastic valve and rubber valve are available.
- ◆ 长度: 400mm ± 10mm Length: 400mm ± 10mm.



泌尿类 UROLOGY SERIES

产品编号 Product No.	规格(Fr) Size(Fr)	色标号 Color Code	球囊容量(ml) Balloon Volume(ml)	备注 Remark
FC1630L3	16	橙色/Orange		
FC1830L3	18	红色/Red		
FC2030L3	20	黄色/Yellow	30ml	不带导丝
FC2230L3	22	紫色/Purple		Without stylet
FC2430L3	24	蓝色/Blue		
FC2630L3	26	粉色/Pink		

直肠管 Rectal Tube

采用医用级PVC生产,管体透明度高。
 Made of medical grade PVC, the tube is transparent and clear.

- ◆ 长度: 400mm±10mm。 Length: 400mm±10mm.
- 尺寸: 24Fr-36Fr。 Size: 24Fr-36Fr.



备注 Remark	球囊容量(ml) Balloon Volume(ml)	色标号 Color Code	规格(Fr) Size(Fr)	产品编号 Product No.
带导丝	3–5ml	黑色/Black	8	FC0803L2
With stylet	5 611	灰色/Grey	10	FC1003L2
		白色/White	12	FC1215L2
	5–15ml	绿色/Green	14	FC1415L2
		橙色/Orange	16	FC1615L2
不带导丝		红色/Red	18	FC1815L2
Without stylet		黄色/Yellow	20	FC2015L2
		紫色/Purple	22	FC2215L2
		蓝色/Blue	24	FC2415L2
		粉红/Pink	26	FC2615L2
		橙色/Orange	16	FC1630L2
	30ml	红色/Red	18	FC1830L2
		黄色/Yellow	20	FC2030L2
不带导丝		紫色/Purple	22	FC2230L2
Without stylet		蓝色/Blue	24	FC2430L2
		粉红/Pink	26	FC2630L2

<mark>泌尿类</mark> UROLOGY SERIES

PVC导尿管 PVC Nelaton Catheter

◎ 采用医用级PVC生产。

Made of medical grade PVC.

◆ 圆润的侧孔和头端处理,减少对患者黏膜的损害。 Smooth side eyes and closed distal tip design reduce patient's pain during surgical.

亲水性涂层型PVC导尿管 Hydrophilic Coated Nelaton Catheter

- 采用医用级PVC生产。 Made of medical grade PVC.
- 圆润的侧孔和头端处理,减少对患者黏膜的损害。
 Smooth side eyes and closed distal tip design reduce patient's pain during surgical.
- 可提供弯头型设计。
 Tiemann tip design is available.
- 导管表面最亲水涂层处理,遇水后极大增加润 滑度,减少患者痛苦的同事更方便置入。
 The surface of catheter is hydrophilic coated. It becomes extremely smooth while using with normal saline. It reduce the patient's pain and is much easier to be inserted.
- 🔷 提供带水包和不带水包型供客户选择。

Both kits with and without water sachet are available.

产品编号 Product No.	规格(Fr) _{Size(Fr)}	色标号 Color Code	长度(mm±10mm) Length(mm±10mm)	备注 Remark
NC08PF	8	浅蓝色/Light blue	200	女用/Female
NC10PF	10	黑色/Black	200	女用/Female
NC12PF	12	白色/White	200	女用/Female
NC14PF	14	绿色/Green	200	女用/Female
NC16PF	16	橙色/Orange	200	女用/Female
NC08PM	8	浅蓝色/Light blue	400	男用/Male
NC10PM	10	黑色/Black	400	男用/Male
NC12PM	12	白色/White	400	男用/Male
NC14PM	14	绿色/Green	400	男用/Male
NC16PM	16	橙色/Orange	400	男用/Male
NC18PM	18	红色/Red	400	男用/Male
NC20PM	20	黄色/Yellow	400	男用/Male

一次性测温导尿管 Disposable Foley Catheter with Temperature Probe

产品供一次性使用。 This catheter is designed for single use only.

- 产品用于在给患者导尿的同时,测量患者膀胱内尿液的温度。 It is used for detecting urine temperature inside patient's bladder while urine drainage.
- ◆ 温度探头置于导尿管尖端和球囊之间。 The detector is indwelled between tip and balloon of catheter.



产品编号 Product No.	产品型号(Fr) ^{Size(Fr)}	球囊容积(ml) Balloon Volume(ml)	产品长度(mm) _{Length(mm)}
FC0803ST	8	3	358
FC1010ST	10	5-10	358
FC1210ST	12	5-10	458
FC1410ST	14	5-10	458
FC1630ST	16	15–30	458
FC1830ST	18	15–30	458
FC2030ST	20	15-30	458
FC2230ST	22	15-30	458
FC2430ST	24	15–30	458



妇科类 GYNECOLOGY SERIES

宫颈扩张球囊 Cervical Ripening Balloon

- ◎ 在宫颈条件不佳时,富善宫颈扩张球囊为患者提供了一种简单促使 宫颈成熟的方法,同时没有其他促宫颈成熟方法的潜在副作用。 The FUSHAN Cervical Ripening Balloon offers a simple method for cervical ripening in unfavourable cervical conditions, without the potential side effects of other ripening methods.
- 对宫颈条件不佳的患者实施引产是一个颇具挑战性的产科治疗措施, 而促宫颈成熟措施的实施可能会使结局更加复杂。某些药物治疗(如 前列腺素)不仅可能会产生副作用和导致引产失败率上升,而且最终 仍可能需借助剖腹产结束妊娠。

It's challenging to induce labour in patients with unfavourable cervical conditions. Such inductions may be complicated by the cervical ripening method that you use. Certain medical therapies, such as the use of prostaglandins, may have side effects and high failed induction rates while still resulting in undesired caesarean section deliveries.

🔹 富善宫颈扩张球囊能自然地,渐进式扩张宫颈并易化引产。宫颈的成熟与扩张通过球囊在宫颈内口和 外口提供温和、持久的扩张力而实现。

The FUSHAN Cervical Ripening Balloon is engineered to naturally and gradually dilate the cervix and facilitate labour induction. Ripening and dilation are accomplished by the balloons gentle and constant pressure at the level of the cervix from both the internal and external ostia.



1、通过宫颈放入宫颈扩张球囊,使两个球囊 均进入宫颈管。

The Cervical Ripening Balloon through the cervix until both balloons have entered the cervical canal.

4、此时可见阴道球囊位于宫颈外口,用20ml 生理盐水充盈阴道球囊。

You should see the vaginal balloon outside the external cervical os. In ate the vaginal balloon with 20 ml of saline.

● 2-3、用40ml生理盐水充盈子宫球囊。充盈 后,将球囊往后拉至子宫球囊紧贴宫颈内口。

Inflate the uterine balloon with 40 ml of saline. Once the uterine balloon is inflated, pull the device back until the balloon abuts the internal cervical os.

● 5、将球囊置于宫颈两侧后,充入生理盐水一每 个球囊最多80ml,球囊放置后应计时,在进入 分娩活跃期前球囊留置不得超过12小时。欲了 解采用富善宫颈扩张球囊进行产程管理的更多 信息,请参阅产品使用说明书。

Once the balloons are situated on either side of the cervix, add saline-a maximum of 80 ml per balloon. Time the placement of the balloon so that it is in place no longer than 12 hours before active labour is induced. For more information on labour management using the FUSHAN Cervical Ripening Balloon, please refer to the Instructions for Use.



Description of the second seco

Colors and the states and the states of the

MARY METERS IN Include and the Automatical Activity Dependencies Activity Dependencies Activity

APPENDENT BEARING STREET

PARTIES DESCRIPTION

47.









22.



2010

CE



and in MALINY

FSC



FSC



FSC

Provident State

