

Carbide Bur



WE MAKE DENTAL CARE EASIER

Manufacturer : Tealth Foshan Medical Equipment Co.,Ltd

Address:The 2nd Floor, No.4, Qiling Road, Lutang Industrial Zone,
Luocun, ShishanTown, Nanhai District, Foshan City, Guangdong
Province, 528226, P. R. China

Website:www.tealthhandpiece.com

E-mail:info@tealthhandpiece.com

Tel:0086-757-86436722



Carbide Bur

Round											
Head Length	0.50	0.60	0.80	1.00	1.20	1.40	1.60	1.80	2.10	2.30	
Full Length	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	
US NO.	FG 1/4	FG 1/2	FG 1	FG 2	FG 3	FG 4	FG 5	FG 6	FG 7	FG 8	
ISO No	001 / 005	001 / 006	001 / 008	001 / 010	001 / 012	001 / 014	001 / 016	001 / 018	001 / 021	001 / 023	
Inverted Cone					Straight Fissure			Taper Fissure			
Head Length	0.70	0.80	1.00	1.20	1.40	4.00	4.00	4.00	5.00	4.00	
Full Length	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	21.00	19.00	
US NO.	FG 33 1/2	FG 34	FG 35	FG 36	FG 37	FG 56	FG 57	FG 58	FGL 169L	FG 170	
ISO No	010 / 006	010 / 008	010 / 010	010 / 012	010 / 014	107 / 009	107 / 010	107 / 012	171 / 009	168 / 010	
Amalgam Prep				Straight Fissure Crosscut							
Head Length	4.00	4.00	3.00	4.00	2.00	2.00	2.00	4.00	4.00	4.00	
Full Length	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	
US NO.	FG 171	FG 172	FG 245	FG 246	FG 330	FG 331	FG 332	FG 556	FG 557	FG 558	
ISO No	168 / 012	168 / 016	237 / 008	234 / 011	232 / 008	232 / 010	232 / 012	107 / 009	107 / 010	107 / 012	
Taper Fissure Crosscut					Straight Fissure Dome						
Head Length	4.00	4.00	4.00	4.00	4.00	4.00	4.00	4.00	4.00	4.00	
Full Length	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	19.00	
US NO.	FG 559	FG 560	FG 699	FG 700	FG 701	FG 702	FG 703	FG 1156	FG 1157	FG 1158	
ISO No	107 / 014	107 / 016	168 / 009	168 / 010	168 / 012	168 / 016	168 / 021	137 / 009	137 / 010	137 / 012	
Remove Silver, Mercury, Natural Teeth				Extraction And Crown Separation							
Head Length	4.00	4.00	4.00	4.00	6.00	6.00	6.00	6.00	6.00	6.00	
Full Length	19.00	19.00	19.00	19.00	32.00	25.00	28.00	32.00	25.00	28.00	
US NO.	FG 1159	FG 1559	FG 1957	FG 1958	FGOS 700L	FGOS 701L	FGOS 701L	FGOS 701L	FGOS 702L	FGOS 702L	
ISO No	137 / 014	137 / 014	137 / 010	137 / 012	171 / 010	171 / 012	171 / 012	171 / 012	171 / 016	171 / 016	
Head Length	6.00										
Full Length	32.00										
US NO.	FGOS 702L										
ISO No	171 / 016										

Precautions For Use

Sterilization Conditions:
 (1)121°C for more than 20 minutes.
 (2)126°C for more than 15 minutes. Applicable handpiece:
 When using a full-length 25mm, 28mm bur, it must be matched with 4 times the speed, 5 times the speed of the handpiece (micro motor) and must strictly comply with the maximum allowable number of revolutions.
 Note for use: Corrosion may occur if the product is immersed in NaClO/EDTA and other solutions for a long time.

The Background Color Represents The Maximum Rotations Of This Bur:

150,000 RPM
200,000 RPM
300,000 RPM