


PRODUCT INTRODUCTION

BIFIX® Surgical Incision Closure Set


A new surgical suturemethod. This method is suitable for the closure of the epidermis of surgical wounds without causing secondary damage to the skin. The skin's cutting edges are aligned by the proper external force of the product and remain until the wound heals naturally. The unique fastening structure design makes the operation simpler. The precise closing force adjustment improves the aesthetics of the w-ound after healing, and it can be peeled off easily by hand.

COMPOSITION OF PROUDUCT STRUCTURE




CASE PICTURE

• ORTHOPEDICS




• OBSTETRICS



PRODUCT CATALOGUE

• 4cm




MODEL:5HWH-4A
INTENDED USE:Disposable non-invasive incision closing device for the epidermis.
MATERIALS:PSA,PA,PA,PET,HYDROCOLLIOD
PRODUCT COMPONENTS:PSA tape, regulating lever, locking mechanism, and release paper.
VALID DATE:30 MONTHS
STROAGE:This product shall be stored indoors at 5~30 °C out of direct su-
nlight, and shall be free from corrosive gas and well ventilated.
APPILCABLE WOUND LENGTH:4cm
STERILIZATION METHOD:Irradiation sterilization
DIMENSIONS:

length of PSA tape	width of tape	thickness	distance between tapes	length of lever
55mm	17mm	5mm	11mm	36mm

PACKING DETAILS:10sets per box, 20 boxes per carton
Total: 200 sets per carton

• 8cm




MODEL:5HWH-8A
INTENDED USE:Disposable non-invasive incision closing device for the epidermis.
MATERIALS:PSA,PA,PA,PET,HYDROCOLLIOD
PRODUCT COMPONENTS:PSA tape, regulating lever, locking mechanism, and release paper.
VALID DATE:30 MONTHS
STROAGE:This product shall be stored indoors at 5~30 °C out of direct su-
nlight, and shall be free from corrosive gas and well ventilated.
APPILCABLE WOUND LENGTH:8cm
STERILIZATION METHOD:Irradiation sterilization
DIMENSIONS:

length of PSA tape	width of tape	thickness	distance between tapes	length of lever
95mm	17mm	5mm	11mm	36mm

PACKING DETAILS:1 set per individule pack,10sets per box, 12 boxes per c-
arton
Total: 120 sets per carton

• 2cm*3,4cm*1





MODEL:5HWH-10A (Combination pack includes 1X 5HWH-4A,3X 5HWH-2A)
INTENDED USE:Disposable non-invasive incision closing device for the epidermis.
MATERIALS:PSA,PA,PA,PET,HYDROCOLLIOD
PRODUCT COMPONENTS:PSA tape, regulating lever, locking mechanism, and release paper.
VALID DATE:30 MONTHS
STROAGE:This product shall be stored indoors at 5~30 °C out of direct su-
nlight, and shall be free from corrosive gas and well ventilated.
APPILCABLE WOUND LENGTH:10-14cm
STERILIZATION METHOD:Irradiation sterilization
DIMENSIONS:

length of PSA tape	width of tape	thickness	distance between tapes	length of lever
5HWH-4A:55mm, 5HWH-2A:35mm	17mm	5mm	11mm	36mm

PACKING DETAILS:1 set per individule pack,10sets per box, 12 boxes per c-
arton
Total: 120 sets per carton
REMARK:5HWH-2A is not sold separately

• 14cm



MODEL:5HWH-14A
INTENDED USE:Disposable non-invasive incision closing device for the epidermis.
MATERIALS:PSA,PA,PA,PET,HYDROCOLLIOD
PRODUCT COMPONENTS:PSA tape, regulating lever, locking mechanism, and release paper.
VALID DATE:30 MONTHS
STROAGE:This product shall be stored indoors at 5~30 °C out of direct su-
nlight, and shall be free from corrosive gas and well ventilated.
APPILCABLE WOUND LENGTH:14cm
STERILIZATION METHOD:Irradiation sterilization
DIMENSIONS:

length of PSA tape	width of tape	thickness	distance between tapes	length of lever
155mm	17mm	5mm	11mm	36mm

PACKING DETAILS:1 set per individule pack,10sets per box, 12 boxes per c-
arton
Total: 120 sets per carton



BIFIX®
Surgical Incision Closure Set

Zhejiang Haichuang Medical Device Co., Ltd.

Product name: Surgical Incision Closure Set
Trade mark: BI-FIX®
Registration number: 20162650103
License number: 20160009
Packing size: Please see catalogue
Scope of application: Obstetrics and gynecology, general surgery, orthopedics, paediatric surgery, thoracic surgery, orthopedics and minimally invasive surgery
Storage: Store at room temperature, sealed, cool and dry place, avoid direct sunlight
Validity period: 30 months
Address: 1st Floor, No. 146, East Chaofeng Road, Yuhang Economic And Technological Development Zone, Yuhang District, Hangzhou, China
Contact: 01186-18616652496
Fax: 01186-571-87178683
Website: www.zjhycl.com

