LDR Delivery Bed Daphne +86-18345168420 daphne@hrbhengwei.com HW-502-C2

Daphne : daphne



- The back and seat plate, the elevation, the elevation of the auxiliary table and trendelenburg are all completed by electric operation of the motor.

- The back plate and the seat plate adopt the linkage structure. While operating the angle of the back plate, the seat plate also reaches a reasonable angle, making the pregnant women more comfortable.

- The split lifting guardrail is adopted. The outer side of the guardrail is pasted for the control of 3.4.

medical staff. The inner side of the guardrail is pasted for pregnant women to operate the bed by themselves, It is a one-step simple operation.

- The surface is equipped with a floor light function. The floor light is set below the head of the dap

- The base has four universal wheels and adopts central control brake to realize the fixed or moving function of the bed. Both sides of the bed are equipped with brake or brake release pedals, and there are obvious signs to indicate the braking status, which is convenient for operation.

- The backplane has CPR quick release function, so as to improve the rescue speed.

- The outer side of the head of the bed is equipped with a push bracket. At the same time, it is equipped with an infusion stand mounting base and a crash wheel design, so as to better prevent unnecessary damage in case of collision.

- The foot pedal leg support is matched with the leg bed lifting operation button to raise of 1834 lower. The hidden structure saves the time of moving had a support in the leg bed lifting operation button to raise of 1834. lower. The hidden structure saves the time of moving back and forth. It can be opened and him of closed automatically left and right according to the maternal state and delivery state, meeting various needs ap

- The concealed high load-bearing auxiliary table can be quickly disassembled to facilitate the conversion from the maternity bed to the hospital bed and facilitate the operation of medical

- The sitting midwifery handle can be hidden under the hip plate for easy retraction.
- The mattress is made of high-density memory cotton, which is waterproof, soft, anti-corrosion, anti odor and easy to clean.
- It is equipped with leg support, infusion stand, dirt basin, sitting midwifery handle and

Hip lifting seat plate facilitates midwifery and operation of Obstetrics and Gynecology hipheng personnel, and accelerates timely delivery of pregnant women.





