



Supine Position



Prone Position



Lateral Position



Lithotomy Position

- The advanced **unique gel material** is formed at one time.
- The position is **comfortable**, the fixation is firm, no movement, the operation field of vision is exposed clearly, and the operation is convenient.
- **Suitable for all kinds of skin, weight evenly distributed on gel, fluidity gel effectively dispersing pressure, providing excellent pressure release.**
- **No heat storage**, can prevent the generation of bedsore.
- **Long service life**, durable.
- No latex, no allergic reaction, non-conductive, non flammable, non adsorption odor, **X-ray transparent**.
- **Easy to disinfect, easy to clean**, with good flexibility and pressure resistance, it is very suitable for long-time, high-intensity surgical needs.



| Donut Head Pad | |
|----------------|-----------|
| Item | Size (CM) |
| H015 | 12*5*2.5 |
| H016 | 15*6*4 |
| H020 | 20*7*5 |
| H021 | 21*9*5 |
| H030 | 25*11*7 |

Support and fixed protection for head and neck during supine surgery

| Horseshoe Head Pad | |
|--------------------|-----------|
| Item | Size (CM) |
| H113 | 12*5*2.5 |
| H115 | 15*5*3 |
| H120 | 20*7*7 |
| H121 | 20*7*5 |



For the support of the head and face during lateral and prone position surgery, convenient for anesthesia cannula.



| Ophthalmic Head Pad | |
|---------------------|-----------|
| Item | Size (CM) |
| H118 | 18*6*4 |
| H119 | 23*22*4 |
| H222 | 28*25*6 |

For the support and fixation of the head and neck during ophthalmic surgery



| Prostrate Head Pad | |
|--------------------|-----------|
| Item | Size (CM) |
| H273 | 26*22*12 |
| H274 | 27*23*13 |
| H275 | 28*24*14 |

Designed for support of the face during prone position surgery, convenient for anesthesia intubation, keeping breathing smooth.

| Prostrate Head Pad | |
|--------------------|-----------|
| Item | Size (CM) |
| H220 | 23*20*11 |
| H221 | 23*20*5.5 |



Designed for support of the face during prone surgery, convenient for anesthesia intubation, keeping breathing smooth.



| Prostrate head pad | |
|--------------------|---------------------|
| Item | Size (CM) |
| H237 | 25*22*14 |
| H238 | 25*22*(15-18) |
| H214 | 31*27*17 Adjustable |
| H218 | 23*20* (14-11) |
| H213 | 23*20*13.5 /14.5 |
| H212 | 24*21*15.5 / 16.5 |

Support and protection of the face during prone surgery, adjustable, convenient for anesthesia intubation, and a mirror on the bottom to observe the patient during the operation.



| Arm/Leg Fixed Pad | |
|-------------------|-----------|
| Item | Size (CM) |
| L040 | 40*15*4 |
| L050 | 50*15*4 |
| L060 | 60*15*4 |
| L052 | 52*16*5.5 |

For the fixation and protection of limbs

| Head Fixation Pad | |
|-------------------|--------------|
| Item | Size (CM) |
| H224 | 24*22*7 |
| H223 | 29*25*8.5 |
| H233 | 20*20φ13*7-6 |
| H232 | 25*20*3/8 |



For the support and fixation of the head and neck during ophthalmic surgery



| Thyroid Pad | |
|-------------|-------------|
| Item | Size (CM) |
| H260 | 48*39*9 |
| H261 | 60*45*10-14 |
| H262 | 45*36*4 |

For thyroid and neck surgery, strengthen the support and protection of the back of the neck.



| Contoured Supine Head Pad | |
|---------------------------|-----------|
| Item | Size (CM) |
| H215 | 16*10*5.5 |
| H225 | 20*15*5 |
| H226 | 20*15*10 |
| H231 | 18*18*8/4 |

Support and fixed protection for head and neck during supine surgery.

| Semicircular Pad | |
|------------------|-----------|
| Item | Size (CM) |
| C130 | 30*10*4 |
| C140 | 30*15*7 |
| C145 | 42*12*10 |
| C150 | 50*15*7 |
| C151 | 50*10*7 |
| C152 | 50*10*4 |
| C160 | 50*15*10 |
| C161 | 50*15*11 |



Used to protect the compression point during surgery, relieve pressure on the shoulders, chest and abdomen, and protect the nerve tissue.



| Column Shaped Pad | |
|-------------------|-----------|
| Item | Size (CM) |
| C060 | 20*8*10 |
| C070 | 20*12*10 |
| C080 | 40*12*10 |
| C090 | 50*12*10 |

It is used to protect the compression point during surgery, relieve pressure on the shoulders, chest and abdomen, and protect the nerve tissue.

| Slanting shaped pad | |
|---------------------|-----------|
| Item | Size (CM) |
| L225 | 25*15*6-2 |
| L230 | 30*18*7-3 |



Suitable for supine position, lateral position and prone position operation.



| Heel Pad | |
|----------|-----------|
| Item | Size (CM) |
| L016 | 16*9*5 |
| L017 | 19*11*7 |

Used to protect the heel and ankle.

| Hip Pad | |
|---------|-----------|
| Item | Size (CM) |
| L003 | 31*32.5*5 |
| L005 | 50*40*1.5 |
| L006 | 40*40*2 |
| L007 | 40*40*3 |
| L010 | 40*40*4 |
| L014 | 50*40*4 |
| L001 | 40*40*5 |
| L002 | 50*40*5 |



Used to protect the buttocks and sacrum during supine surgery, which can reduce the pressure to the compression point and prevent the formation of the pressure ulcers.



| Strap Pad | |
|-----------|-----------|
| Item | Size (CM) |
| L125 | 25*5*1 |
| L132 | 32*5*1 |
| L145 | 45*7*1 |
| L160 | 60*7*1 |
| L161 | 60*11*1 |
| L162 | 60*14*1 |

Protect the fixed torso and joints for easy surgery and avoid injury.



| Prone Position Pad | |
|--------------------|-------------|
| Item | Size (CM) |
| C260 | 63*46*16 |
| C261 | 62*50*13 |
| C262 | 62*50*16/13 |

It is designed to relieve chest pressure during prone position and maintain a smooth breathing cycle.

| Prone Position Pad | |
|--------------------|-----------|
| Item | Size (CM) |
| C263 | 63*46*18 |



It is designed to relieve chest pressure during prone position and maintain a smooth breathing cycle.



| Prone Position Pad | |
|--------------------|-----------|
| Item | Size (CM) |
| C253 | 80*50*11 |
| C254 | 60*40*11 |

It is designed to relieve chest pressure during prone position and maintain a smooth breathing cycle.



| Prone Position Pad | |
|--------------------|-------------|
| Item | Size (CM) |
| C266 | 75*50*25-33 |

The stent is made of carbon fiber and can pass X-rays. It is mostly used for prone surgery under perspective.

| Prone Position Pad | |
|--------------------|-------------------|
| Item | Size (CM) |
| C264 | 66* (33-40) *22 |



It is designed to relieve chest pressure during prone position and maintain a smooth breathing cycle.



| Prone Position Pad Set | |
|------------------------|------------|
| Item | Size (CM) |
| C257 | 30*12*10/4 |
| C258 | 40*12*10/4 |
| C259 | 50*12*10/4 |

A special protective pad for prone surgery, reduces the pressure on the chest and abdomen and provides protection for the compression points of the nerve tissue.



| Prone Position pad | |
|--------------------|-----------|
| Item | Size (CM) |
| C281 | 56*20*15 |
| C283 | 65*20*12 |
| C284 | 54*20*12 |
| C282 | 67*20*15 |

A special protective pad for prone surgery, reduces the pressure on the chest and abdomen and provides protection for the compression points of the nerve tissue.

| Tunnel Pad | |
|------------|-----------|
| Item | Size (CM) |
| C273 | 64*46*20 |
| C275 | 64*46*25 |
| C276 | 63*40*20 |



It is used for the protection of lower limbs and the prevention and treatment of pressure sore during the operation of lateral position. Suitable for people with large weight.



| Lateral Position Pad | |
|----------------------|-----------|
| Item | Size (CM) |
| C270 | 70*50*17 |
| C271 | 70*50*15 |
| C272 | 72*50*10 |

It is designed for use in lateral position surgery to relieve pressure on the shoulders and arms and protect nerve tissue.



It has strong versatility and is widely used. It can reduce the pressure of the pressure points and prevent the formation of the pressure ulcers.

| Universal square pad | |
|----------------------|-----------|
| Item | Size (CM) |
| L301 | 10.5*62*1 |
| L302 | 16*15*2.5 |
| L303 | 16*60*1 |
| L304 | 27*20*1.5 |
| L305 | 40*24*1.5 |
| L306 | 50*23*1.5 |
| L307 | 50*30*1.5 |
| L308 | 50*40*1.5 |
| L309 | 50*50*1.5 |
| L310 | 53*45*1 |
| L311 | 57*35*1 |
| L312 | 63*45*1.5 |
| L313 | 30*20*2.5 |
| L314 | 18*18*1.5 |
| L315 | 59*12*1.5 |
| L316 | 40*40*1.5 |
| L317 | 75*50*2 |
| L318 | 45*7*1.2 |
| L319 | 25*5*1.2 |
| L320 | 32*5*1.2 |

| Universal square pad | |
|----------------------|-----------|
| Item | Size (CM) |
| L321 | 15*5*1 |
| L322 | 60*7*1.2 |
| L323 | 35*35*1.5 |
| L357 | 15*12*1.8 |
| L324 | 15*12*1.5 |
| L325 | 30*20*1.5 |
| L326 | 40*30*1.5 |
| L327 | 63*40*1.3 |
| L328 | 40*20*4 |
| L329 | 40*40*2.5 |
| L346 | 40*40*2 |
| L347 | 40*40*3 |
| L330 | 40*40*4 |
| L331 | 50*20*2 |
| L332 | 50*20*3 |
| L333 | 50*20*4 |
| L348 | 50*30*2 |
| L334 | 50*30*3 |
| L335 | 50*30*4 |
| L336 | 50*40*2 |