## ELECTIC HOSPITAL BED TEST REPORT

## FOSHAN DONGFANG MEDICAL EQUIPMENT MANUFACTORY (LTD.)

NO: 20250616

PRODUCT NAME		ELECTRIC HOSPITAL BED		NO.	2505023	QTY	6PCS	TEST QTY	2PCS	STANDAR	YZB/	粤 2371-2014	
MODEL		FS3230WF52		CHECKER	GUIJUAN ZENG	DATE	2025-06-16	EXAMINATOR	YIZHU ZHANG	YIZHU ZHANG DATE 2		025-06-16	
N o	Standar terms No.	Name of item	Technology specification							Result		Single conclusion	
1	4.2	Appearance structure requirement	2) ch obvious 3) app scratch, 4) the	somed surface shous s scratch. pearance of powder no dehiscense and seam should be syn	orderliness, no irregula ald be smooth and flatted coating parts should be no obvious repair mark mmetrical, no crack, no which write on the hos	no	Pass		Pass				
	4.3	Dimension	a) L:185cm~215cm; W:85cm~105cm; H:40cm~80cm							Pass		Pass	
			b) MAX angle of backrest is no less than 70°							Pass		Pass	
2			c) MAX angle of footrest is no less than 40°							Pass		Pass	
			d) when legrest up to $40^\circ$ and footrest up to the highest angle, the angle between them is no less than $130^\circ$ .						nan	Pass		Pass	
			e) backrest, legrest and footrest could adjust freely in the adjustable angle.							Pass		Pass	
3	4.4	Structure requiremnt	<ol> <li>bed board and bed frame could be replaced so that easy to assemble and disassemble.</li> <li>it's firm between bed board and bed frame, no loosened</li> <li>after connecting all parts of the bed, all operational parts could work flexible without block and loosened.</li> </ol>						and	Pass		Pass	
	4.5	Operational capability	After loading 175kg dead load, the operational system should work exactly and reliably, can stop on any time. When it elevates to any middle place, the speed of falling by itselft is no less than 10mm/h.					ny	Pass		Pass		
4			The operational spare part should design on the place which can't be easy touched.  When bed operates, it should not pause or jump suddently.  There is at least one emergency stop button on the bed. All damage action will stop when press the button.					n.	Pass Pass Pass		Pass		
5	4.7	Stability		Assembled bed should be stable, the height difference of four angles is no more than 10mm.						Pass		Pass	
	4.8	Wheels	Diameter of wheel is no less than φ75mm.						Pass 11 1		Pass		
			Whee	Wheels work flexibly and stably, locomotivity is no more than 50N when it loads nothing.					/	Pass		Pass	
6			Whee Ther	els are stable assem	bled under the bed.  ng wheels on any one p			he locking system shou		Pass			
7	A. 2. 34	Continuous leakage electric	Moto User:	or: normal condition	<b>J</b>		condition:≤0.5m		one mulfunction co User: normal con	ndition: ≤0.1mA; ndition:≤0.5mA dition:≤0.1mA; just o condition:≤0.5mA		Pass	
8	A. 2. 36	Electrolyte strength The part A-a <sub>2</sub> S. B-a of the hospital bed could bear 400Vand 50Hz testing pressure for 1mim, without breakdown and flashover.											
	Test result	According to the test, all qualifications reach the manufacturer's standard request of YZB/學 2371-2014 <medical bed="" electric="" hospital="">. The goods pass the test, allowed to leave the factory.  Quality managing department: Mindan Zhang Date: 2025-06-16</medical>											