



Shenzhen Rainbow Technology Co., Ltd.

Add: Building C, Fenda Technology Innovation Park, Gushu Area,
Xixiang, Baoan District, Shenzhen, Guang dong province, China.

Web: www.rainbowdo.com

Tel: (86)-755-2301-0360

Phone: (86) 176-2036-3135

Email: rainbow10@rainbowdo.com



Wireless Breast massager

深圳市瑞音博科技有限公司

Shenzhen Rainbow Technology Co., Ltd.

Helpful
Effectively improve milk
flow relieve clogged ducts

RB-019 Warming Lactation Massager



Breast massager

Quiet comfort / Thin and invisible



- 01 | Content
- 02 | Introduction
- 03 | Parameter
- 04 | Product features
- 05 | Comparison
- 06 | Operation
- 07 | Package includes
- 08 | OEM&ODM
- 09 | Contact Us

Lactation Massager Can Help You

- Relieve Pumping Discomfort
- Reduce Breastfeeding Issue
- Reduce Pumping Time



- *Are you often low on milk when breastfeeding?*
- *Are you often blocked while breastfeeding?*
- *Do you often experience breast tenderness?*

Warming Lactation Massager can help you effectively improve these problems!

- ✓ **Unique Design**


Unique wing design which is different from the market

- ✓ **Material Safety**

Silicone on the contact surface is BPA Free

- ✓ **Good Concealment**

Thin and light design for easy concealment



Reverse--Silicone Material



Front - ABS Plastic



Ergonomic desian



Combining Heat & Vibration





>70
MINS

Longer massage time

The scraping board design-More contact points for massage precisely

2-in-1 Lactation Massager

Our Warming Lactation Massage Pads bring together heat and vibration to make your pumping and nursing experience better



The warming lactation massager helps with breastfeeding ailments

- | | | | |
|---|-------|---|-------|
| 1 | ~~~~~ | 4 | ~~~~~ |
| 2 | ~~~~~ | 5 | ~~~~~ |
| 3 | ~~~~~ | 6 | ~~~~~ |

Effective lactation tool with 6 different vibration modes

WARMING LACTATION MASSAGER

6 Vibration Modes Control



FOOD GRADE
SILICONE



RELIEVE CLOGGED
DUCTS - MASTITIS



HEATING MODE



LATCHING



LETDOWN

Lactation Massager Can Help You

- Relieve Pumping Discomfort
- Reduce Breastfeeding Issue
- Reduce Pumping Time



BLOCKAGE



PAIN




SWELLING



LESS MILK



Reduce Pumping Discomfort Save Time



Vibration mode
Heating mode

Before




After



3 heating levels



6 vibration modes



Parameter




Model	RB-019
Product net weight	50g
Power Source	5V 1A
Motor Vibration	9000rpm \pm 10%
Working temperature	38°C-48°C
Working time	Around 1.5 hours
Battery capacity	500mah
Cable	Type-C
Waterproof	IPX5
CTN Size	370*250*330MM

FUNCTIONAL COMPARISON


		
		
OTHER MASSAGER	2022 NEW	HEATING PAD
3 Levels	Adjustable Vibration 6 Levels	NO
NO	Warming Fast	Warming
Rechargeable	Rechargeable	NO
NO		NO

Ways to Use


Unblock plugged ducts , ergonomic desian




Use the massage pad and rub breast towards the nipple



Use the massage pad and rub breast towards the nipple



Gently massage breast using the narrow part of massager



Use the narrow part of the massager and apply light pressure to the clogged ducts


Monther and Baby Care-make a set



RB-018



RB-019











COMPANY PROFILE

Shenzhen Rainbow are an original manufacturer:

- Audited by Kmart & Costco
- With BSCI /ISO9001
- Professional in OEM & ODM, with our ID department
- Free ODM service, with free samples for evaluation



Contact US




Joyce

 rainbow10@rainbowdo.com

 (86)-130-7787-1959

 (86)-755-2301-0360

 Building C, Fenda Technology Innovation Park, Xixiang, Baoan District, Shenzhen,
Guang dong province, China.



Shenzhenrainbow



Shenzhen Rainbow



Shenzhen Rainbow Technology Co.,Ltd



Company website



Facebook