

Nucleic Acid Extraction Kit

Product Introduction

It's used for extraction of viral DNA/RNA from Samples in a rapid, highly sensitive and efficient way. The purified nucleic acid is suitable for reverse transcription, PCR, RT-PCR and other common downstream experiments.

Product Features

Nontoxic

Without using expensive equipment; Non-exposure to toxic reagents, like phenol and chloroform; can be operated indoor and outdoor.

Wide


Be suitable for DNA/RNA extraction of serum, plasma, urine and other cell-free body fluids, as well as virus culture fluid, cell culture fluid, virus preservation fluid and other samples.

Simple

Simple to operate and automate, Suitable for high-throughput extraction.

Technology Roadmap

Magnetic bead method for nucleic acid extraction



Product Information

Product	Specification	Remark
Nucleic Acid Extraction Kit (Manual)	50 Tests	
Nucleic Acid Extraction Kit (Automatic)	48 Tests	Using with nucleic acid extractor

