

## Rapid Test Kit(Colloidal Gold)



### Introduction:

Hospitals, customs, health and epidemic prevention departments, pharmacies, personal and family use, etc.

### Parameters:

| No. | Name                                                                | Packing Parameters                                                                                         |
|-----|---------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|
| 1   | Chikungunya IgG/IgM Combo Rapid Test Kit (Colloidal Gold)           | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 2   | Combined Detection Kit For Influenza A/B (Colloidal Gold)           | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 3   | Dengue Rapid Test Kit (Colloidal Gold)                              | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 4   | Helicobacter Pylori Antibody Rapid Test Kit (Colloidal Gold)        | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 5   | Malaria P.Falciparum/P.Vivax Rapid Test Kit (Colloidal Gold)        | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 6   | Typhoid Rapid Test Kit (Colloidal Gold)                             | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 7   | 6-Monoacetylmorphine Test Kit (Colloidal Gold)                      | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 8   | Amphetamine Test Kit (Colloidal Gold)                               | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 9   | Benzodiazepine Test Kit (Colloidal Gold)                            | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 10  | Buprenorphine Test Kit (Colloidal Gold)                             | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 11  | MDMA Test Kit (Colloidal Gold)                                      | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 12  | Cocaine Test Kit (Colloidal Gold)                                   | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 13  | Ketamine Test Kit (Colloidal Gold)                                  | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 14  | Morphine Test Kit (Colloidal Gold)                                  | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 15  | Tetrahydrocannabinolic Acid Test Kit (Colloidal Gold)               | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 16  | Luteinizing Hormone (LH) Test Kit (Colloidal Gold)                  | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 17  | Follicle-Stimulating Hormone (FSH) Test Kit (Colloidal Gold)        | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |
| 18  | Beta-Human Chorionic Gonadotropin (β-HCG) Test Kit (Colloidal Gold) | Packaging information: 25 pieces/box, 36 boxes/carton<br>Packaging size: 630*460*250mm; Gross weight: 12kg |

## Disposable Virus Sampling Tube Kit



### Introduction:

In the current epidemic situation, virus sample collection is an important part of virus detection. The single-use virus sampling tube can collect, transport, Inactivate and store virus samples from specific parts of the human body.

### Advantage:

- ①. Easy to operate and use.
- ②. Adding virus stabilizing ingredients can maintain the activity of the virus in a wide temperature range and reduce the decomposition rate of the virus (Active type).
- ③. Contains virus cleavage and virus nucleic acid preservation solution, which can quickly cleave the virus to release nucleic acid and store the nucleic acid stably (Inactive type).

### Test Procedures:




1. According to different sampling requirements, use a sampling swab to sample at the corresponding part.




2. Place the swab after collecting the sample in a sterile sampling tube.



3. Put the swab into the preservation solution and break it from the breaking point.



4. Tighten the sterile sampling tube cover.



5. Mark the relevant sample information on the label of the sampling tube.

### Product Components:

| Item                                   | Packing Information                                                 |
|----------------------------------------|---------------------------------------------------------------------|
| 10ml tube with 2ml medium; 50pcs/board | Total of 1200pcs, Packing size: 605*485*255mm, gross weight: 11.5kg |
| 10ml tube with 3ml medium; 50pcs/board | Total of 1200pcs, Packing size: 605*485*255mm gross weight:12.5 kg  |
| Remark                                 | The swab can be optional                                            |