
BIOBASE®

BIOBASE

E-mail: export@biobase.com www.biobase.com

Mixer Multi-Tube Vortex Mixer



Features:

- * Up to 96 samples can be processed at one time, and the rotation speed can reach 2500rpm.
- * Run time and pause time can be set.
- * Configure various types of sponge test tube racks.
- * Button-type jog mode, double switch design.
- * Power-off memory function, saving setup time.
- * Accelerates uniformly, effectively avoiding sample splashing.
- * Three working modes of jog, continuous and intermittent pulse can be selected at will.
- * LED display speed and time, simple operation panel, precise control by microprocessor.

Technical Parameters:

Model	BHY-2500	
Speed Range	500~2500rpm	
Speed Regulation Accuracy	±1rpm	
Amplitude	3.6mm	
Time Setting Range	1s~99h59min	
Pulse Interval Timing Range	1s~99s	
Pulse Operation Timing Range	1s~99min59s	
Consumption	120W	
Power Supply	100~230V, 50/60Hz	
Standard Accessory	I, II, III, IV	
External Size (L*W*H)	438*301*429mm	
Package Size (W*D*H)	530*380*510mm	
Gross Weight	27kg	

Accessories:

Model	Pore Diameter	Compatible Type	Size
1	48-Well Virus Sampling Tube Rack	20ml or 30ml virus sampling tube	272*205*50mm
П	96-Well Virus Sampling Tube Rack	ck 7ml or 10ml virus sampling tube 10ml or 15ml centrifuge tube	296*202*50mm
III	Replaceable Tray Pad	1	302*212*25mm
IV	Ф25~Ф29mm	50ml centrifuge tube	300*210*45mm

113