


A Individual Package


Package for 1 Test:



Packing Quantity	Package Size	Gross Weight
1 boxes/carton	300*250*200mm	1kg
12 boxes/carton	650*440*430mm	11kg
16 boxes/carton	840*560*260mm	14kg

B Normal Package

What In the Box?




Packing Quantity	Package Size	Gross Weight
1 boxes/carton	300*200*150mm	0.5kg
22 boxes/carton	460*350*350mm	7.5kg
50 boxes/carton	650*440*430mm	15kg

HOW TO USE

Sample collection:



Sample testing:



The Results of Judgment:

Positive: T-line and C-line present within 15-20 minutes after sample addition.

Negative: only C-line, no T-line, occurs within 15-20 minutes after dropping the sample.

Invalid: As long as there is no line C, this test card is invalid. Please use another test card to retest.

Do not check the results after 20minutes.

