

Cocaine Test Kit (Colloidal Gold)



Applications:

Drug control system departments, traffic control system departments, professional institutions, hospitals, health and epidemic prevention departments, etc.

Technical Parameters:

| Methodology | Colloidal gold method |
|---------------------|---------------------------|
| Sample Type | Urine |
| Storage Condition | 4-30°C |
| Period of Validity | 24 months |
| Sample Volume | 3 drops, about 80µL |
| Detection Duration | 5-10 minutes |
| Detection Limit | 300ng/mL |
| Package Information | 20pcs/box, 30boxes/carton |
| Packing Size | 310*330*360mm |
| Gross Weight | 8kg |