

## Dengue Rapid Test Kit (Colloidal Gold)



### Applications:

Hospitals, customs, health and epidemic prevention departments, drugstore and personal home usage, etc.

### Technical Parameters:

Methodology	Colloidal gold method
Sample Type	Whole blood/serum/plasma
Storage Condition	4-30°C
Period of Validity	24 months
Sample Volume	5µL serum/plasma; 10ul Whole blood
Detection Duration	15-20 minutes
Performance Index	The lowest detection limit for IgG is 1:40; the lowest detection limit for IgM is 1:20.
Package Information	20pcs/box, 30boxes/carton
Packing Size	310*330*360mm
Gross Weight	8kg