

## Follicle-Stimulating Hormone (FSH) Test Kit (Colloidal Gold)



### Applications:

Professional institutions, hospitals, clinics, drugstore and personal home usage, etc.

### Technical Parameters:

Methodology	Colloidal gold method
Sample Type	Urine
Storage Condition	4-30°C
Period of Validity	24 months
Sample Volume	3 drops, about 80µL
Detection Duration	5-10 minutes
Detection Limit	25mIU/mL
Package Information	20pcs/box, 30boxes/carton
Packing Size	310*330*360mm
Gross Weight	8kg