

Malaria P.Falciparum/P.Vivax Rapid Test Kit (Colloidal Gold)



Applications:

Hospitals, customs, health and epidemic prevention departments, drugstore and personal home usage, etc.

Technical Parameters:

Methodology	Colloidal gold method
Sample Type	Whole blood
Storage Condition	4-30°C
Period of Validity	24 months
Sample Volume	5ul Whole blood
Detection Duration	15-20 minutes
Performance Index	100 Plasmodium/uL
Package Information	20pcs/box, 30boxes/carton
Packing Size	310*330*360mm
Gross Weight	8kg