



## Multi-Tube Vortex Mixer




I



II



III



IV

### Features:

- \* Up to 96 samples can be processed at one time, and the rotation speed can reach 2500rpm.
- \* Run time and pause time can be set.
- \* Configure various types of sponge test tube racks.
- \* Button-type jog mode, double switch design.
- \* Power-off memory function, saving setup time.
- \* Accelerates uniformly, effectively avoiding sample splashing.
- \* Three working modes of jog, continuous and intermittent pulse can be selected at will.
- \* LED display speed and time, simple operation panel, precise control by microprocessor.

### Technical Parameters:

Model	BHY-2500
Speed Range	500~2500rpm
Speed Regulation Accuracy	±1rpm
Amplitude	3.6mm
Time Setting Range	1s~99h59min
Pulse Interval Timing Range	1s~99s
Pulse Operation Timing Range	1s~99min59s
Consumption	120W
Power Supply	100~230V, 50/60Hz
Standard Accessory	I, II, III, IV
External Size (L*W*H)	438*301*429mm
Package Size (W*D*H)	530*380*510mm
Gross Weight	27kg

### Accessories:

Model	Pore Diameter	Compatible Type	Size
I	48-Well Virus Sampling Tube Rack	20ml or 30ml virus sampling tube	272*205*50mm
II	96-Well Virus Sampling Tube Rack	7ml or 10ml virus sampling tube 10ml or 15ml centrifuge tube	296*202*50mm
III	Replaceable Tray Pad	/	302*212*25mm
IV	Φ25~Φ29mm	50ml centrifuge tube	300*210*45mm