

Roller Mixer



Introduction:

Roller mixer can be used for all kinds of common test tubes. It is especially suitable for mixing blood, viscous liquids and liquid-solid suspensions, prevention of blood coagulation and immune precipitation etc.

Features:

- * LCD display real-time indicates time and speed.
- * Time setting range of 1~99h59min.
- * Adjustable speed range of 20~100rpm, DC brushless motor, stable performance, low noise, maintenance free.
- * Gentle rolling process, continuous operation mode, safe and reliable.
- * Compatible with the different sizes of test tubes.

Technical Parameters:

Model	MX-T6-D
Mixing Mode	Rocking & Rolling
Speed Range	20~100rpm
Rocking Amplitude	24
Number of Roller	6
Roller Length	280mm
Roller Interval	5mm
Max. Loading Weight	4kg
Timing Range	1~99h59min
Operation Mode	Timing or continuous
Speed Display	LCD
Timer	Yes
Timer display	LCD
Motor Type	DC brushless motor
Ambient Temp.	5~40°C
Relative Humidity	80%
Power Consumption	15W
Power Supply	100~240V, 50/60Hz
External Size(W*D*H)	475*265*110mm
Package Size(W*D*H)	540*475*215mm
Net Weight	5.5kg
Gross Weight	7.2kg