

Tetrahydrocannabinolic Acid Test Kit (Colloidal Gold)



Applications:

Drug control system departments, traffic control system departments, professional institutions, hospitals, health and epidemic prevention departments, etc.

Technical Parameters:

Methodology	Colloidal gold method
Sample Type	Urine
Storage Condition	4-30°C
Period of Validity	24 months
Sample Volume	3 drops, about 80µL
Detection Duration	5-10 minutes
Detection Limit	50ng/mL
Package Information	20pcs/box, 30boxes/carton
Packing Size	310*330*360mm
Gross Weight	8kg