

## Transdermal Diffusion Meter



### Applications:

The unit can objectively reflect the speed and degree of penetration of drug preparations through the skin of living animals in the specified solvent, and scientifically screen effective drugs with transdermal absorption conditions.

It is one of the internationally recognized standard detection methods for drug transdermal release.

### Technical Parameters:

Model	BK-TP6
Transdermal Cup	6pcs
Capacity of Cell	15ml(standard)
Temperature Control Range	30-40°C
Temperature Control Precision	±0.1°C
Stirring Speed Range	100-800r/min
Speed Precision	± 5%
Power Supply	Standard: 220V, 50/60Hz; optional: 110V, 50/60Hz
Power	360W
Packing Size	420*420*380mm
Gross Weight	12kg