

Negative Pressure Wound Therapy System

The details about NPWT Device:

- **Weight:** 0.62 kg/pcs. approx.
- **Size:** L13.5cm*H12cm*W10.2cm—5.7cm
- **Mobility:** Portable, optional carrying canister : Volume 140ml./450ml.
- **Therapy Modes:** Continuous Therapy; Dynamic Therapy; Intermittent Therapy.
- **Pressure Options:** -50 mmHg to -225 mmHg.
- **Charging Time:** 3 hours. for full charge.
- **Battery Life:** 25 hours, depending on configuration.
- **Accessory:** a power adapter and a use manual.
- **Specification of adapter:** Input AC 100-240VAC, Output DC18V, 2A. Voltage: DC18V
- **Package:** individually wrapped in a box with a machine.
- **Operation Language:** English or Others language can be customized such as Spanish.



Indications:

The Rainhome System for Negative Pressure Therapy is a portable, battery-powered pump designed to generate negative pressure or suction to remove exudates, infectious material, and tissue debris from the wound bed to aid the healing process. .

- Full thickness wounds.
- Chronic wounds that DO NOT heal despite receiving proper treatment.
- Venous Ulcers
- Diabetic foot ulcers
- Neuropathic Ulcers
- Bedsores
- Acute wounds such as dehiscent surgical wounds, traumatic wounds, and other open ulcers that may improve when treated with NPWT. Our products have been application in world hospitals.



Rainhome NPWT Dressing Kit



It is including black/white PU foam, tubes with the pad, adhesive film, I/Y connector and clip. Each accessory is packed individually and sterilized.

1; Black/white PU foam

The size : 5 × 5 × 3cm ; 15 × 5 × 3cm ; 10 × 7 × 3cm ; 15 × 10 × 3cm ; 20 × 10 × 3cm , these sizes are in stock for sale.

And extra large sizes: 25 × 15 × 3cm and 30 × 20 × 3cm for order. MOQ: 200sets

2; Adhesive film

The size: 30cm × 20cm

3; Tubes with the pad

Length of tube: 30cm/60cm; Diameter of pad: 10cm

4; I/Y connector


Length: 5.6cm, Dia:0.4~1.0cm

5; Clip

It is used to absorb moderate to heavy exudate associated with moderate-deep wounds and is customizable to any wound size. Rainhome dressing kit is used in conjunction with NPWT device to promote wound healing by removal of wound exudates, debris and infectious material.

All components are packaged together and have been sterilized by gamma irradiation or ethylene oxide (EO).

Package for NPWT device, canister and dressing kit:



Size/carton: 46×35×37cm; 4.5kg/carton.

3pcs NPWT device/carton, 24pcs canisters/carton, 28sets dressing kits/carton.

Here are another option for foam dressing



NPWT dressing white PVA foam

Rainhome white foam is an open cell polyvinyl alcohol (PVA) foam, with hydrophilic properties and a good level of absorption, allowing a reliable removal of exudate from the wound.

Rainhome white foam is available in two sizes:

- 10 × 7 × 1 cm
- 15 × 10 × 1cm
- 20 × 10 × 1 cm

The products have been sterilized by gamma irradiation.

Package:

White PVA foam with tubes is packed individually and I/Y connector & clip are packed in a plastic bag. They are packed together in a box. 10boxes/carton (46*34*36cm/carton).

Order Process: Receiving your PO, offer PI for you, Receiving payment and arrange shipment for you during 2-3days.

Preparing time for the products: 2-3days

Delivery time: 5-7days

Shipping terms: DHL, UPS, FedEx, Special Air Line or By Sea.

Payment terms: Pay it online or T/T by bank or others.

We are the manufacturer for these items and there are enough stock for order and ready to ship.



Manufacturer: Guangzhou Rainhome Pharm&Tech Co., Ltd

Address: On 5th Floor, Building 6, No.10 Yongsheng Road, Huangpu District, Guangzhou City of Guangdong Province, China, 511356

Email: ds03@rhkj.com.cn

Tel: +86-13632492562 (WhatsApp/Wechat)

http://www.rainhomepharm.en.alibaba.com