



DYNAMIC PRECISION
EMBEDDING TRUST



Jin Ming Medical Device Co., Ltd.

Address: Room 601, 602, 6th Floor, Building 4, No. 12 Zhicheng Road, Nexus Valley, Songshan

Lake Park, Dongguan City, Guangdong Province, China

Tel: +86 769 22803823

Email: info@jinmingmed.com

Website: www.jinmingmed.com

Copyright © 2025 Jin Ming Medical Device Co., Ltd. All rights reserved.

Trauma Plate Series

Company profile



We are a China-based company dedicated to the research, development, manufacturing, and sales of specialized orthopedic surgical instruments and implants. Our product portfolio encompasses nine major categories: Trauma Plate Series, Locking Plate Series, Cranio-Maxillofacial Plate Series, intramedullary Nail Series, Spine Series, External Fixation System Series, Joint Series, Sports Medicine Solutions, and Power Tools.

We cooperate with world-class medical device companies and customize production to produce a variety of products that meet the laws and regulations and requirements of different countries.

"Insisting on quality first".Pioneering,pragmatic, enterprising and innovative, and strive to move forward.



CONTENTS

- 1** **SMALL FRAGMENT SERIES**
- 5** **LARGE FRAGMENT SERIES**
- 11** **SCREW SERIES**
- 14** **INSTRUMENTS SET**

Proximal Humeral Plate I



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12010001S	JM12010001T	4 holes L	87
JM12010002S	JM12010002T	5 holes L	99
JM12010003S	JM12010003T	6 holes L	111
JM12010004S	JM12010004T	7 holes L	123
JM12010005S	JM12010005T	8 holes L	135
JM12010006S	JM12010006T	9 holes L	147
JM12010007S	JM12010007T	10 holes L	159
JM12010008S	JM12010008T	4 holes R	87
JM12010009S	JM12010009T	5 holes R	99
JM12010010S	JM12010010T	6 holes R	111
JM12010011S	JM12010011T	7 holes R	123
JM12010012S	JM12010012T	8 holes R	135
JM12010013S	JM12010013T	9 holes R	147
JM12010014S	JM12010014T	10 holes R	159

Use 4.0mm cancellous screws,3.5mm cortical screws.

Distal Humeral Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12020001S	JM12020001T	3 holes L	88.5
JM12020002S	JM12020002T	4 holes L	104.5
JM12020003S	JM12020003T	5 holes L	120.5
JM12020004S	JM12020004T	6 holes L	136.5
JM12020005S	JM12020005T	3 holes R	88.5
JM12020006S	JM12020006T	4 holes R	104.5
JM12020007S	JM12020007T	5 holes R	120.5
JM12020008S	JM12020008T	6 holes R	136.5

Use 4.0mm cancellous screws,3.5mm cortical screws.

Clavicular Hook Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12030001S	JM12030001T	4 holes L	92.5
JM12030002S	JM12030002T	5 holes L	101.5
JM12030003S	JM12030003T	6 holes L	110
JM12030004S	JM12030004T	7 holes L	119
JM12030005S	JM12030005T	8 holes L	129
JM12030006S	JM12030006T	4 holes R	92.5
JM12030007S	JM12030007T	5 holes R	101.5
JM12030008S	JM12030008T	6 holes R	110
JM12030009S	JM12030009T	7 holes R	119
JM12030010S	JM12030010T	8 holes R	129

Use 4.0mm cancellous screws,3.5mm cortical screws.

Clover Leaf Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12040001S	JM12040001T	3 holes	74
JM12040002S	JM12040002T	4 holes	90
JM12040003S	JM12040003T	5 holes	106
JM12040004S	JM12040004T	6 holes	122
JM12040005S	JM12040005T	7 holes	138
JM12040006S	JM12040006T	8 holes	154

Use 4.0mm cancellous screws,3.5mm cortical screws.

Obliqued T-Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12050001S	JM12050001T	3 holes L	53
JM12050002S	JM12050002T	4 holes L	64
JM12050003S	JM12050003T	5 holes L	75
JM12050004S	JM12050004T	6 holes L	86
JM12050005S	JM12050005T	3 holes R	53
JM12050006S	JM12050006T	4 holes R	64
JM12050007S	JM12050007T	5 holes R	75
JM12050008S	JM12050008T	6 holes R	86

Use 4.0mm cancellous screws,3.5mm cortical screws.

T-Plate II



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12060001S	JM12060001T	3 holes	46
JM12060002S	JM12060002T	4 holes	57
JM12060003S	JM12060003T	5 holes	68
JM12060004S	JM12060004T	6 holes	79

Use 4.0mm cancellous screws,3.5mm cortical screws.


Distal Radius Plate I



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12070001S	JM12070001T	4 holes L	75.7
JM12070002S	JM12070002T	6 holes L	96.7
JM12070003S	JM12070003T	8 holes L	117.7
JM12070004S	JM12070004T	10 holes L	138.7
JM12070005S	JM12070005T	4 holes R	75.7
JM12070006S	JM12070006T	6 holes R	96.7
JM12070007S	JM12070007T	8 holes R	117.7
JM12070008S	JM12070008T	10 holes R	138.7

Use 4.0mm cancellous screws,3.5mm cortical screws.

Distal Radius Plate II



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12080001S	JM12080001T	3 holes L	48.5
JM12080002S	JM12080002T	4 holes L	56.5
JM12080003S	JM12080003T	5 holes L	63.5
JM12080004S	JM12080004T	6 holes L	70.5
JM12080005S	JM12080005T	3 holes R	48.5
JM12080006S	JM12080006T	4 holes R	56.5
JM12080007S	JM12080007T	5 holes R	63.5
JM12080008S	JM12080008T	6 holes R	70.5

Use 4.0mm cancellous screws,3.5mm cortical screws.

Clavicle S-Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12090001S	JM12090001T	6 holes L	94
JM12090002S	JM12090002T	7 holes L	108
JM12090003S	JM12090003T	8 holes L	121
JM12090004S	JM12090004T	6 holes R	94
JM12090005S	JM12090005T	7 holes R	108
JM12090006S	JM12090006T	8 holes R	121

Use 4.0mm cancellous screws,3.5mm cortical screws.

Straight Reconstruction Plate II



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12100001S	JM12100001T	4 holes	48
JM12100002S	JM12100002T	5 holes	60
JM12100003S	JM12100003T	6 holes	72
JM12100004S	JM12100004T	7 holes	84
JM12100005S	JM12100005T	8 holes	96
JM12100006S	JM12100006T	9 holes	108
JM12100007S	JM12100007T	10 holes	120
JM12100008S	JM12100008T	11 holes	132
JM12100009S	JM12100009T	12 holes	144
JM12100010S	JM12100010T	14 holes	168

Use 4.0mm cancellous screws,3.5mm cortical screws.

Straight Reconstruction Plate I



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12110001S	JM12110001T	4 holes	45.5
JM12110002S	JM12110002T	6 holes	57.5
JM12110003S	JM12110003T	7 holes	69.5
JM12110004S	JM12110004T	8 holes	81.5
JM12110005S	JM12110005T	9 holes	93.5
JM12110006S	JM12110006T	10 holes	105.5
JM12110007S	JM12110007T	12 holes	117.5
JM12110008S	JM12110008T	14 holes	141.5
JM12110009S	JM12110009T	16 holes	189.5

Use 4.0mm cancellous screws,3.5mm cortical screws.

Reconstruction S-Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12120001S	JM12120001T	6 holes	68
JM12120002S	JM12120002T	7 holes	79
JM12120003S	JM12120003T	8 holes	89.4
JM12120004S	JM12120004T	9 holes	100
JM12120005S	JM12120005T	10 holes	112


Use 4.0mm cancellous screws,3.5mm cortical screws.

Encircling Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	W/L(mm)
JM12130001S	JM12130001T	4 holes	19/30
JM12130002S	JM12130002T	4 holes	19/35
JM12130003S	JM12130003T	4 holes	19/45
JM12130004S	JM12130004T	4 holes	19/55
JM12130005S	JM12130005T	4 holes	19/60
JM12130006S	JM12130006T	4 holes	19/65
JM12130007S	JM12130007T	4 holes	17/30
JM12130008S	JM12130008T	4 holes	17/35
JM12130009S	JM12130009T	4 holes	17/45
JM12130010S	JM12130010T	4 holes	17/55
JM12130011S	JM12130011T	4 holes	17/60
JM12130012S	JM12130012T	4 holes	17/65
JM12130013S	JM12130013T	4 holes	14/30
JM12130014S	JM12130014T	4 holes	14/35
JM12130015S	JM12130015T	4 holes	14/45
JM12130016S	JM12130016T	4 holes	14/55
JM12130017S	JM12130017T	4 holes	14/60
JM12130018S	JM12130018T	4 holes	14/65

Tubular Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12140001S	JM12140001T	5 holes	61
JM12140002S	JM12140002T	6 holes	73
JM12140003S	JM12140003T	7 holes	85
JM12140004S	JM12140004T	8 holes	97
JM12140005S	JM12140005T	9 holes	109
JM12140006S	JM12140006T	10 holes	121

Use 4.0mm cancellous screws,3.5mm cortical screws.


Ulna & Radius Limited Contact
Dynamic Compression Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12150001S	JM12150001S	4 holes	51
JM12150002S	JM12150002S	5 holes	64
JM12150003S	JM12150003S	6 holes	77
JM12150004S	JM12150004S	7 holes	90
JM12150005S	JM12150005S	8 holes	103
JM12150006S	JM12150006S	9 holes	116
JM12150007S	JM12150007S	10 holes	129
JM12150008S	JM12150008S	11 holes	142
JM12150009S	JM12150009S	12 holes	155

Use 4.0mm cancellous screws,3.5mm cortical screws.

Humeral Limited Contact Dynamic
Compression Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12160001S	JM12160001T	4 holes	56
JM12160002S	JM12160002T	5 holes	70
JM12160003S	JM12160003T	6 holes	84
JM12160004S	JM12160004T	7 holes	98
JM12160005S	JM12160005T	8 holes	112
JM12160006S	JM12160006T	9 holes	126
JM12160007S	JM12160007T	10 holes	140
JM12160008S	JM12160008T	12 holes	168
Use 4.0mm cancellous screws,3.5mm cortical screws.			

Proximal Femoral Plate I



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12170001S	JM12170001T	5 holes L	116.7
JM12170002S	JM12170002T	7 holes L	150.7
JM12170003S	JM12170003T	9 holes L	184.7
JM12170004S	JM12170004T	11 holes L	218.7
JM12170005S	JM12170005T	13 holes L	252.7
JM12170006S	JM12170006T	5 holes R	116.7
JM12170007S	JM12170007T	7 holes R	150.7
JM12170008S	JM12170008T	9 holes R	184.7
JM12170009S	JM12170009T	11 holes R	218.7
JM12170010S	JM12170010T	13 holes R	252.7
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Distal Lateral Femoral Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12180001S	JM12180001T	5 holes L	134
JM12180002S	JM12180002T	7 holes L	166
JM12180003S	JM12180003T	9 holes L	198
JM12180004S	JM12180004T	11 holes L	230
JM12180005S	JM12180005T	13 holes L	262
JM12180006S	JM12180006T	5 holes R	134
JM12180007S	JM12180007T	7 holes R	166
JM12180008S	JM12180008T	9 holes R	198
JM12180009S	JM12180009T	11 holes R	230
JM12180010S	JM12180010T	13 holes R	262
Use 6.5mm cancellous screws,3.5mm cortical screws.			

Distal Medial Femoral Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12190001S	JM12190001S	4 holes L	107.5
JM12190002S	JM12190002S	6 holes L	139.5
JM12190003S	JM12190003S	8 holes L	171.5
JM12190004S	JM12190004S	10 holes L	203.5
JM12190005S	JM12190005S	4 holes R	107.5
JM12190006S	JM12190006S	6 holes R	139.5
JM12190007S	JM12190007S	8 holes R	171.5
JM12190008S	JM12190008S	10 holes R	203.5
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Femur Condyle Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12200001S	JM12200001T	5 holes L	122.3
JM12200002S	JM12200002T	7 holes L	156.3
JM12200003S	JM12200003T	9 holes L	190.3
JM12200004S	JM12200004T	11 holes L	224.3
JM12200005S	JM12200005T	13 holes L	258.3
JM12200006S	JM12200006T	5 holes R	122.3
JM12200007S	JM12200007T	7 holes R	156.3
JM12200008S	JM12200008T	9 holes R	190.3
JM12200009S	JM12200009T	11 holes R	224.3
JM12200010S	JM12200010T	13 holes R	258.3
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Proximal Dorsal Tibia Plate I




REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12210001S	JM12210001T	5 holes L	128
JM12210002S	JM12210002T	7 holes L	160
JM12210003S	JM12210003T	9 holes L	192
JM12210004S	JM12210004T	11 holes L	224
JM12210005S	JM12210005T	13 holes L	256
JM12210006S	JM12210006T	5 holes R	128
JM12210007S	JM12210007T	7 holes R	160
JM12210008S	JM12210008T	9 holes R	192
JM12210009S	JM12210009T	11 holes R	224
JM12210010S	JM12210010T	13 holes R	256
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Proximal Medial Tibia Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12220001S	JM12220001T	5 holes L	125
JM12220002S	JM12220002T	7 holes L	157
JM12220003S	JM12220003T	9 holes L	189
JM12220004S	JM12220004T	11 holes L	221
JM12220005S	JM12220005T	13 holes L	253
JM12220006S	JM12220006T	5 holes R	125
JM12220007S	JM12220007T	7 holes R	157
JM12220008S	JM12220008T	9 holes R	189
JM12220009S	JM12220009T	11 holes R	221
JM12220010S	JM12220010T	13 holes R	253
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Distal Dorsal Tibia Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12230001S	JM12230001T	5 holes L	108
JM12230002S	JM12230002T	7 holes L	140
JM12230003S	JM12230003T	9 holes L	172
JM12230004S	JM12230004T	11 holes L	204
JM12230005S	JM12230005T	13 holes L	236
JM12230006S	JM12230006T	5 holes R	108
JM12230007S	JM12230007T	7 holes R	140
JM12230008S	JM12230008T	9 holes R	172
JM12230009S	JM12230009T	11 holes R	204
JM12230010S	JM12230010T	13 holes R	236
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Distal Medial Tibia Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12240001S	JM12240001T	5 holes L	108
JM12240002S	JM12240002T	7 holes L	140
JM12240003S	JM12240003T	9 holes L	172
JM12240004S	JM12240004T	11 holes L	204
JM12240005S	JM12240005T	13 holes L	236
JM12240006S	JM12240006T	5 holes R	108
JM12240007S	JM12240007T	7 holes R	140
JM12240008S	JM12240008T	9 holes R	172
JM12240009S	JM12240009T	11 holes R	204
JM12240010S	JM12240010T	13 holes R	236
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Distal Lateral Tibia L-Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12250001S	JM12250001T	6 holes L	97
JM12250002S	JM12250002T	8 holes L	121
JM12250003S	JM12250003T	10 holes L	145
JM12250004S	JM12250004T	12 holes L	169
JM12250005S	JM12250005T	14 holes L	193
JM12250006S	JM12250006T	16 holes L	217
JM12250007S	JM12250007T	18 holes L	241
JM12250008S	JM12250008T	20 holes L	265
JM12250009S	JM12250009T	22 holes L	289
JM12250010S	JM12250010T	6 holes R	97
JM12250011S	JM12250011T	8 holes R	121
JM12250012S	JM12250012T	10 holes R	145
JM12250013S	JM12250013T	12 holes R	169
JM12250014S	JM12250014T	14 holes R	193
JM12250015S	JM12250015T	16 holes R	217
JM12250016S	JM12250016T	18 holes R	241
JM12250017S	JM12250017T	20 holes R	265
JM12250018S	JM12250018T	22 holes R	289
Use 4.0mm cancellous screws,3.5mm cortical screws.			

L-Plate III



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12260001S	JM12260001T	4 holes L	76.5
JM12260002S	JM12260002T	5 holes L	92.5
JM12260003S	JM12260003T	6 holes L	108.5
JM12260004S	JM12260004T	8 holes L	140.5
JM12260005S	JM12260005T	10 holes L	172.5
JM12260006S	JM12260006T	4 holes R	76.5
JM12260007S	JM12260007T	5 holes R	92.5
JM12260008S	JM12260008T	6 holes R	108.5
JM12260009S	JM12260009T	8 holes R	140.5
JM12260010S	JM12260010T	10 holes R	172.5
Use 6.5mm cancellous screws,4.5mm cortical screws.			

T-Plate IV



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12270001S	JM12270001T	4 holes L	81
JM12270002S	JM12270002T	6 holes L	113
JM12270003S	JM12270003T	8 holes L	145
JM12270004S	JM12270004T	4 holes R	81
JM12270005S	JM12270005T	6 holes R	113
JM12270006S	JM12270006T	8 holes R	145
Use 6.5mm cancellous screws,4.5mm cortical screws.			

Patella Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12280001S	JM12280001T	S	47
JM12280002S	JM12280002T	M	52
JM12280003S	JM12280003T	L	57.5
Use 4.0mm cancellous screws,3.5mm cortical screws.			

Distal Fibular Plate III



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12290001S	JM12290001T	2 holes L	53
JM12290002S	JM12290002T	3 holes L	63.5
JM12290003S	JM12290003T	4 holes L	74
JM12290004S	JM12290004T	5 holes L	84.5
JM12290005S	JM12290005T	6 holes L	95
JM12290006S	JM12290006T	7 holes L	105.5
JM12290007S	JM12290007T	8 holes L	116
JM12290008S	JM12290008T	9 holes L	126.5
JM12290009S	JM12290009T	10 holes L	137
JM12290010S	JM12290010T	11 holes L	147.5
JM12290011S	JM12290011T	12 holes L	158
JM12290012S	JM12290012T	2 holes R	53
JM12290013S	JM12290013T	3 holes R	63.5
JM12290014S	JM12290014T	4 holes R	74
JM12290015S	JM12290015T	5 holes R	84.5
JM12290016S	JM12290016T	6 holes R	95
JM12290017S	JM12290017T	7 holes R	105.5
JM12290018S	JM12290018T	8 holes R	116
JM12290019S	JM12290019T	9 holes R	126.5
JM12290020S	JM12290020T	10 holes R	137
JM12290021S	JM12290021T	11 holes R	147.5
JM12290022S	JM12290022T	12 holes R	158

Use 4.0mm cancellous screws,3.5mm cortical screws.

Calcaneal Plate I



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12300001S	JM12300001T	13 holes L	80
JM12300002S	JM12300002T	13 holes R	80
JM12300003S	JM12300003T	15 holes L	90
JM12300004S	JM12300004T	15 holes R	90

Use 4.0mm cancellous screws,3.5mm cortical screws.

Calcaneal Plate II



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12310001S	JM12310001T	8 holes L	44
JM12310002S	JM12310002T	8 holes R	44

Use 4.0mm cancellous screws,3.5mm cortical screws.

Calcaneal Plate III



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12320001S	JM12320001T	11 holes L	60
JM12320002S	JM12320002T	13 holes L	70
JM12320003S	JM12320003T	11 holes R	60
JM12320004S	JM12320004T	13 holes R	70

Use 4.0mm cancellous screws,3.5mm cortical screws.

Calcaneal Plate VII



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12330001T	JM12330001T	12 holes L	59
JM12330002T	JM12330002T	12 holes R	59

Use 4.0mm cancellous screws,3.5mm cortical screws.

Curved Reconstruction Plate I



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12340001S	JM12340001T	4 holes	48.2
JM12340002S	JM12340002T	5 holes	59.5
JM12340003S	JM12340003T	6 holes	70.5
JM12340004S	JM12340004T	7 holes	80.9
JM12340005S	JM12340005T	8 holes	90.7
JM12340006S	JM12340006T	9 holes	99.8
JM12340007S	JM12340007T	10 holes	108.2
JM12340008S	JM12340008T	12 holes	122.3
JM12340009S	JM12340009T	14 holes	132.7

Use 4.0mm cancellous screws,3.5mm cortical screws.

Tibia Limited Contact
Dynamic Compression Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12350001S	JM12350001T	6 holes	106
JM12350002S	JM12350002T	7 holes	124
JM12350003S	JM12350003T	8 holes	142
JM12350004S	JM12350004T	9 holes	160
JM12350005S	JM12350005T	10 holes	178
JM12350006S	JM12350006T	12 holes	214
JM12350007S	JM12350007T	14 holes	250
JM12350008S	JM12350008T	16 holes	286

Use 4.0mm cancellous screws,3.5mm cortical screws.

Femoral Limited Contact
Dynamic Compression Plate



REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12360001S	JM12360001T	6 holes	106
JM12360002S	JM12360002T	7 holes	124
JM12360003S	JM12360003T	8 holes	142
JM12360004S	JM12360004T	9 holes	160
JM12360005S	JM12360005T	10 holes	178
JM12360006S	JM12360006T	12 holes	214
JM12360007S	JM12360007T	14 holes	250
JM12360008S	JM12360008T	16 holes	286

Use 4.0mm cancellous screws,3.5mm cortical screws.

3.5mm Cortical Screw,Fully Threaded



REF (Stainless Steel)	REF (Titanium)	Specifications
JM12373512S	JM12373512T	Φ3.5×12mm
JM12373514S	JM12373514T	Φ3.5×14mm
JM12373516S	JM12373516T	Φ3.5×16mm
JM12373518S	JM12373518T	Φ3.5×18mm
JM12373520S	JM12373520T	Φ3.5×20mm
JM12373522S	JM12373522T	Φ3.5×22mm
JM12373524S	JM12373524T	Φ3.5×24mm
JM12373526S	JM12373526T	Φ3.5×26mm
JM12373528S	JM12373528T	Φ3.5×28mm
JM12373530S	JM12373530T	Φ3.5×30mm
JM12373532S	JM12373532T	Φ3.5×32mm
JM12373534S	JM12373534T	Φ3.5×34mm
JM12373536S	JM12373536T	Φ3.5×36mm
JM12373538S	JM12373538T	Φ3.5×38mm
JM12373540S	JM12373540T	Φ3.5×40mm
JM12373542S	JM12373542T	Φ3.5×42mm
JM12373544S	JM12373544T	Φ3.5×44mm
JM12373546S	JM12373546T	Φ3.5×46mm
JM12373548S	JM12373548T	Φ3.5×48mm
JM12373550S	JM12373550T	Φ3.5×50mm
JM12373555S	JM12373555T	Φ3.5×55mm
JM12373560S	JM12373560T	Φ3.5×60mm

4.5mm Cortical Screw,Fully Threaded



REF (Stainless Steel)	REF (Titanium)	Specifications
JM12384512S	JM12384512T	Φ4.5×12mm
JM12384514S	JM12384514T	Φ4.5×14mm
JM12384516S	JM12384516T	Φ4.5×16mm
JM12384518S	JM12384518T	Φ4.5×18mm
JM12384520S	JM12384520T	Φ4.5×20mm
JM12384522S	JM12384522T	Φ4.5×22mm
JM12384524S	JM12384524T	Φ4.5×24mm
JM12384526S	JM12384526T	Φ4.5×26mm
JM12384528S	JM12384528T	Φ4.5×28mm
JM12384530S	JM12384530T	Φ4.5×30mm
JM12384532S	JM12384532T	Φ4.5×32mm
JM12384534S	JM12384534T	Φ4.5×34mm
JM12384536S	JM12384536T	Φ4.5×36mm
JM12384538S	JM12384538T	Φ4.5×38mm
JM12384540S	JM12384540T	Φ4.5×40mm
JM12384542S	JM12384542T	Φ4.5×42mm
JM12384544S	JM12384544T	Φ4.5×44mm
JM12384546S	JM12384546T	Φ4.5×46mm
JM12384548S	JM12384548T	Φ4.5×48mm
JM12384550S	JM12384550T	Φ4.5×50mm
JM12384552S	JM12384552T	Φ4.5×52mm
JM12384554S	JM12384554T	Φ4.5×54mm
JM12384556S	JM12384556T	Φ4.5×56mm
JM12384558S	JM12384558T	Φ4.5×58mm
JM12384560S	JM12384560T	Φ4.5×60mm

4.0mm Cancellous Screw,
Partially Threaded



REF (Stainless Steel)	REF (Titanium)	Specifications
JM12390412S	JM12390412T	Φ4.0×12mm
JM12390414S	JM12390414T	Φ4.0×14mm
JM12390416S	JM12390416T	Φ4.0×16mm
JM12390418S	JM12390418T	Φ4.0×18mm
JM12390420S	JM12390420T	Φ4.0×20mm
JM12390422S	JM12390422T	Φ4.0×22mm
JM12390424S	JM12390424T	Φ4.0×24mm
JM12390426S	JM12390426T	Φ4.0×26mm
JM12390428S	JM12390428T	Φ4.0×28mm
JM12390430S	JM12390430T	Φ4.0×30mm
JM12390432S	JM12390432T	Φ4.0×32mm
JM12390434S	JM12390434T	Φ4.0×34mm
JM12390436S	JM12390436T	Φ4.0×36mm
JM12390438S	JM12390438T	Φ4.0×38mm
JM12390440S	JM12390440T	Φ4.0×40mm
JM12390442S	JM12390442T	Φ4.0×42mm
JM12390444S	JM12390444T	Φ4.0×44mm
JM12390446S	JM12390446T	Φ4.0×46mm
JM12390448S	JM12390448T	Φ4.0×48mm
JM12390450S	JM12390450T	Φ4.0×50mm

4.0mm Cancellous Screw, Fully Threaded



REF (Stainless Steel)	REF (Titanium)	Specifications
JM12400414S	JM12400414T	Φ4.0×14mm
JM12400416S	JM12400416T	Φ4.0×16mm
JM12400418S	JM12400418T	Φ4.0×18mm
JM12400420S	JM12400420T	Φ4.0×20mm
JM12400422S	JM12400422T	Φ4.0×22mm
JM12400424S	JM12400424T	Φ4.0×24mm
JM12400426S	JM12400426T	Φ4.0×26mm
JM12400428S	JM12400428T	Φ4.0×28mm
JM12400430S	JM12400430T	Φ4.0×30mm
JM12400432S	JM12400432T	Φ4.0×32mm
JM12400434S	JM12400434T	Φ4.0×34mm
JM12400436S	JM12400436T	Φ4.0×36mm
JM12400438S	JM12400438T	Φ4.0×38mm
JM12400440S	JM12400440T	Φ4.0×40mm
JM12400442S	JM12400442T	Φ4.0×42mm
JM12400444S	JM12400444T	Φ4.0×44mm
JM12400446S	JM12400446T	Φ4.0×46mm
JM12400448S	JM12400448T	Φ4.0×48mm
JM12400450S	JM12400450T	Φ4.0×50mm
JM12400452S	JM12400452T	Φ4.0×52mm
JM12400454S	JM12400454T	Φ4.0×54mm
JM12400456S	JM12400456T	Φ4.0×56mm
JM12400458S	JM12400458T	Φ4.0×58mm
JM12400460S	JM12400460T	Φ4.0×60mm

6.5mm Cancellous Screw, Partially Threaded




REF (Stainless Steel)	REF (Titanium)	Specifications (mm)	Length (mm)
JM12416535S	JM12416535T	Φ6.5×35mm	(Thread L=16mm)
JM12416540S	JM12416540T	Φ6.5×40mm	(Thread L=16mm)
JM12416545S	JM12416545T	Φ6.5×45mm	(Thread L=16mm)
JM12416550S	JM12416550T	Φ6.5×50mm	(Thread L=16mm)
JM12416555S	JM12416555T	Φ6.5×55mm	(Thread L=16mm)
JM12416560S	JM12416560T	Φ6.5×60mm	(Thread L=16mm)
JM12416565S	JM12416565T	Φ6.5×65mm	(Thread L=16mm)
JM12416570S	JM12416570T	Φ6.5×70mm	(Thread L=16mm)
JM12416575S	JM12416575T	Φ6.5×75mm	(Thread L=16mm)
JM12416580S	JM12416580T	Φ6.5×80mm	(Thread L=16mm)
JM12416585S	JM12416585T	Φ6.5×85mm	(Thread L=16mm)
JM12416590S	JM12416590T	Φ6.5×90mm	(Thread L=16mm)
JM12416595S	JM12416595T	Φ6.5×95mm	(Thread L=16mm)
JM12416500S	JM12416500T	Φ6.5×100mm	(Thread L=16mm)
JM12416505S	JM12416505T	Φ6.5×105mm	(Thread L=16mm)
JM12416510S	JM12416510T	Φ6.5×110mm	(Thread L=16mm)
JM12416515S	JM12416515T	Φ6.5×115mm	(Thread L=16mm)
JM12416520S	JM12416520T	Φ6.5×120mm	(Thread L=16mm)
JM12426545S	JM12426545T	Φ6.5×45mm	(Thread L=32mm)
JM12426550S	JM12426550T	Φ6.5×50mm	(Thread L=32mm)
JM12426555S	JM12426555T	Φ6.5×55mm	(Thread L=32mm)
JM12426560S	JM12426560T	Φ6.5×60mm	(Thread L=32mm)
JM12426565S	JM12426565T	Φ6.5×65mm	(Thread L=32mm)
JM12426570S	JM12426570T	Φ6.5×70mm	(Thread L=32mm)
JM12426575S	JM12426575T	Φ6.5×75mm	(Thread L=32mm)
JM12426580S	JM12426580T	Φ6.5×80mm	(Thread L=32mm)
JM12426585S	JM12426585T	Φ6.5×85mm	(Thread L=32mm)
JM12426590S	JM12426590T	Φ6.5×90mm	(Thread L=32mm)
JM12426595S	JM12426595T	Φ6.5×95mm	(Thread L=32mm)
JM12426500S	JM12426500T	Φ6.5×100mm	(Thread L=32mm)
JM12426505S	JM12426505T	Φ6.5×105mm	(Thread L=32mm)
JM12426510S	JM12426510T	Φ6.5×110mm	(Thread L=32mm)
JM12426515S	JM12426515T	Φ6.5×115mm	(Thread L=32mm)
JM12426520S	JM12426520T	Φ6.5×120mm	(Thread L=32mm)

6.5mm Cancellous Screw, Fully Threaded



REF (Stainless Steel)	REF (Titanium)	Specifications
JM12416535S	JM12416535T	Φ6.5×35mm
JM12416540S	JM12416540T	Φ6.5×40mm
JM12416545S	JM12416545T	Φ6.5×45mm
JM12416550S	JM12416550T	Φ6.5×50mm
JM12416555S	JM12416555T	Φ6.5×55mm
JM12416560S	JM12416560T	Φ6.5×60mm
JM12416565S	JM12416565T	Φ6.5×65mm
JM12416570S	JM12416570T	Φ6.5×70mm
JM12416575S	JM12416575T	Φ6.5×75mm
JM12416580S	JM12416580T	Φ6.5×80mm
JM12416585S	JM12416585T	Φ6.5×85mm
JM12416590S	JM12416590T	Φ6.5×90mm
JM12416595S	JM12416595T	Φ6.5×95mm
JM12416500S	JM12416500T	Φ6.5×100mm
JM12416505S	JM12416505T	Φ6.5×105mm
JM12416510S	JM12416510T	Φ6.5×110mm
JM12416515S	JM12416515T	Φ6.5×115mm
JM12416520S	JM12416520T	Φ6.5×120mm

Small Fragment Instrument Set



Product NO	Products	QTY.
JM71300100	Guide Pin,Φ2.0×150mm	4
JM71300200	Drill Bit,Φ2.0mm	2
JM71300300	Drill Bit,Φ2.5mm	2
JM71300400	Drill Bit,Φ3.5mm	2
JM71300500	Tap,HA3.5mm	1
JM71300600	Tap,HB4.0mm	1
JM71300700	Drill Sleeve,Φ2.5mm	1
JM71300800	Double Drill Guide,Φ2.0 - 4.0mm	1
JM71300900	Double Drill Guide,Φ2.5 - 3.5mm	1
JM71301000	LCP Drill Sleeve,Φ2.5mm	1
JM71301100	Screw Extraction,SW2.5	1
JM71301200	Countersink Drill,Φ6.5mm	1
JM71301300	Hollow Reamer,Φ4.0mm	1
JM71301400	Depth Gauge,0 - 60mm	1
JM71301500	Hex Screwdriver With Q,C,SW2.5	1
JM71301600	Screw Holding Forceps	1
JM71301700	Hex Screwdriver,SW2.5	1
JM71301800	T Handle With Quick Coupling	1
JM71301900	Plate Bender,A	1
JM71302000	Plate Bender,B	1
JM71302100	Curved Reduction Forceps	1
JM71302200	Self - centering Bone Holding Forceps	2
JM71302300	Reduction Forceps With Points	1
JM71302400	Sharp Hook	1
JM71302500	Retractor ,15.5	1
JM71302600	Retractor ,8	1
JM71302700	Periosteum Elevator,6.5mm	1
JM71302800	Template,Long	1
JM71302900	Template,Short	1
JM71303000	Double Drill Guide,Φ2.5 - 3.5mm	1
JM71303100	Screw Case	1
JM71303200	Instrument Case	1

1 Guide Pin,Φ2.0×150mm



Product NO JM71300100

2 Drill Bit,Φ2.0mm



Product NO JM71300200

3 Drill Bit,Φ2.5mm



Product NO JM71300300

4 Drill Bit,Φ3.5mm



Product NO JM71300400

5 Tap,HA3.5mm



Product NO JM71300500

6 Tap,HB4.0mm



Product NO JM71300600

7 Drill Sleeve,Φ2.5mm



Product NO JM71300700

8 Double Drill Guide,
Φ2.0 - 4.0mm



Product NO JM71300800

9 Double Drill Guide,
Φ2.5 - 3.5mm



Product NO JM71300900

10 LCP Drill Sleeve,Φ2.5mm



Product NO JM71301000

11 Screw Extraction,SW2.5



Product NO JM71301100

12 Countersink Drill,Φ6.5mm



Product NO JM71301200

13 Hollow Reamer,Φ4.0mm



Product NO JM71301300

14 Depth Gauge,0 - 60mm



Product NO JM71301400

15 Hex Screwdriver With Q.C,
SW2.5



Product NO JM71301500

16 Screw Holding Forceps



Product NO JM71301600

17 Hex Screwdriver,SW2.5



Product NO JM71301700

18 T Handle With Quick Coupling



Product NO JM71301800

19 Plate Bender,A



Product NO JM71301900

20 Plate Bender,B



Product NO JM71302000

21 Curved Reduction Forceps



Product NO JM71302100

22 Self-centering Bone Holding Forceps



Product NO JM71302200

23 Reduction Forceps With Points




Product NO JM71302300

24 Sharp Hook



Product NO JM71302400

31 Screw Case



Product NO JM71303100

25 Retractor ,15.5



Product NO JM71302500

26 Retractor ,8




Product NO JM71302600

27 Periosteum Elevator,6.5mm



Product NO JM71302700

28 Template,Long



Product NO JM71302800

29 Template,Short



Product NO JM71302900

30 Double Drill Guide, Φ2.5 - 3.5mm



Product NO JM71303000

Encircling Plate Instrument Set



Product NO	Products	QTY.
JM71310100	Compressing Clamp	1
JM71310200	One - Side Compressing Clamp	1
JM71310300	Curved Compressing Clamp	1
JM71310400	Cutter	1
JM71310500	Instrument Case	1

1 Compressing Clamp



Product NO	JM71310100
------------	------------

2 One - Side Compressing Clamp



Product NO	JM71310200
------------	------------

3 Curved Compressing Clamp



Product NO	JM71310300
------------	------------

4 Cutter



Product NO	JM71310400
------------	------------

Large Fragment Instrument Set



Product NO	Products	QTY.
JM71301800	T Handle With Quick Coupling	1
JM71302400	Sharp Hook	1
JM71320300	Periosteum Elevator, 15mm	1
JM71320400	Double Drill Guide, Φ3.2 - 4.5mm	1
JM71320500	Retractor , Small	1
JM71320600	Retractor , Large	1
JM71320700	Guide Pin, Φ2.0×250mm	2
JM71320800	Guide Pin, Φ2.0×300mm	2
JM71320900	Periosteum Eevator, 8mm	1
JM71321000	Template, Short	1
JM71321100	Template, Middle	1
JM71321200	Template, Long	1
JM71321300	Drill Bit, Φ3.2mm	2
JM71321400	Drill Bit, Φ4.5mm	2
JM71321500	Drill Sleeve, Φ3.2mm	1
JM71321600	Tap, HA4.5mm	1
JM71321700	Tap, HB6.5mm	1
JM71321800	Double Drill Guide, Φ3.2 - 4.5mm	1
JM71321900	Double Drill Guide, Φ3.2 - 6.5mm	1
JM71322000	LCP - DCP Guide, Φ3.2mm	1
JM71322100	Countersink Drill, Φ8.8mm	1
JM71322200	Depth Gauge, 0 - 100mm	1
JM71322300	Plate Bender, A	1
JM71322400	Plate Bender, B	1
JM71322500	Screw Holding Forceps	1
JM71322600	Hex Screwdriver With Q.C.SW3.5	1
JM71322700	Screw Extraction, SW3.5	1
JM71322800	Hex Screwdriver, SW3.5	1
JM71322900	Hollow Reamer, Φ5.0mm	1
JM71323000	Self - centering Bone Holding Forceps	2
JM71323100	Curved Reduction Forceps	1
JM71323200	Reduction Forceps, with Points	1
JM71323300	Screw Case	1
JM71323400	Instrument Case	1

1 T Handle With Quick Coupling



Product NO	JM71301800
------------	------------

2 Sharp Hook



Product NO	JM71302400
------------	------------

3 Periosteum Elevator, 15mm



Product NO	JM71320300
------------	------------

4 Double Drill Guide,
Φ3.2 - 4.5mm




Product NO JM71320400

5 Retractor ,Small



Product NO JM71320500

6 Retractor , Large



Product NO JM71320600

7 Guide Pin, Φ2.0×250mm



Product NO JM71320700

8 Guide Pin, Φ2.0×300mm




Product NO JM71320800

9 Periosteum Eevator, 8mm



Product NO JM71320900

10 Template, Short



Product NO JM71321000

11 Template, Middle



Product NO JM71321100

12 Template, Long



Product NO JM71321200

13 Drill Bit, Φ3.2mm



Product NO JM71321300

14 Drill Bit, Φ4.5mm



Product NO JM71321400

15 Drill Sleeve, Φ3.2mm



Product NO JM71321500

16 Tap, HA4.5mm



Product NO JM71321600

17 Tap, HB6.5mm



Product NO JM71321700

18 Double Drill Guide,
Φ3.2 - 4.5mm



Product NO JM71321800

19 Double Drill Guide,
Φ3.2 - 6.5mm



Product NO JM71321900

20 LCP - DCP Guide, Φ3.2mm



Product NO JM71322000

21 Countersink Drill, Φ8.8mm



Product NO JM71322100

22 Depth Gauge, 0 - 100mm



Product NO JM71322200

23 Plate Bender, A



Product NO JM71322300

24 Plate Bender, B



Product NO JM71322400

25 Screw Holding Forceps



Product NO JM71322500

26 Hex Screwdriver With Q.C.SW3.5



Product NO JM71322600

27 Screw Extraction, SW3.5



Product NO JM71322700

28 Hex Screwdriver, SW3.5



Product NO JM71322800

29 Hollow Reamer, Ø5.0mm



Product NO JM71322900

30 Self - centering Bone Holding Forceps




Product NO JM71323000

31 Curved Reduction Forceps




Product NO JM71323100

32 Reduction Forceps, with Points



Product NO JM71323200

33 Screw Case



Product NO JM71323300

Reconstruction Plate Instrument Set



Product NO	Products	QTY.
JM71300500	Tap, HA3.5mm	1
JM71300600	Tap, HB4.0mm	1
JM71300800	Double Drill Guide, Φ2.0 - 4.0mm	1
JM71300900	Double Drill Guide, Φ2.5 - 3.5mm	1
JM71301400	Depth Gauge, 0 - 60mm	1
JM71301500	Hex Screwdriver, SW2.5	1
JM71330700	Ball Tips, Curved	1
JM71330800	Ball Tips, Flat	1
JM71330900	Ball Spike	1
JM71331000	Side - bending Reduction Forceps, Long	1
JM71331100	Side - bending Reduction Forceps, Short	1
JM71331200	Reduction Forceps, Straight	1
JM71331300	Small Bone Holding Forceps	1
JM71331400	Large Bone Holding Forceps	1
JM71331500	Double Claw Reduction Forceps, Different Claw	1
JM71331600	Double Claw Reduction Forceps	1
JM71331700	Acetabular Reduction Forceps	1
JM71331800	Plate Bending Forceps	1
JM71331900	Hook	1
JM71332000	T Hook	1
JM71332100	Sciatic Nerve Retractor	1
JM71332200	Universal Hex Screwdriver, SW2.5	1
JM71332300	Flexible Drill Bit, Φ2.5mm	1
JM71332400	Drill Guide, Φ2.5mm	1
JM71332500	Template, Straight	1
JM71332600	Template, Curved	1
JM71332700	Three Claw Reduction Forceps	1
JM71332800	Plate Bender, A	1
JM71332900	Plate Bender, B	1
JM71333000	Sleeve for Hex Screwdriver	1
JM71333100	Drill Bit, Φ2.5×125mm	1
JM71333200	Screw Holding Forceps	1
JM71333300	Instrument Case	1

1 Tap, HA3.5mm




Product NO JM71300500

2 Tap, HB4.0mm



Product NO JM71300600

3 Double Drill Guide, Φ2.0 - 4.0mm



Product NO JM71300800

4 Double Drill Guide, Φ2.5 - 3.5mm



Product NO JM71300900

5 Depth Gauge, 0 - 60mm




Product NO JM71301400

6 Hex Screwdriver, SW2.5




Product NO JM71301500

7 Ball Tips, Curved



Product NO JM71330700

8 Ball Tips, Flat



Product NO JM71330800

9 Ball Spike




Product NO JM71330900

10 Side - bending Reduction Forceps, Long



Product NO JM71331000

11 Side - bending Reduction Forceps, Short



Product NO JM71331100

12 Reduction Forceps, Straight



Product NO JM71331200

13 Small Bone Holding Forceps




Product NO JM71331300

14 Large Bone Holding Forceps



Product NO JM71331400

15 Double Claw Reduction Forceps, Different Claw



Product NO JM71331500

16 Double Claw Reduction Forceps



Product NO JM71331600

17 Acetabular Reduction Forceps



Product NO JM71331700

18 Plate Bending Forceps




Product NO JM71331800

19 Hook



Product NO JM71331900

20 T Hook



Product NO JM71332000

21 Sciatic Nerve Retractor



Product NO JM71332100

22 Universal Hex Screwdriver, SW2.5



Product NO JM71332200

23 Flexible Drill Bit, Ø2.5mm



Product NO JM71332300

24 Drill Guide, Ø2.5mm



Product NO JM71332400

25 Template, Straight



Product NO JM71332500

26 Template, Curved



Product NO JM71332600

27 Three Claw Reduction Forceps




Product NO JM71332700

28 Plate Bender, A



Product NO JM71332800

29 Plate Bender, B



Product NO JM71332900

30 Sleeve for Hex Screwdriver



Product NO JM71333000

31 Drill Bit, Ø2.5×125mm



Product NO	JM71333100
------------	------------

32 Screw Holding Forceps



Product NO	JM71333200
------------	------------