



**Jin Ming Medical Device Co., Ltd.**

Address: Room 601, 602, 6th Floor, Building 4, No. 12 Zhicheng Road, Nexus Valley, Songshan Lake Park, Dongguan City, Guangdong Province, China  
Tel: +86 769 22803823  
Email: [info@jinmingmed.com](mailto:info@jinmingmed.com)  
Website: [www.jinmingmed.com](http://www.jinmingmed.com)  
Copyright © 2025 Jin Ming Medical Device Co., Ltd. All rights reserved.

# POWER TOOLS

## PRODUCT CATALOG

N series multifunctional power system



**JM81011000**



**Multifunctional host**  
Speed/Frequency 0-18000 Times/min

N series multifunctional power system



**Swing saw**  
Swing saw assembly  
Sawing frequency 0-16700 Times/min



**Reciprocating saw**  
Reciprocating saw assembly  
Sawing frequency 0-11000Times/min  
Reciprocating stroke 4.2mm



**Sternum saw**  
Sternum saw assemb  
Sawing frequency 0-11000Times/min  
Reciprocating stroke 4.2mm



**Craniotomy drill**  
Craniotomy drill assembly  
Speed 0-800Times/min  
Clamping range 6-12mm



**Slow drilling**  
Slow drill assembly  
Speed 0-290 Times/min  
Clamping range 0.8-8mm



**Hollow drill**  
Hollow drill assembly  
Speed 0-800 Times/min  
Clamping range 0.8-4.2mm



**Craniotomy milling**  
Craniotomy milling assembly  
Speed 0-37000Times/min  
Milling cutter 3mm



**Kirschner wire**  
Kirschner wire assembly  
Speed 0-1200Times/min  
0.8-1.8mm  
Clamping range 1.5-2.5mm  
2.0-3.2mm

D series multifunction



Multiple functions in one, wide application range. Quick replacement of components, convenient and fast, suitable for closing section, trauma, neurosurgery, intracardiac surgery, fitness. The components used include hollow drill assembly and Kirschner wire set parts, bone drill assembly, slow drill assembly, swing saw assembly, craniotomy milling assembly, reciprocating saw assembly, sternum saw group 1, craniotomy drill assembly.

D series multifunction



Swing saw

Swing saw assembly

Sawing frequency 0-16700 Times/min



Reciprocating saw

Reciprocating saw assembly

Sawing frequency 0-11000 Times/min

Reciprocating stroke 4.2mm



Sternum saw

Sternum saw assembly

Sawing frequency 0-11000 Times/min

Reciprocating stroke 4.2mm



Craniotomy drill

Craniotomy drill assembly

Speed 0-800 Times/min

Clamping range 6-12mm



Slow drilling

Slow drill assembly

Speed 0-290 Times/min

Clamping range 0.8-8mm



Hollow drill

Hollow drill assembly

Speed 0-800 Times/min

Clamping range 0.8-4.2mm



Craniotomy milling

Craniotomy milling assembly

Speed 0-37000 Times/min

Milling cutter 3mm



Kirschner wire

Kirschner wire assembly

Speed 0-1200 Times/min  
0.8-1.8mm

Clamping range 1.5-2.5mm  
2.0-3.2mm

S series power

S series power



JM81010200

Speed	0-1200Times/min
Range	0.8-8mm



JM81010300

Speed	0-800 Times/min
Range	0.8-8mm



JM81010400

Speed	0-290Times/min
Interface	stryker



JM81010500

Sawing frequency	16700 Times/min
Interface	AO

S series power

S series power



JM81010600

Sawing frequency	11000 Times/min
Reciprocating stroke	4.2mm



JM81010700

Speed	0-800 Times/min
-------	-----------------



JM81010800

Sawing frequency	11000 Times/min
Reciprocating stroke	4.2mm



JM81010900

Speed	0-37000 Times/min
Milling cutter	3mm

G series power

G series power

JM81011200

Speed	0-1200 Times/min
Range	0.8-8mm



JM81011300

Speed	0-800 Times/min
Range	0.8-8mm



JM81011400

Speed	0-290 Times/min
Range	0.8-8mm
Interface	stryker



JM81011500

Sawing frequency	16700 Times/min
Interface	AO



G series power



JM81011600

Sawing frequency	11000 Times/min
Reciprocating stroke	4.2mm

G series power



JM81011800

Speed	0-800 Times/min
-------	-----------------



JM81011700

Sawing frequency	11000 Times/min
Reciprocating stroke	4.2mm



JM81011900

Speed	0-37000 Times/min
Milling cutter	3mm

G series power



JM81012000

Speed	0-1200 Times/min
Range	0.8-8mm

G series power




JM81012100

Speed	0-800 Times/min
Range	0.8-8mm



JM81012200

Speed	0-290 Times/min
Range	0.8-8mm
Interface	stryker




JM81012300

Sawing frequency	16700 Times/min
Interface	AO


G series power

G series power



JM81012400

Sawing frequency	11000 Times/min
Reciprocating stroke	4.2mm



JM81012600

Speed	0-800 Times/min
-------	-----------------



JM81012500

Sawing frequency	11000 Times/min
Reciprocating stroke	4.2mm



JM81012700

Speed	0-37000 Times/min
Miling cutter	3mm

H series power

E-Series power



**JM81012800**


Hand External Bone Drill

Speed	0-1200 Times/min
Range	0.6-4mm

**JM81012900**

Hand External Hollow drill

Speed	0-1200 Times/min
Range	0.6-3.2mm



**JM81013000**

Hand External Swing saw

Sawing frequency	0-1800 Times/min
------------------	------------------

**JM81013100**


Multifunctional host

Speed/Frequency	0-18000 Times/min
-----------------	-------------------



Swing saw assembly

Frequency	0-18000 Times/min
No	JM81013101



Hollow drill assembly

Seed	0-1200 Times/min
Range	0.6-3.2mm
NO	JM81013102



Gram needle assembly

Seed	0-1200 Times/min
Range	0.8-2mm
NO	JM81013104



Bone drill assembly

Seed	0-1200 Times/min
Range	0.6-4mm
NO	JM81013103



Ao bone drill assembly

Seed	0-1200 Times/min
NO	JM81013105

Gypsum saw

Gypsum saw

TPLO circular saw

Speed	0-16000 Times/min
Voltage	12V

Accessories



Name

Battery  
JM81013201

Charger  
JM81013202

Disinfection channel  
JM81013203


Wrench  
JM81013204

Size	8mm	10mm	12mm	15mm	18mm	
	20mm	21mm	24mm	27mm	30mm	33mm



Saw blades

Equipment disinfection box



Gypsum saw

Voltage	Speed
110V	0-13000cpm
220V	0-13000cpm

Size

φ52

Size

φ43

Size

φ64



6 gears adjustable

Accessories series



**G series battery B**  
GB-DC13.2V 1800mAh  
GB-DC14.4V 1800mAh



**G series battery S**  
GB-DC13.2V 1800mAh  
GB-DC14.4V 1800mAh



**G series H Series  
E series C charger**  
G-CDQ 7.2V-14.4V



**E series H series battery**  
HS-DC7.2V 1800mAh



**G series H series N series Battery  
disinfection channel**  
GXDTD



**G Series N series D Series s series  
Chuck key**



**Drill chuck**  
0.8-8mm  
0.6-4mm

Accessories series



**D series s series charger**



**D series s series battery disinfection channel**  
SXDTD



**H Series E series collet key**



**D&S Series battery**  
S-DC13.2V 1800mAh  
S-DC14.4V 1800mAh



**Disinfection plug**











**Disinfection cover**

STERILIZATION CASE




Saw blade series

Saw blade



Bit



Transfer interface

The transfer interface can be customized according to customer needs

Main joint / Adapter	Stryker	Johnson	Big Ao	Hudson	Zhengtiar	Aikang	Chun Li
Stryker	●		●	●	●	●	●
Johnson	●		●	●	●	●	●
Big Ao	●	●		●	●	●	●
Hudson	●	●	●		●	●	●
Zhengtian	●	●	●	●		●	●
Aikang	●	●	●	●	●		●
Chun Li	●	●	●	●	●	●	