

Company profile



We are a China-based company dedicated to the research, development, manufacturing, and sales of specialized orthopedic surgical instruments and implants. Our product portfolio encompasses nine major categories: Trauma Plate Series, Locking Plate Series, Cranio-Maxillofacial Plate Series, intramedullary Nail Series, Spine Series, External Fixation System Series, Joint Series, Sports Medicine Solutions, and Power Tools.

We cooperate with world-class medical device companies and customize production to produce a variety of products that meet the laws and regulations and requirements of different countries.

"Insisting on quality first". Pioneering, pragmatic, enterprising and innovative, and strive to move forward.



CONTENTS

1

ANNULAR EXTERNAL FIXATOR SYSTEM

14

COMBINATIONAL EXTERNAL FIXATOR SYSTEM

39

SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

49

COMMON EXTERNAL FIXATOR SYSTEM

59

ACCESSORIES OF EXTERNAL FIXATOR SYSTEM

Transparent Ring Frame Combination (PEEK)

Product Number:JM61010001



Transparent Ring Frame Combination(PEEK)

Product Number:JM61010002

Tibial Fixator Transparant (PEEK)

Product Number:JM61010003



Calcaneal Fixator Transparant (PEEK)

Product Number:JM61010004



Ring Connection

Product Number:JM61010005



2/3 Rings

Product Number:JM61010007

Size: $\Phi 125 \setminus \Phi 155 \setminus \Phi 175$
 $\Phi 195 \setminus \Phi 215$



1/3 Rings

Product Number:JM61010006

Size: $\Phi 125 \setminus \Phi 155 \setminus \Phi 175 \setminus \Phi 195 \setminus \Phi 215$



Foot Rings

Product Number:JM61010008



Rod To Rod Coupling I

Product Number:JM61010009

Size: $\Phi 6/6$



Rod To Rod Coupling II

Product Number:JM61010010

Size: $\Phi 6/6$



Pin To Ring Coupling

Product Number:JM61010011

Size: $\Phi 1.5/3$



Two-hole K-wire Pin

Product Number:JM61010012



90° Extension Plate

Product Number:JM61010013



Single Screw Clamp

Product Number:JM61010014



Diagonal Tie

Product Number:JM61010015



Single Direction joint With Displacement

Product Number:JM61010016



Ring Bar

Product Number:JM61010017
Size: 50\100\150



K-wire Pin

Product Number:JM61010018



Ring External Fixator

Product Number:JM61010027
Used for: Tibia, Femur
Configuration:
Half rings 6 pcs
M6×300 threaded rod 3 pcs
M6×20 thread 6 pcs
M6×20 nut 6 pcs
Pin clamp 16 pcs



Quick Connecting Rod

Product Number:JM61010019



Single Direction Joint

Product Number:JM61010020



Micro-motion Part For Rings

Product Number:JM61010021



Ring Clamp

Product Number:JM61010022



T Wrench

Product Number:JM61010023



Kirschner Wire Tensioning Pliers

Product Number:JM61010024



Kirschner Wire Guider

Product Number:JM61010025



Open Spanner

Product Number:JM61010026



Lower Limb Lengthening External Fixator

Product Number:JM61010029
Used for: Tibia, Femur
Configuration:
Rings 4 pcs
Inner threaded connecting rod 4 pcs
M6×180 threaded rod 4 pcs
M6×120 threaded rod 4 pcs
Pin clamp 16 pcs



2/3 Ring External Fixator

Product Number:JM61010028
Used for: Tibia, Femur
Configuration:
Rings with 6 groove 3 pcs
M6×400 threaded rod 3 pcs
Pin clamp 12 pcs



Ring External Fixator

Product Number: JM61010030
 Used for: Tibia, Femur
 Configuration:
 Rings 4 pcs
 M6×400 threaded rod 4 pcs
 Adjusting nut 4 pcs
 Pin clamp 16 pcs
 Annular tubes 4 pcs



Knee Joint External Fixator

Product Number: JM61010034
 Used for: Knee joint
 Configuration:
 Rings 4 pcs Foot ring 1 pc
 Half rings 1 pc
 Annular diagonal tie 1 pc
 Post with 4 edges 2 pcs Single joint 2 pcs
 Pin clamp 16 pcs
 M6 threaded rod 12 pcs



Elbow Joint External Fixator

Product Number: JM61010031
 Used for: Elbow joint fixator
 Configuration:
 Rings 4 pcs
 Annular diagonal tie 1 pc
 Single joint 2 pcs
 M6×200 threaded rod 8 pcs
 Pin clamp 16 pcs
 2-hole connecting plate 1 pc



Ankle Joint External Fixator

Product Number: JM61010035
 Used for: Equinus
 Configuration: Rings 2 pcs
 Foot ring 1 pc L connecting 2 pcs
 Half rings 1 pc
 Annular diagonal tie
 Single joint 2 pcs
 M6×200 threaded rod 3 pcs
 Pin clamp 14 pcs



Ankle External Fixator

Product Number: JM61010032
 Used for: Equinus
 Configuration: Rings 4 pcs
 L connecting 2 pcs
 Half rings 1 pc Annular diagonal tie 1 pc
 Inner threaded post 4 pcs Single joint 2 pcs
 M6 threaded rod 8 pcs
 Pin clamp 24 pcs

Parts & Instruments Of Ring External Fixator

▶ Rings

Product Number:JM61010068 Φ100mm
 Product Number:JM61010069 Φ120mm
 Product Number:JM61010070 Φ140mm
 Product Number:JM61010071 Φ160mm
 Product Number:JM61010072 Φ180mm
 Product Number:JM61010073 Φ200mm
 Product Number:JM61010074 Φ220mm
 Material: Aluminum




▶ Half Rings

Product Number: JM61010075 Φ100mm
 Product Number:JM61010076 Φ120mm
 Product Number:JM61010077 Φ140mm
 Product Number:JM61010078 Φ160mm
 Product Number:JM61010079 Φ180mm
 Product Number:JM61010080 Φ200mm
 Product Number:JM61010081 Φ220mm
 Material: Aluminum



▶ 2/3 Rings

Product Number: JM61010082 Φ100mm
 Product Number:JM61010083 Φ120mm
 Product Number:JM61010084 Φ140mm
 Product Number:JM61010085 Φ160mm
 Product Number:JM61010086 Φ180mm
 Product Number:JM61010087 Φ200mm
 Product Number:JM61010088 Φ220mm
 Material: Aluminum




▶ Foot Ring

Product Number:JM61010089 L
 Product Number:JM61010090 M
 Product Number:JM61010091 S
 Material: Stainless steel



▶ Annular Diagonal Tie

Product Number:JM61010092



▶ Threaded Rod M6

Product Number:JM61010093
 40-400mm



▶ Adjusting Nut

Product Number:JM61010094



▶ Annular Tubes

Product Number:JM61010095



▶ Spring

Product Number:JM61010096



▶ Single Joint

Product Number:JM61010097
 Material: Stainless steel



▶ Single Joint

Product Number:JM61010098
 Material: Stainless steel



▶ Universal Joint

Product Number:JM61010099
 Material: Stainless steel



Parts & Instruments Of Ring External Fixator

▶ Threaded Column

Product Number:JM61010100 1 holes
 Product Number:JM61010101 2 holes
 Product Number:JM61010102 3 holes
 Product Number:JM61010103 4 holes
 Product Number:JM61010104 5holes
 Product Number:JM61010105 6holes
 Product Number:JM61010106 7 holes
 Product Number:JM61010107 8 holes
 Material: Stainless steel



▶ Nut Column

Product Number:JM61010108 1 holes
 Product Number:JM61010109 2 holes
 Product Number:JM61010110 3 holes
 Product Number:JM61010111 4 holes
 Product Number:JM61010112 5 holes
 Product Number:JM61010113 6 holes
 Product Number:JM61010114 7 holes
 Product Number:JM61010115 8 holes
 Material: Stainless steel



▶ L Connecting Plate

Product Number:JM61010116 5 holes
 Product Number:JM61010117 6 holes
 Product Number:JM61010118 7 holes
 Product Number:JM61010119 8 holes
 Product Number:JM61010120 9 holes
 Product Number:JM61010121 10 holes
 Product Number:JM61010122 11 holes
 Product Number:JM61010123 12 holes
 Material: Stainless steel



▶ Connecting Plate

Product Number:JM61010124 2 holes
 Product Number:JM61010125 3 holes
 Product Number:JM61010126 4 holes
 Product Number:JM61010127 5 holes
 Product Number:JM61010128 6 holes
 Product Number:JM61010129 7 holes
 Product Number:JM61010130 8 holes
 Material: Stainless steel



▶ Cross Connecting Plate

Product Number:JM61010131 2 holes
 Product Number:JM61010132 3 holes
 Product Number:JM61010133 4 holes
 Product Number:JM61010134 5 holes
 Product Number:JM61010135 6 holes
 Product Number:JM61010136 7 holes
 Product Number:JM61010137 8 holes
 Material: Stainless steel



▶ Inner Threaded Connecting Rod

Product Number:JM61010138 80mm
 Product Number:JM61010139 120mm
 Material: Stainless steel



▶ M6 Bolt

Product Number:JM61010140
 Size:M6x10\12\16\20\25\30
 Material: Stainless steel



▶ Pin Clamp

Product Number:JM61010141 Φ4
 Product Number:JM61010142 Φ5-6
 Material: Stainless steel



▶ Pin Clamp Thick

Product Number:JM61010143 Φ1.5mm
 Product Number:JM61010144 Φ2.5mm
 Material: Stainless steel



▶ Drill Bit Stop Unit

Product Number:JM61010145
 Material: Stainless steel



▶ Hex Nut M6

Product Number:JM61010146
 Material: Stainless steel



Ring External Fixator


REF: JM61010033
 Used for: Femur, Tibia
 Configuration: Rings with 1 ear, rings with six ears, standard post 6 pcs



Parts & Instruments Of Ring External Fixator

Ring With Six Ears

Product Number: JM61010036 L
 Product Number: JM61010037 M
 Product Number: JM61010038 S
 Material: Aluminum



Ring With 3 Ears

Product Number: JM61010039 L
 Product Number: JM61010040 M
 Product Number: JM61010041 S
 Material: Aluminum



Half Ring

Product Number: JM61010042 L
 Product Number: JM61010043 M
 Product Number: JM61010044 S
 Material: Aluminum



Foot Ring

Product Number: JM61010045 L
 Material: Aluminum



Post With 4 Edges

Product Number: JM61010046 1 holes
 Product Number: JM61010047 2 holes
 Product Number: JM61010048 3 holes
 Product Number: JM61010049 4 holes
 Product Number: JM61010050 5 holes
 Product Number: JM61010051 6 holes
 Material: Stainless steel



Post With 6 Edges

Product Number: JM61010052
 Material: Stainless steel



Rod To Ring Bolt

Product Number: JM61010053
 Material: Stainless steel



Bush

Product Number: JM61010054 Φ4
 Product Number: JM61010055 Φ5
 Product Number: JM61010056 Φ6
 Material: Aluminum



Pin to Ring Coupler

Product Number: JM61010057 Φ1.5
 Product Number: JM61010058 Φ2
 Product Number: JM61010059 Φ3
 Material: Stainless steel



Pin To Post Coupling

Product Number: JM61010060 Φ4mm
 Product Number: JM61010061 Φ5mm
 Material: Stainless steel



Standard Rod

Product Number: JM61010062 L
 Product Number: JM61010063 M
 Product Number: JM61010064 S



Packaged Rod

Product Number: JM61010065 L
 Product Number: JM61010066 M
 Product Number: JM61010067 S



Wrist Joint External Fixator

Product Number: JM61010147
 Used for: wrist
 Configuration:
 2pcs Φ160 full ring
 3pcs Φ160 half ring
 2pcs hinge
 1pc retractor
 7pcs M6×250 threaded rod



Ankle Joint External Fixator

Product Number: JM61010148
 Used for: Ankle joint
 Configuration:
 4pcs extender
 6pcs M6×150 threaded rod
 3pcs porous plate with 1 hole
 3pcs universal joint
 2pcs Φ160 ring
 2pcs Φ160U ring



Parts & Instruments

▶ **Ring**

Product Number:JM61010149
Size: Φ120\Φ140\Φ150\Φ160\Φ180
Φ200\Φ220\Φ240
Material: Aluminum




▶ **Half Ring**

Product Number:JM61010150
Size: Φ120\Φ140\Φ150\Φ160\Φ180
Φ200\Φ220\Φ240
Material: Aluminum



▶ **U Ring**

Product Number:JM61010151
Size: Φ120\Φ140\Φ150\Φ160\Φ180
Φ200
Material: Aluminum



▶ **90° Ring**

Product Number:JM61010152
Size: R90 Long\R90 Short\R120 Long\R120 Short
Material: Aluminum



▶ **Porous Plate**

Product Number:JM61010153
Size: 1\2\3\4\5\6\7\8 Holes
Material: Stainless steel



▶ **Six Corner Connection**

Product Number:JM61010154
Size: 15mm\30mm
Material: Stainless steel



▶ **Connecting Plate**

Product Number:JM61010155
Size: 20\30\40\50
Material: Stainless steel



▶ **U Type Ring Extension Rod**

Product Number:JM61010156
Size: 3\4\5\6\7\8 Holes
Material: Aluminum



▶ **Ball Head Screw**

Product Number:JM61010157
Size: Φ100\Φ150\Φ200\Φ280
Material: Stainless steel



▶ **90° Conversion**

Product Number:JM61010158
Material: Stainless steel



▶ **Square Nuts**

Product Number:JM61010159
Material: Stainless steel



▶ **Nylon Square Nut**

Product Number:JM61010160
Material: Nylon



▶ **Cone Washer**

Product Number:JM61010161
Material: Stainless steel



▶ **Cone Washer**

Product Number:JM61010162
Material: Stainless steel



▶ **Annular Tubes**

Product Number:JM61010163
Size: 10.20.30.40.50



Parts & Instruments

▶ **Connecting Rod**

Product Number:JM61010164
Material: Aluminum



▶ **Quick Coupling Rod**

Product Number:JM61010165
Material: Stainless steel



▶ **Universal Joint**

Product Number:JM61010166
Material: Stainless steel



▶ **Retractor**

Product Number:JM61010167
Size: 100\150\200
Material: Stainless steel



▶ **Retractor II**

Product Number:JM61010168
Material: Stainless steel



▶ **Threaded Rod**

Product Number:JM61010169



▶ **Needle Clamp**

Product Number:JM61010170
Material: Stainless steel



▶ **Square Needle Clamp**

Product Number:JM61010171
Material: Stainless steel



▶ **Pin To Ring Coupling**

Product Number:JM61010172
Material: Stainless steel



▶ **Hinge**

Product Number:JM61010173
Material: Stainless steel



▶ **L Type Hinge**

Product Number:JM61010174
Material: Stainless steel



▶ **Extender**

Product Number:JM61010175
Size: 70\100\150\200
Material: Aluminum



▶ **Tight Wire Clamp**

Product Number:JM61010176



▶ **Socket Wrench #10**

Product Number:JM61010177



▶ **Straight Screwdriver #3**

Product Number:JM61010178



▶ **Screw Driver**

Product Number:JM61010179
Size: 6mm



▶ **Double Ended Spanner**

Product Number:JM61010180
Size: #10




▶ **Socket Wrench**

Product Number:JM61010181
Size: #10



ANNULAR EXTERNAL FIXATOR SYSTEM

ANNULAR EXTERNAL FIXATOR SYSTEM




▶ **Lateral Limb Reconstruction External Fixator**

Product Number: JM61010182
 Size: $\Phi 4.0 \times 100$
 $\Phi 4.0 \times 160$

▶ **Large Skin Retractor**

Product Number: JM61010183
 Size: $\Phi 2.0 \text{mm} \times 250$
 $\Phi 2.5 \text{mm} \times 250$



▶ **Dual-purpose Skin Retractor**




Product Number: JM61010184

▶ **Small Skin Retractor**

Product Number: JM61010185
 Size: $\Phi 1.5 \text{mm} \times 250$
 $\Phi 2.0 \text{mm} \times 250$



Small Fragment External Fixator Instrument Set



Product NO.	NO.	Product Name	Specification	Qty
JM61020100	1	Rod To Rod Coupling	$\Phi 5/5$	12
JM61020200	2	Pin To Rod Coupling	$\Phi 5/3-4$	12
JM61020300	3	4-hole Pin Clamp	$\Phi 5/3-4$	3
JM61020400	4	Peri-articular Pin Clamp	$\Phi 3/4$	1
JM61020500	5	Post	$\Phi 5/0^\circ$	2
JM61020600	6	Post	$\Phi 5/30^\circ$	4
JM61020700	7	Apex Pin Self-drilling	HA 3.0×80	4
JM61020800			HA 3.0×100	4
JM61020900			HA 4.0×80	4
JM61021000			HA 4.0×100	4
JM61021100	8	Drill Sleeves	HA 4.0×120	4
JM61021200			$\Phi 3$	1
JM61021300			$\Phi 4$	1
JM61021400	9	Carbon Fiber Connecting Rod	$\Phi 5 \times 120$	4
JM61021500			$\Phi 5 \times 150$	4
JM61021600			$\Phi 5 \times 180$	2
JM61021700			$\Phi 5 \times 200$	4
JM61021800	10	T Wrench	#5	1
JM61021900	11	Stabilization/Reduction Wrench	#5	1
JM61022000	12	Elbow Joint Movement Device	$\Phi 5$	1
JM61022100	13	Apex Pin Chuck	$\Phi 3$	1
JM61022200			$\Phi 4$	1
JM61022300	14	T-shaped Pin Wrench	$\Phi 3$	1
JM61022400			$\Phi 4$	1
JM61022500	15	Hand Drill		1
JM61022600	16	Thumb Wheel	#5/7	1
JM61022700	17	Aluminum Box		1

Parts & Instruments Of Small Fragment External Fixator Instrument Set

▶ **Pin To Rod Coupling**

Product Number:JM61020200
Size: Φ5/3-4
Material:Aluminum



▶ **Rod To Rod Coupling**

Product Number:JM61020100
Size: Φ5/5
Material:Aluminum



▶ **Pin To Rod Coupling**

Product Number:JM61022800
Size: Φ5/3-4
Material:PEEK



▶ **Rod To Rod Coupling**

Product Number:JM61022900
Size: Φ5/5
Material:PEEK



▶ **Pin To Rod Integrated Coupling**

Product Number:JM61023000
Size: Φ3/4/5
Material:Aluminum



▶ **Peri-articular Pin Clamp**

Product Number:JM61020400
Size: Φ5/3-4
Material:Aluminum



▶ **4-Hole Pin Clamp**

Product Number:JM61020300
Size: Φ5/3-4
Material:Aluminum



▶ **Pin To Rod Coupling**

Product Number:JM61023100
Size: Φ3/2
Material:Aluminum



▶ **Rod To Rod Coupling**

Product Number:JM61023200
Size:Φ3/3
Material:Aluminum



Parts & Instruments Of Small Fragment External Fixator Instrument Set

▶ **Tube To Rod Coupling**

Product Number:JM61023300
Size: Φ15/5
Material:Aluminum



▶ **Connecting Rod**

Product Number:JM61023400/JM61021400
Material:Aluminum/Carbon fiber



▶ **Straight Post**

Product Number:JM61020500
Size: Φ5/0°
Material:Stainless steel



▶ **Elbow Joint Movement Device**

Product Number:JM61022000
Size: Φ5
Material:Stainless steel



▶ **30° Post**

Product Number:JM61020600
Size: Φ5/30°
Material:Stainless steel



▶ **Wrist Joint Movement Device**

Product Number:JM61023500
Size: Φ5
Material:Stainless steel



▶ **Compression/distraction Tube**

Product Number:JM61023600
Size: 15mm
Material:Aluminum



Parts & Instruments Of Small Fragment External Fixator Instrument Set

Stabilization/reduction Wrench

Product Number:JM61021900
Size: #15



Hand Drill

Product Number:JM61022500



Wrench

Product Number:JM61023700
Size: #3/#4/#5




T Wrench

Product Number:JM61021800
Size: #5



T-shaped Pin Wrench

Product Number:JM61022300/JM61022400
Size: Φ3/Φ4



Drill Sleeves

Product Number:JM61021200
Size: Φ3/Φ4



Apex Pin Chuck

Product Number:JM61022100/JM61022200
Size: Φ3/Φ4



Thumb Wheel

Product Number:JM61020200
Size:#5/7



Ulna & Radius External Fixator

Product Number:JM61030000
Used for: Ulna, Radius

REF	Product Name	Specification	Qty
JM61030001	Rod to rod coupling	Φ5/5	4
JM61030002	4-hole pin clamp	Φ5/3-4	2
JM61030003	Post	Φ5/30°	4
JM61030004	Connecting rod	Φ5×200mm	2
JM61030005	Apex Pin Self-drilling	Φ3×120mm	4



Forearm External Fixator

Product Number:JM61040000
Used for: Forearm

REF	Product Name	Specification	Qty
JM61040001	Rod to rod coupling	Φ5/5	2
JM61040002	4-hole pin clamp	Φ5/3-4	2
JM61040003	Post	0°	2
JM61040004	Connecting rod	Φ5×200mm	1
JM61040005	Apex Pin Self-drilling	Φ3×120mm	4



Distal Radius External Fixator

Product Number:JM61050000
Used for: Ulna, Radius

REF	Product Name	Specification	Qty
JM61050001	Rod to rod coupling	Φ5/5	1
JM61050002	Pin to rod coupling	Φ5/3-4	2
JM61050003	Peri-articular pin clamp	Φ5/3-4	1
JM61050004	Connecting rod	Φ5×200mm	1
JM61050005	Apex Pin Self-drilling	Φ3×120mm	4



► **Wrist Joint External Fixator**

Product Number: JM61060000
Used for: Wrist joint

REF	Product Name	Specification	Qty
JM61060001	Wrist joint movement device	Φ5	1
JM61060002	Pin to rod coupling	Φ5/3-4	4
JM61060003	Apex Pin Self-drilling	Φ3×120mm	4



► **Elbow Joint External Fixator**

Product Number: JM61070000
Used for: Elbow joint

REF	Product Name	Specification	Qty
JM61070001	Elbow joint movement device	Φ5	1
JM61070002	Pin to rod coupling	Φ5/3-4	4
JM61070003	Apex Pin Self-drilling	4×120mm	4



► **Finger & Palm External Fixator**

Product Number: JM61080000
Used for: Finger, Palm

REF	Product Name	Specification	Qty
JM61080001	Pin to rod coupling	Φ3/2	2
JM61080002	Connecting rod	Φ3×80mm	1
JM61080003	Apex Pin Self-drilling	Φ2×70mm	4



► **Humeral Shaft External Fixator**

Product Number: JM61090000
Used for: Humeral shaft

REF	Product Name	Specification	Qty
JM61090001	Rod to rod coupling	Φ5/5	2
JM61090002	4-hole pin clamp	Φ5/3-4	2
JM61090003	Post	Φ5/30°	2
JM61090004	Connecting rod	Φ5×200mm	1
JM61090005	Apex Pin Self-drilling	Φ4×120mm	4



► **Proximal Humerus External Fixator**

Product Number: JM61100000
Used for: Proximal humerus, Condyle neck

REF	Product Name	Specification	Qty
JM61100001	Rod to rod coupling	Φ5/5	4
JM61100002	Pin to rod coupling	Φ5/3-4	2
JM61100003	4-hole pin clamp	Φ5/3-4	1
JM61100004	Connecting rod	Φ5×250mm	2
JM61100005	Connecting rod	Φ5×200mm	1
JM61100006	Post	Φ5/30°	2
JM61100007	Apex Pin Self-drilling	Φ4×120mm	4






► **Humerus External Fixator**

Product Number: JM61110000
Used for: Bone shaft, distal & proximal bone

REF	Product Name	Specification	Qty
JM61110001	Compression/distraction tube	Φ15	1
JM61110002	Tube to rod coupling	Φ15/5	2
JM61110003	4-hole pin clamp	Φ5/3-4	2
JM61110004	Post	Φ5/30°	2
JM61110005	Apex Pin Self-drilling	Φ4×120mm	4

Large Fragment External Fixator Instrument Set

Parts & Instruments Of Large Fragment External Fixator Instrument Set



Product NO.	NO.	Product Name	Specification	Qty
JM61200100	1	Rod To Rod Coupling	Φ8/8	12
JM61200200	2	Pin To Rod coupling	Φ8/5-6	12
JM61200300	3	5-hole Pin Clamp	Φ8/5-6	3
JM61200400	4	Post	Φ8/0°	2
JM61200500	5	Post	Φ8/30°	4
JM61200600			HA 5.0×120	4
JM61200700	6	Apex Pin Self-drilling	HA 5.0×150	8
JM61200800			HA 5.0×180	4
JM61200900	7		Drill Sleeve	Φ5
JM61201000			Φ8×200	4
JM61201100	8	Carbon Fiber Connecting Rod	Φ8×250	4
JM61201200			Φ8×300	4
JM61201300			Φ8/180	2
JM61201400	9	Semi-circular Curved Rod	Φ8/ 200	2
JM61201500	10	Stabilization/Reduction Wrench	#7	1
JM61201600	11	Torque Wrench	Φ7	1
JM61201700	12	Thumb Wheel	#5/7	1
JM61201800	13	Apex Pin Chuck	Φ5	1
JM61201900	14	Hand Drill		1
JM61202000	15	T-shaped Pin Wrench	Φ5	1
JM61202100	16	T Wrench	#7	1
JM61202200	17	Aluminum Box		1

Pin To Rod Coupling

Product Number:JM61200100
Size:Φ8/5-6
Material:Aluminum



Rod To Rod Coupling

Product Number:JM61200200
Size:Φ8/8
Material:Aluminum



Pin To Rod Coupling

Product Number:JM61022800
Size:Φ8/5-6
Material:Aluminum



Rod To Rod Coupling

Product Number:JM61022900
Size:Φ8/8
Material:PEEK



Pin To Rod Integrated Coupling

Product Number:JM61023000
Size: Φ5/6/8
Material:Aluminum



5-hole Pin Clamp

Product Number:JM61200100
Size:Φ8/5-6
Material:Aluminum



Tube To Rod Coupling

Product Number:JM61203500
Size: Φ20/8
Material:Aluminum



Semi-circular Curved Rod

Product Number:JM61203300/JM61201300/JM61201400
Size: Φ8×160/180/200
Material:Aluminum



Straight Post

Product Number:JM61200500
Size: Φ8/0°
Material:Stainless steel



Parts & Instruments Of Large Fragment External Fixator Instrument Set

▶ **30°post**

Product Number:JM61200400
Size: Φ8/30°
Material:Stainless steel



▶ **Connecting Rod**

Product Number:JM61203400/JM61201000
Material:Aluminum/Carbon fiber



▶ **Compression/Distracton Tube**

Product Number:JM61203600
Size:Φ20
Material:Aluminum



▶ **Compression/Distracton Single Rod**

Product Number:JM61202800
Size: Φ8
Material:Aluminum



▶ **Compression/Distracton Double Rods**

Product Number:JM61202900
Size: Φ8
Material:Aluminum



▶ **Ankle Joint Movement Device**

Product Number:JM61203700
Size: Φ8
Material:Stainless steel



▶ **Knee Joint Movement Device**

Product Number:JM61203800
Size: Φ8 L/R
Material:Stainless steel



Parts & Instruments Of Large Fragment External Fixator Instrument Set

▶ **Torque Wrench**

Product Number:JM61201600
Size: #7



▶ **Stabilization/Reduction Wrench**

Product Number:JM61201500
Size: #22



▶ **Wrench**

Product Number:JM61203000
Size: #5/6/7



▶ **Hand Drill**

Product Number:JM61201900




▶ **T Wrench**

Product Number:JM61202200
Size: #7



▶ **T-shaped Pin Wrench**

Product Number:JM61202100
Size: Φ5/Φ6



▶ **Drill Sleeve**

Product Number:JM61200900/JM61203100
Size: Φ5/Φ6



▶ **Apex Pin Chuck**

Product Number:JM61201800
Size: Φ5/Φ6



▶ **Thumb Wheel**

Product Number:JM61201700
Size: Φ5/Φ7



► **Proximal Tibia External Fixator**

Product Number: JM61120000
Used for: Proximal tibia


REF	Product Name	Specification	Qty
JM61120001	Rod to rod coupling	Φ8/8	4
JM61120002	Pin to rod coupling	Φ8/5-6	3
JM61120003	5-hole pin clamp	Φ8/5-6	1
JM61120004	Post	Φ8/30°	2
JM61120005	Semi-circular Curved Rod	Φ8 180	1
JM61120006	Connecting rod	Φ8×250mm	2
JM61120007	Apex Pin Self-drilling	Φ5×150mm	5



► **Proximal Femur External Fixator**

Product Number: JM61150000
Used for: Petrochanteric

REF	Product Name	Specification	Qty
JM61150001	Rod to rod coupling	Φ8/8	4
JM61150002	S-hole pin clamp	Φ8/5-6	2
JM61150003	Post	Φ8/30°	4
JM61150004	Connecting rod	Φ8×300mm	2
JM61150005	Apex Pin Self-drilling	Φ5×150mm	2
JM61150006	Apex Pin Self-drilling	Φ5×180mm	2



► **Pelvic External Fixator**

Product Number: JM61160000
Used for: Pelvic

REF	Product Name	Specification	Qty
JM61160001	Rod to rod coupling	Φ8/8	6
JM61160002	5-hole pin clamp	Φ8/5-6	2
JM61160003	Post	Φ8/30°	2
JM61160004	Connecting rod	Φ8×300mm	2
JM61160005	Connecting rod	Φ8×200mm	2
JM61160006	Apex Pin Self-drilling	Φ5×150mm	4

► **Femoral Shaft External Fixator**

Product Number: JM61130000
Used for: Femoral shaft

REF	Product Name	Specification	Qty
JM61130001	Rod to rod coupling	Φ8/8	2
JM61130002	Pin to rod coupling	Φ8/5-6	4
JM61130003	Connecting rod	Φ8×200mm	1
JM61130004	Connecting rod	Φ8×250mm	2
JM61130005	Apex Pin Self-drilling	Φ5×150mm	4



► **Half Ring External Fixator**

Product Number: JM61140000
Used for: Tibia, perone

REF	Product Name	Specification	Qty
JM61140001	Rod to rod coupling	Φ8/8	6
JM61140002	5-hole pin clamp	Φ8/5-6	2
JM61140003	Post	Φ8/30°	4
JM61140004	Connecting rod	Φ8×300mm	2
JM61140005	Connecting rod	Φ8×250mm	1
JM61140006	Apex Pin Self-drilling	Φ5×150mm	4

► **Ankle Joint External Fixator**

Product Number: JM61170000
Used for: Ankle joint

REF	Product Name	Specification	Qty
JM61170001	Rod to rod coupling	Φ8/8	3
JM61170002	Pin to rod coupling	Φ8/5-6	4
JM61170003	Connecting rod	Φ8×250mm	2
JM61170004	Connecting rod	Φ8×200mm	1
JM61170005	Apex Pin Self-drilling	Φ5×150mm	4



► Knee Joint External Fixator

Product Number:JM61180000
Used for: Knee joint

REF	Product Name	Specification	Qty
JM61180001	Rod to rod coupling	Φ8/8	2
JM61180002	Pin to rod coupling	Φ8/5-6	4
JM61180003	Knee joint movement device	Φ8	1
JM61180004	5-hole pin clamp	Φ8/5-6	1
JM61180005	Post	Φ8/30°	1
JM61180006	Connecting rod	Φ8×150mm	1
JM61180007	Apex Pin Self-drilling	Φ5×150mm	6



► Ankle Fusion Joint External Fixator

Product Number:JM61190000
Used for: Ankle joint

REF	Product Name	Specification	Qty
JM61190001	Pin to rod coupling	Φ5/3-4	1
JM61190002	Pin to rod coupling	Φ8/5-6	4
JM61190003	Ankle joint movement device	Φ8	1
JM61190004	Apex Pin Self-drilling	Φ5×150mm	4
JM61190005	Apex Pin Self-drilling	Φ3×120mm	1



► Micro-motion Type External Fixator

Product Number:JM61210000

REF	Product Name	Specification	Qty
JM61210001	Tube to rod coupling	Φ20/8	2
JM61210002	Compression/distraction tube	Φ20	1
JM61210003	5-hole pin clamp	Φ8/5-6	2
JM61210004	Post	Φ8/30°	2
JM61210005	Apex Pin Self-drilling	Φ5×150mm	4




► Micro-motion Type External Fixator

Product Number:JM61220000

REF	Product Name	Specification	Qty
JM61220001	Tube-Rod Coupling	Φ20/8	2
JM61220002	Compression Distraction Tube	Φ20	1
JM61220003	S-Hole Pin Clamp	Φ8/5-6	2
JM61220004	Post	Φ8/30°	2
JM61220005	Apex Pin Self-drilling	Φ6×150mm	4



ZKK External Fixator Instrument Set-I





Product NO.	NO.	Product Name	Specification	Qty
JM61360100	1	Ringent Rod To Rod Coupling	Φ8/8	4
JM61360200	2	Ringent Pin To Rod Coupling	Φ8/4-5	8
JM61360300	3	Ringent Rod To Rod Coupling	Φ11/11	4
JM61360400	4	Ringent Pin To Rod Coupling	Φ11/5-6	12
JM61360500	5	T Clamp	Φ11/5	2
JM61360600	6	Shanz Pin	HA 4.0×90	4
JM61360700		Shanz Pin	HA 5.0×110	4
JM61360800		Shanz Pin	HA 5.0×130	6
JM61360900		Shanz Pin	HA 6.0×130	6
JM61361000	8	Shanz Pin	HB 5.0×150	2
JM61361100		Shanz Pin	HB 6.0×150	2
JM61361200	7	Drill Bits	Φ2.8×200	2
JM61361300		Drill Bits	Φ3.8×200	2
JM61361400		Drill Bits	Φ4.8×200	2
JM61361500	8	Screw Driver	Φ4	1
JM61361600		Screw Driver	Φ5	1
JM61361700	9	Tapered Trocar & Guide	Φ4	1
JM61361800		Tapered Trocar & Guide	Φ5	1
JM61361900		Tapered Trocar & Guide	Φ6	1
JM61362000	10	Open Wrench	#10	2
JM61362100	11	Pressurizer		1
JM61362200	12	Connecting Rod	Φ8×150	2
JM61362300		Connecting Rod	Φ8×200	2
JM61362400		Connecting Rod	Φ8×250	2
JM61362500		Connecting Rod	Φ8×300	2
JM61362600		Connecting Rod	Φ11×150	2
JM61362700		Connecting Rod	Φ11×200	2
JM61362800		Connecting Rod	Φ11×250	2
JM61362900	Connecting Rod	Φ11×300	2	
JM61363000	Connecting Rod	Φ11×350	2	
JM61363100	13	Aluminum Box		1

ZAS External Fixator Instrument Set-II




Product NO.	NO.	Product Name	Specification	Qty
JM61370100	1	Pin To Rod Coupling	Φ11/5-6	12
JM61370200	2	Ringent Rod To Rod Coupling	Φ11/11	9
JM61370300	3	T Clamp	Φ11/5	2
JM61370400	4	Shanz Pin	HA 5.0×110	4
JM61370500		Shanz Pin	HA 5.0×130	6
JM61370600		Shanz Pin	HB 5.0×150	2
JM61370700		Shanz Pin	HA 6.0×130	6
JM61370800		Shanz Pin	HA 6.0×150	2
JM61370900	5	Drill Bits	Φ3.8×200	2
JM61371000		Drill Bits	Φ4.8×200	2
JM61371100	6	Screw Driver	Φ5	1
JM61371200		Screw Driver	Φ6	1
JM61371300	7	Tapered Trocar & Guide	Φ4	1
JM61371400		Tapered Trocar & Guide	Φ5	1
JM61371500		Tapered Trocar & Guide	Φ6	1
JM61371600	8	Open Wrench	#10	2
JM61371700	9	Pressurizer		1
JM61371800	10	Connecting Rod	Φ11×150	4
JM61371900		Connecting Rod	Φ11×200	4
JM61372000		Connecting Rod	Φ11×250	4
JM61372100		Connecting Rod	Φ11×300	4
JM61372200		Connecting Rod	Φ11×350	4
JM61372300	11	Aluminum Box		1



ZJY External Fixator Instrument Set-III



Product NO.	NO.	Product Name	Specification	Qty
JM61380100	1	Pin To Rod Coupling	Φ8/3.5-4	8
JM61380200	2	Rod To Rod Coupling	Φ8/8	4
JM61380300	3	Rod To Rod Coupling	Φ11/11	4
JM61380400	4	Pin To Rod Coupling	Φ11/5-6	12
JM61380500	5	T Clamp	Φ11/5	2
JM61380600	6	Shanz Pin	HA 4.0×90	4
JM61380700		Shanz Pin	HA 5.0×110	4
JM61380800		Shanz Pin	HA 5.0×130	6
JM61380900		Shanz Pin	HA 6.0×130	6
JM61381000		Shanz Pin	HB 5.0×150	2
JM61381100		Shanz Pin	HB 6.0×150	2
JM61381200	7	Drill Bits	Φ2.8×200	2
JM61381300		Drill Bits	Φ3.8×200	2
JM61381400		Drill Bits	Φ4.8×200	2
JM61381500	8	Screw Driver	Φ4	1
JM61381600		Screw Driver	Φ5	1
JM61381700	9	Tapered Trocar & Guide	Φ4	1
JM61381800		Tapered Trocar & Guide	Φ5	1
JM61381900		Tapered Trocar & Guide	Φ6	1
JM61382000	10	Pressurizer		1
JM61382100	11	Open Wrench	#10	2
JM61382200	12	Connecting Rod	Φ8×150	2
JM61382300		Connecting Rod	Φ8×200	2
JM61382400		Connecting Rod	Φ8×250	2
JM61382500		Connecting Rod	Φ8×300	2
JM61382600		Connecting Rod	Φ11×150	2
JM61382700	13	Connecting Rod	Φ11×200	2
JM61382800		Connecting Rod	Φ11×250	2
JM61382900		Connecting Rod	Φ11×300	2
JM61383000	13	Connecting Rod	Φ11×350	2
JM61383100		Aluminum Box		1



Parts & Instruments Of ZKK,ZAS,ZJY External Fixator Instrument Set

Pin To Rod Coupling

Product Number:JM61390001
Size: Φ11/5
Material:Stainless steel



Pin To Rod Coupling

Product Number:JM61390002
Size: Φ11/5-6, Φ8/3.5-4
Material:Aluminum



T Clamp

Product Number:JM61390003
Size: Φ11/5
Material:Stainless steel



Rod To Rod Coupling

Product Number:JM61390004
Size: Φ11/11
Material:Stainless steel



Rod To Rod Coupling

Product Number:JM61390005
Size: Φ11/11, Φ8/8
Material:Aluminum



T Clamp

Product Number:JM61390006
Size: Φ8/4 Φ11/5
Material:Aluminum



Ringent Rod To Rod Coupling

Product Number:JM61390007
Size: Φ11/11, Φ11/8, Φ8/8
Material:Aluminum



Ringent Pin To Rod Coupling

Product Number:JM61390008
Size: Φ8/2.5-3.5, Φ8/4-5
Material:Aluminum



Ringent Pin To Rod Coupling


Product Number:JM61390009
Size: Φ11/5-6, Φ11/3-4
Material:Aluminum



Parts & Instruments Of ZKK,ZAS,ZJY External Fixator Instrument Set

Integrated Pin To Rod Coupling

Product Number:JM61390010
Size: Φ8/3-4, Φ11/5-6
Material:Aluminum



Integrated Rod To Rod Coupling

Product Number:JM61390011
Size: Φ8/8, Φ11/11
Material:Aluminum



Ringent Rod To Rod Coupling

Product Number:JM61390012
Size: Φ11/11
Material:Aluminum



Pin To Rod Coupling

Product Number:JM61390013
Size: Φ8/3-4, Φ11/5-6
Material:Aluminum



Semi-circular Curved Rod

Product Number:JM61390014
Size: Φ8mm
Material:Aluminum



Trapezoidal Curved Rod

Product Number:JM61390015
Size: Φ8mm
Material:Aluminum



Connecting Rod

Product Number:JM61390016
Size: 8×150/200/250/300
11×150/200/250/300/350



Parts & Instruments Of ZKK,ZAS,ZJY External Fixator Instrument Set

Curved Rod

Product Number:JM61390017
Size: Φ8mm
Material:Aluminum



Pressurizer

Product Number:JM61390018



Compression/distraction Rod

Product Number:JM61390019
Size: Φ8
Material:Stainless steel



Open Spanner

Product Number:JM61390020
Size: 10#



Tapered Trocar & Guide

Product Number:JM61390021
Size: Φ3mm, Φ4mm, Φ5mm, Φ6mm
Note: Single or Triple is OK




Hex Wrench

Product Number:JM61390022
Size: Φ2mm, Φ3mm, Φ4mm, Φ5mm, Φ6mm



Screw Driver

Product Number:JM61390023
Size: Φ2mm, Φ2.5mm, Φ3mm, Φ3.5mm, Φ4mm, Φ5mm, Φ6mm



Half Ring Type External Fixator

Product Number:JM61230000
Used for: Tibia & perone open & closed fracture

REF	Product Name	Specification	Qty
JM61230001	Ringent rod to rod coupling	Φ8/8	6
JM61230002	Ringent pin to rod coupling	Φ8/4-5	6
JM61230003	Semi-circular curved rod	Φ8	2
JM61230004	Connecting rod	Φ8×300	1
JM61230005	Compression/distraction rod	Φ8	2
JM61230006	Shanz Pin	HA 4.0×120	6



Condylar External Fixator

Product Number:JM61240000
Used for: Humeral epicondyle, tibial epicondyle, distal fracture

REF	Product Name	Specification	Qty
JM61240001	Ringent rod to rod coupling	Φ8/8	3
JM61240002	Ringent pin to rod coupling	Φ8/4-5	6
JM61240003	Semi-circular curved rod	Φ8	1
JM61240004	Trapezoidal curved rod	Φ8	1
JM61240005	Connecting rod	Φ8×250	1
JM61240006	Shanz Pin	HA 4.0×120	6



Tibial Joint Side External Fixator

Product Number:JM61250000
Used for: Distal tibial joint, proximal fracture

REF	Product Name	Specification	Qty
JM61250001	Ringent pin to rod coupling	Φ8/8	3
JM61250002	Ringent pin to rod coupling	Φ8/4-5	7
JM61250003	Semi-circular curved rod	Φ8	1
JM61250004	Curved rod	Φ8	1
JM61250005	Connecting rod	Φ8×300	1
JM61250006	Shanz Pin	HA 4.0×120	7

Z Frame(T Type) External Fixator

Product Number:JM61260000
Used for: Femur comminuted fracture

REF	Product Name	Specification	Qty
JM61260001	Ringent rod to rod coupling	Φ8/8	2
JM61260002	Ringent pin to rod coupling	Φ8/4-5	4
JM61260003	T clamp	Φ8/4	1
JM61260004	Connecting rod	Φ8×100	1
JM61260005	Connecting rod	Φ8×150	2
JM61260006	Shanz Pin	HA 4.0×120	6



T Frame External Fixator

Product Number:JM61270000
Used for: Femur, Tibia

REF	Product Name	Specification	Qty
JM61270001	Pin to rod coupling	Φ11/5-6	4
JM61270002	T clamp	Φ11/5	1
JM61270003	Connecting rod	Φ11×300	1
JM61270004	Shanz Pin	Φ5×130	4
JM61270005	Shanz Pin	HBΦ5×150	2



Straight Type External Fixator

Product Number:JM61280000
Used for: Femur, Tibia

REF	Product Name	Specification	Qty
JM61280001	Ringent rod to rod coupling	Φ11/11	1
JM61280002	Pin to rod coupling	Φ11/5-6	4
JM61280003	Connecting rod	Φ11×200	2
JM61280004	Shanz Pin	Φ5×130	4



Z Type External Fixator

Product Number:JM61290000
Used for: Femur, Tibia

REF	Product Name	Specification	Qty
JM61290001	Rod to rod coupling	Φ11/11	2
JM61290002	Pin to rod coupling	Φ11/5-6	4
JM61290003	Connecting rod	Φ11×200	2
JM61290004	Connecting rod	Φ11×150	1
JM61290005	Shanz Pin	Φ5×130	4



T Type External Fixator

Product Number:JM61300000
Used for: Tibia comminuted fracture

REF	Product Name	Specification	Qty
JM61300001	Pin to rod coupling	Φ11/5-6	4
JM61300002	T clamp	Φ11/5	1
JM61300003	Connecting rod	Φ11×300	1
JM61300004	Shanz Pin	Φ5×130	4
JM61300005	Shanz Pin	HBΦ5×150	2



Z Type(Enhancement Mode) External Fixator

Product Number:JM61310000
Used for: Femur, Tibia

REF	Product Name	Specification	Qty
JM61310001	Rod to rod coupling	Φ11/11	2
JM61310002	Pin to rod coupling	Φ11/5-6	4
JM61310003	Connecting rod	Φ11×200	2
JM61310004	Connecting rod	Φ11×150	1
JM61310005	Shanz Pin	Φ5×130	4



Straight Type External Fixator

Product Number:JM61320000
Used for: Tibia

REF	Product Name	Specification	Qty
JM61320001	Pin to rod coupling	Φ11×5-6	4
JM61320002	Connecting rod	Φ11×250	1
JM61320003	Shanz Pin	Φ5×130	4

Pelvic Fixator External Fixator

Product Number:JM61330000
Used for: Pelvic

REF	Product Name	Specification	Qty
JM61330001	Ringent rod to rod coupling	Φ11/11	2
JM61330002	Ringent pin to rod coupling	Φ11/5-6	4
JM61330003	Ringent rod to rod coupling	Φ11/8	2
JM61330004	Connecting rod	Φ11×200	2
JM61330005	Connecting rod	Φ11×300	1
JM61330006	Compression/distraction rod	Φ8	1
JM61330007	Shanz Pin	Φ5×150	4



Universal Pin/Rod To Rod Coupling



Z Type External Fixator

Product Number:JM61340000

REF	Product Name	Specification	Qty
JM61340001	Universal rod to rod coupling	Φ11/11	2
JM61340002	Universal pin to rod coupling	Φ6/16	4
JM61340003	Connecting rod	Φ11/200	2
JM61340004	Connecting rod	Φ11×150	1
JM61340005	Threaded rod	Φ6×130	4

Universal Rod To Rod Coupling

Product Number:JM61340006
Size: Φ11/11.Φ8/8



Universal Pin To Rod Coupling

Product Number:JM61340007
Size: Φ11/6.Φ8/4



Carbon Fiber Composite Ring



Three-dimensional Type External Fixator

Product Number:JM61350000

REF	Product Name	Specification	Qty
JM61350001	Rings	L	1
JM61350002	Connecting rod	Φ11/250	2
JM61350003	Rod to ring coupling	Φ11/11	2
JM61350004	Pin to ring coupling	Φ2.0	4
JM61350005	Ringlet pin to rod coupling	Φ11/5-6	2
JM61350006	Ringlet rod to rod coupling	Φ11/11	1
JM61350007	Threaded rod	Φ6×130	2
JM61350008	K-wire	Φ2.0×250	2

Tight Wire Clamp

Product Number:JM61350009



Rings

Product Number:JM61350010



Rod To Ring Coupling

Product Number:JM61350011
Size: Φ11/11



Pin To Ring Coupling

Product Number:JM61350012
Size: Φ2-4



◆ **Mini Fragment External Fixator-Type A**

Product Number:JM61400001
 Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 100mm



◆ **Mini Fragment External Fixator-Type B**

Product Number:JM61400002
 Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 90mm



◆ **Mini Fragment External Fixator-Type E**

Product Number:JM61400005
 Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 60mm



◆ **Mini Fragment External Fixator-Type F**

Product Number:JM61400006
 Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 70mm



◆ **Mini Fragment External Fixator-Type C**

Product Number:JM61400003
 Matching bone screws: $\Phi 2.5 \times 50\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 74mm



◆ **Mini Fragment External Fixator-Type D**

Product Number:JM61400004
 Matching bone screws: $\Phi 2 \times 50\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 45mm



◆ **Mini Fragment External Fixator-Type G**

Product Number:JM61400007
 Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 90mm



◆ **Mini Fragment External Fixator-Type H**

Product Number:JM61400008
 Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 70mm



▶ **Wrist External Fixator-Type A**

Product Number: JM6140009
 Matching bone screws: $\Phi 3.5 \times 80$ mm 4 pcs
 Matching instruments: 3.5mm screwdriver, 4mm hex wrench




▶ **Wrist External Fixator-Type E**

Product Number: JM6140012
 Matching bone screws: $\Phi 2.5 \times 60$ mm 2 pcs, $\Phi 3.5 \times 80$ mm 2 pcs
 Matching instruments: 3.5mm screwdriver, 4mm hex wrench



▶ **Wrist External Fixator-Type B**

Product Number: JM6140010
 Matching bone screws: $\Phi 2.5 \times 60$ mm 2 pcs, $\Phi 3.5 \times 80$ mm 2 pcs
 Matching instruments: 10 mm open spanner, 2.5&3.5mm screwdriver, 4mm hex wrench



▶ **Wrist External Fixator-Type F**

Product Number: JM6140013
 Matching bone screws: $\Phi 3.5 \times 80$ mm 4 pcs
 Matching instruments: 3.5mm screwdriver, 5mm hex wrench, 10mm open spanner



▶ **Wrist External Fixator-Type D**

Product Number: JM6140011
 Matching bone screws: $\Phi 2.5 \times 60$ mm 2 pcs, $\Phi 3.5 \times 80$ mm 2 pcs
 Matching instruments: 3.5mm screwdriver, 4mm hex wrench



▶ **Wrist External Fixator (Radiolucent Type)**

Product Number: JM6140014
 Matching bone screws: $\Phi 2.5 \times 60$ mm 2 pcs, $\Phi 3.5 \times 80$ mm 2 pcs
 Matching instruments: 3.5mm screwdriver, 4mm hex wrench



▶ **Wrist External Fixator (Radiolucent Type)**

Product Number: JM6140015
 Matching bone screws: $\Phi 2.5 \times 60$ mm 2 pcs, $\Phi 3.5 \times 80$ mm 2 pcs
 Matching instruments: 3.5mm screwdriver, 5mm hex wrench



Dynamic Axial External Fixator(Straight Type)

XL
Product Number:JM61400016
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
Product Number:JM61400017
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

M
Product Number:JM61400018
Matching bone screws: $\Phi 6 \times 130 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

S
Product Number:JM61400019
Matching bone screws: $\Phi 5 \times 120 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

XS
Product Number:JM61400020
Matching bone screws: $\Phi 4 \times 100 \text{mm}$ 4 pcs
Matching instruments: 4mm hex wrench, 4mm screwdriver



Dynamic Axial External Fixator(T Type)

XL
Product Number:JM61400021
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150 \text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
Product Number:JM61400022
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150 \text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

M
Product Number:JM61400023
Matching bone screws: $\Phi 6 \times 130 \text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150 \text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

S
Product Number:JM61400024
Matching bone screws: $\Phi 5 \times 120 \text{mm}$ 2 pcs, HB $\Phi 5.0 \times 130 \text{mm}$ 2 pcs
Matching instruments: 5mm hex wrench, 5mm screwdriver

XS
Product Number:JM61400025
Matching bone screws: $\Phi 4 \times 100 \text{mm}$ 4 pcs
Matching instruments: 4mm hex wrench, 4mm screwdriver



Dynamic Axial External Fixator(Ankle Joint Type)

XL
Product Number:JM61400026
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

L
Product Number:JM61400027
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

M
Product Number:JM61400028
Matching bone screws: $\Phi 6 \times 130 \text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver



Dynamic Axial External Fixator(Hip Joint Type)

XL
Product Number:JM61400029
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150 \text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
Product Number:JM61400030
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150 \text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

Petrochanteric External Fixator

Product Number:JM61400032
Matching bone screws: $\Phi 6 \times 150 \text{mm}$ 2 pcs
HB $\Phi 6.0 \times 180 \text{mm}$ 2 pcs
Matching instruments: 5mm hex wrench, 6mm screwdriver



◆ **Elbow Joint External Fixator**

Product Number: JM61400032
 Matching bone screws: $\Phi 5 \times 110\text{mm}$ 4 pcs
 Matching instruments: 5mm hex wrench, 3mm hex wrench, 6mm screw driver



◆ **Knee Joint External Fixator**

Product Number: JM61400033
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs, $\Phi 6 \times 130\text{mm}$ 2 pcs
 Matching instruments: 6mm hex wrench, 3mm hex wrench, 6mm screwdriver



◆ **Pelvic External Fixator-Type A**

Product Number: JM61400034
 Matching bone screws: HB $\Phi 6 \times 150\text{mm}$ 4 pcs
 Matching instruments: 6mm hex wrench, 6mm screwdriver
 Option: 5mm Pin



◆ **Pelvic External Fixator-Type B**

Product Number: JM61400035
 Matching bone screws: HB $\Phi 6 \times 150\text{mm}$ 4 pcs
 Matching instruments: 6mm hex wrench, 6mm screwdriver
 Option: 5mm Pin



◆ **Pelvic External Fixator-Type C**

Product Number: JM61400036
 Matching bone screws: HB $\Phi 6 \times 150\text{mm}$ 4 pcs
 Matching instruments: 6mm hex wrench, 6mm screwdriver
 Option: 5mm Pin



► **Bone Lengthening External Fixator-Type A**

Product Number:JM61400037 (250mm) For children
 Matching bone screws: $\Phi 4.0 \times 120\text{mm}$ 6 pcs
 Matching instruments: 4mm hex wrench, 4mm screwdriver

Product Number:JM61400038 (300mm)
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs
 Matching instruments: 6mm hex wrench, 6mm screwdriver

Product Number:JM61400039 (350mm)
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs
 Matching instruments: 6mm hex wrench, 6mm screwdriver

Product Number:JM61400040 (400mm)
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs
 Matching instruments: 6mm hex wrench, 6mm screwdriver



► **Guide Rail**

Product Number:JM61400044
 Size: 200, 250, 300, 400mm



► **Intermediate Clamp**

Product Number:JM61400045



► **Outer End Clamp**

Product Number:JM61400046



► **Adjustment Device**

Product Number:JM61400047
 Size: 50, 70, 110mm



► **T-shaped correction clamp**

Product Number:JM61400048
 Product Description : Suitable for correcting foot eversion Or varus



► **Bone Lengthening External Fixator-Type B**

Product Number:JM61400041 (250mm)
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs (Select 5mm bone screw)
 Matching instruments: 6mm hex wrench, 6mm screwdriver

Product Number:JM61400042 (300mm)
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs (Select 5mm bone screw)
 Matching instruments: 6mm hex wrench, 6mm screwdriver



► **T-clamp**

Product Number:JM61400049



► **Elastic Push Block**

Product Number: JM61400050
 Product Description : Apply pressure to the fracture line



► **Double Ring Bridge Block**

Product Number:JM61400051
 Product Description : For bridging with ring frame



► **Bone Lengthening External Fixator-Type C**

Product Number:JM61400043 (300mm)
 Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs,
 HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs
 Matching instruments:6mm hex wrench,
 6mm screwdriver

► **Fine Tuning Rotary Clamp**

Product Number:JM61400052
 Product Description : For correction of angular deformity



► **Multi-directional clamp**

Product Number:JM61400053
 Product Description : Suitable for correction of deformed rotation and movement



► **Universal T-clamp**

Product Number:JM61400054
 Product Description : For joint end



◆ **Straight Type External Fixator**

L
Product Number: JM61410001
Bone screws: $\Phi 6.0 \times 130\text{mm}$ 4PCS
Instruments: 6mm hex wrench, 6mm screw driver

M
Product Number: JM61410002
Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver (Select 5mm bone screw)

S
Product Number: JM61410003
Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 5mm screw driver

XS
Product Number: JM61410004
Matching bone screws: $\Phi 4.0 \times 90\text{mm}$ 4 pcs
Matching instruments: 4mm hex wrench, 4mm screw driver



◆ **T-Type External Fixator**

L
Product Number: JM61410005
Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

M
Product Number: JM61410006
Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 130\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

S
Product Number: JM61410007
Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 2 pcs, HB $\Phi 5.0 \times 120\text{mm}$ 2 pcs
Matching instruments: 5mm hex wrench, 5mm screw driver

XS
Product Number: JM61410008
Matching bone screws: $\Phi 4.0 \times 90\text{mm}$ 4 pcs
Matching instruments: 4mm hex wrench, 4mm screw driver



◆ **Ankle Joint External Fixator**

L
Product Number: JM61410009
Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

M
Product Number: JM61410010
Matching bone screws: $\Phi 5.0 \times 130\text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 5mm screw driver

◆ **Straight Type With Middle Pins External Fixator**

L
Product Number: JM61410011
Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

M
Product Number: JM61410012
Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver (Select 5mm bone screw)

S
Product Number: JM61410013
Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench, 5mm screw driver



◆ **T Type With Middle Pins External Fixator**

L
Product Number: JM61410014
Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

M
Product Number: JM61410015
Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 4 pcs, HB $\Phi 6.0 \times 130\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screw driver

S
Product Number: JM61410016
Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs, HB $\Phi 5.0 \times 120\text{mm}$ 2 pcs
Matching instruments: 5mm hex wrench, 5mm screw driver



◆ **3D Type External Fixator**

L
Product Number: JM61410017
Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 10-14mm open spanner

M
Product Number: JM61410018
Matching bone screws: $\Phi 5.0 \times 120\text{mm}$ 4 pcs
Matching instruments: 6mm hex wrench, 10-14mm open spanner

S
Product Number: JM61410019
Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs
Matching instruments: 5mm hex wrench, 10-14mm open spanner

XS
Product Number: JM61410020
Matching bone screws: $\Phi 4.0 \times 90\text{mm}$ 4 pcs
Matching instruments: 4mm hex wrench, 10mm open spanner



Three-section Straight Type External Fixator

Product Number: JM61410021
 Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 6 pcs
 Matching Instruments: 6mm hex wrench, 6mm screw driver



Three-section T-shaped Type External Fixator

Product Number: JM61410023
 Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs, $\text{HB}\Phi 6.0 \times 150\text{mm}$ 2 pcs
 Matching instruments: 6mm hex wrench, 6mm screw driver



Joint Type External Fixator

Product Number: JM61410022
 Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 2 pcs, $\text{HB}\Phi 6.0 \times 180\text{mm}$ 2 pcs
 Matching instruments: 6mm hex wrench, 6mm screw driver




Oblique Type External Fixator

Product Number: JM61410024
 Matching bone screws: $\Phi 6.0 \times 150\text{mm}$ 2 pcs, $\text{HB}\Phi 6.0 \times 180\text{mm}$ 2 pcs
 Matching instruments: 6mm hex wrench, 6mm screw driver



Dynamic External Fixator

Product Number: JM61410025
 Matching bone screws: $\Phi 4.0 \times 250\text{mm}$ 2 pcs, $\text{HB}\Phi 5.0 \times 250\text{mm}$ 3 pcs
 Matching instruments: 10mm open spanner, 4&5 mm screw driver



Stainless Steel Universal Type External Fixator- I

L
 Product Number: JM61410026
 Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs
 Matching instruments: 13-16mm open spanner, 4&5 mm hex wrench, 6mm screw driver



M
 Product Number: JM61410027
 Matching bone screws: $\Phi 5.0 \times 120\text{mm}$ 4 pcs
 Matching instruments: 10-14mm open spanner, 4&5 mm hex wrench, 5mm screw driver



S
 Product Number: JM61410028
 Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs
 Matching instruments: 10-14mm open spanner, 4&5 mm hex wrench, 5 mm screw driver

Stainless Steel Straight Type External Fixator-II

L
 Product Number: JM61410029
 Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs
 Matching instruments: 10-16mm open spanner, 4 mm hex wrench, 6 mm screw driver



M
 Product Number: JM61410030
 Matching bone screws: $\Phi 5.0 \times 120\text{mm}$ 4 pcs
 Matching instruments: 10-14mm open spanner, 4 mm hex wrench, 5 mm screw driver



S
 Product Number: JM61410031
 Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs
 Matching instruments: 10-14mm open spanner, 4 mm hex wrench, 5 mm screw driver

Mini-External Fixator(Finger & Palm)

Product Number: JM61410032
Used for: Finger, Palm



Mini-External Fixator(Finger & Palm)

Product Number:JM61410033
Used for: Finger, Palm



Wrist Joint External Fixator

Product Number:JM61410034
Used for: Wrist joint



Semi-circular Type(Standard) External Fixator

Product Number:JM61410035
Used for: Tibia & perone open/closed fracture

REF	Product Name	Specification	Qty
JM61410036	Rod to rod coupling	Φ8/8	6
JM61410037	Pin to rod coupling	Φ8/4-5	6
JM61410038	Semi-Circular Curved Rod	Φ8mm L	2
JM61410039	Rod	Φ8×300	1
JM61410040	Compression/distraction rod	Φ8	2
JM61410041	Shanz Pin	Φ4×120	2
JM61410042	K-wire	Φ4×250	2



Pelvic Fixator

Product Number:JM61410043
Used for: Pelvic


REF	Product Name	Specification	Qty
JM61410044	Rod to rod coupling	Φ8/8	6
JM61410045	Pin to rod coupling	Φ8/4-5	6
JM61410046	Semi-circular curved rod	Φ8mm	2
JM61410047	Compression/distraction rod	Φ8	1
JM61410048	Rod	Φ8×180	2
JM61410049	Shanz Pin	Φ5×130	6



Double Rods Tuberosity External Fixator

Product Number:JM61410050
Used for: Peritrochanteric

REF	Product Name	Specification	Qty
JM61410051	Double rods pin clamp	Φ8/4-5	5
JM61410052	Rod	Φ8×300	2
JM61410053	Shanz Pin	Φ5×130	2
JM61410054	Shanz Pin	Φ5×180	3



Double Rod Single Side External Fixator

Product Number:JM61410055
Used for: Femur, Tibia

REF	Product Name	Specification	Qty
JM61410056	Pin to double rod coupling	Φ8/4-5	4
JM61410057	Rod	Φ8×250	2
JM61410058	Shanz Pin	Φ5×130	4

Humeral Lateral Condylar Neck External Fixator

Product Number:JM61410059
Used for: Humeral lateral condylar neck, proximal & distal humer

REF	Product Name	Specification	Qty
JM61410060	Pin to rod coupling	Φ8/4-5	4
JM61410061	Rod to rod coupling	Φ8/8	1
JM61410062	Rod	Φ8×200	1
JM61410063	Rod	Φ8×150	1
JM61410064	Shanz Pin	Φ5×130	4



Humeral Frame External Fixator

Product Number:JM61410065
Used for: Humer

REF	Product Name	Specification	Qty
JM61410066	Rod to rod coupling	Φ8/8	2
JM61410067	Pin to rod coupling	Φ8/4-5	6
JM61410068	Rod	Φ8×250	2
JM61410069	Rod	Φ8×200	1
JM61410070	Shanz Pin	Φ4×120	6

◉ **Femoral Condylar External Fixator**

Product Number:JM61410071
Indication: Femoral supracondylar & intercondylar fracture

REF	Product Name	Specification	Qty
JM61410072	Rod to rod coupling	Φ8/8	3
JM61410073	Pin to double rod coupling	Φ8/4-5	2
JM61410074	Pin to rod coupling	Φ8/4-5	5
JM61410075	Semi-circular curved rod	Φ8mm L	1
JM61410076	Rod	Φ8×300	2
JM61410077	Curved rod	Long	1
JM61410078	Shanz Pin	Φ5×130	7



◉ **Tibial Plateau Frame External Fixator**

Product Number:JM61410095
Used for: Proximal tibia

REF	Product Name	Specification	Qty
JM61410096	Rod to rod coupling	Φ8/8	3
JM61410097	Pin to rod coupling	Φ8/4-5	8
JM61410098	Rod	Φ8×300	2
JM61410099	Semi-circular curved rod	Φ8mm L	1
JM61410100	Curved rod	Long	1
JM61410101	Shanz Pin	Φ5×130	6
JM61410102	K-wire	Φ4×250	1



◉ **Femoral Shaft External Fixator**

Product Number:JM61410079
Used for: Femoral shaft 1/3 to the joint comminuted fracture

REF	Product Name	Specification	Qty
JM61410080	Rod to rod coupling	Φ8/8	3
JM61410081	Pin to double rod coupling	Φ8/4-5	4
JM61410082	Pin to rod coupling	Φ8/4-5	3
JM61410083	Semi-circular curved rod	Φ8mm Long	2
JM61410084	Rod	Φ8×300	2
JM61410085	Shanz Pin	Φ5×330	4
JM61410086	K-wire	Φ4×250	1



◉ **Humeral Condylar External Fixator**

Product Number:JM61410103
Used for: Humeral & tibial lateral epicondyle

REF	Product Name	Specification	Qty
JM61410104	Rod to rod coupling	Φ8/8	3
JM61410105	Pin to rod coupling	Φ8/4-5	5
JM61410106	Rod	Φ8×250	1
JM61410107	Semi-circular curved rod	Φ8mm L	1
JM61410108	Curved rod	Φ8mm Short	1
JM61410109	Shanz Pin	Φ4×120	3
JM61410110	K-wire	Φ4×250	1



◉ **Ankle Joint External Fixator**

Product Number:JM61410087
Used for: Distal tibia, Ankle joint

REF	Product Name	Specification	Qty
JM61410088	Rod to rod coupling	Φ8/8	6
JM61410089	Pin to rod coupling	Φ8-4/5	8
JM61410090	Compression/distraction rod	Φ8	2
JM61410091	Rod	Φ8×350	2
JM61410092	Rod	Φ8×250	2
JM61410093	Shanz Pin	Φ4×120	4
JM61410094	K-wire	Φ4×250	2



◉ **Knee Joint External Fixator**

Product Number:JM61410111
Used for: Spanning joint, distal femur & proximal tibia

REF	Product Name	Specification	Qty
JM61410112	Rod to rod coupling	Φ8/8	8
JM61410113	Pin to rod coupling	Φ8/4-5	14
JM61410114	Semi-circular curved rod	Φ8mm L	2
JM61410115	Curved rod	Φ8mm Long	2
JM61410116	Spanning joint connecting rod	Φ8	2
JM61410117	Shanz Pin	Φ5×130	14

Parts And Tools For Combination External Fixator ZZL

Single Rod Clamp

Product Number: JM61410118
Size: $\Phi 6/2.5-3.5\text{mm}$, $\Phi 8/3-4\text{mm}$, $\Phi 8/4-5\text{mm}$, $\Phi 8/5-6\text{mm}$



Double Rods Clamp

Product Number: JM61410119
Size: $\Phi 6/6\text{mm}$, $\Phi 8/8\text{mm}$



Double Rods Pin clamp

Product Number: JM61410120
Size: $\Phi 8/3-4\text{mm}$, $\Phi 8/4.5\text{mm}$, $\Phi 8/5-6\text{mm}$



Semi-circular Curved Connecting Rod

Product Number: JM61410121
Size: $\Phi 6\text{mm L}$, $\Phi 6\text{mm S}$, $\Phi 8\text{mm L}$, $\Phi 8\text{mm S}$
Material: Stainless steel



Double Rod To Rod clamp

Product Number: JM61410122
Size: $\Phi 6/6\text{mm}$, $\Phi 8/8\text{mm}$



Trapezoidal Curved Connecting Rod

Product Number: JM61410123
Size: $\Phi 6\text{mm}$, $\Phi 8\text{mm}$
Material: Stainless steel



Curved Connecting Rod

Product Number: JM61410124
Size: $\Phi 6\text{mm L}$, $\Phi 6\text{mm S}$, $\Phi 8\text{mm L}$, $\Phi 8\text{mm S}$
Material: Stainless steel



Rod

Product Number: JM61410125



Parts And Tools For Combination External Fixator ZZL

Compression/distraction Rod

Product Number: JM61410126
Size: $\Phi 8$
Material: Stainless steel



Hexagonal



Knurling

Spanning Joint Connecting Rod

Product Number: JM61410127
Size: $\Phi 8$
Material: Stainless steel



Open Spanner

Product Number: JM61410128
Size: 10#



Open Spanner

Product Number: JM61410129
Size: 10-13#



Tapered Trocar & Guide

Product Number: JM61410130
Size: $\Phi 3\text{mm}$, $\Phi 4\text{mm}$, $\Phi 5\text{mm}$, $\Phi 6\text{mm}$
Note: Single or Triple is OK



Hex Wrench

Product Number: JM61410131
Size: $\Phi 2\text{mm}$, $\Phi 3\text{mm}$, $\Phi 4\text{mm}$, $\Phi 5\text{mm}$, $\Phi 6\text{mm}$



Screw Driver

Product Number: JM61410126
Size: $\Phi 2\text{mm}$, $\Phi 2.5\text{mm}$, $\Phi 3\text{mm}$, $\Phi 3.5\text{mm}$, $\Phi 4\text{mm}$, $\Phi 5\text{mm}$, $\Phi 6\text{mm}$



Shanz Pin

Product Number: JM74010001/JM74010084
 Diameter: 2.0/2.5/3.0/3.5/4.0/5.0/6.0mm
 Length: 50-250mm



Apex Pin Self-drilling

Product Number: JM74020001/JM7402125
 Diameter: 3.0/4.0/5.0/6.0mm
 Length: 60-250mm



Threaded Guide Wire

Product Number: JM74030001/JM7403125
 Diameter: 1.2/1.5/2.0/2.5/3.0mm
 Length: 100/120/150/250/300mm



Wire with Olive

Product Number: JM74040001/JM74040002/JM74040003
 Specification: 1.5x350, 2.0x350, 2.5x350mm



K-wire

Product Number: JM74050001/JM74050036
 Diameter: 0.8/1.0/1.2/1.5/1.6/1.8/2.0/2.5/3.0/3.5/4.0/4.5mm
 Length: 100/120/150/250/300mm



Connecting Rod

Product Number: JM74060001/JM74060036
 Length: 40-400mm
 Diameter: 3/5/6/8/11mm
 Material: Carbon fiber/Aluminum/Stainless steel

