


01 BIFIX® ADVANTAGE


- Safety**
Low sensitization in the hydrocolloid tape, good gas permeability and easy to stretch; Implantable polymer material, good biocompatibility, no skin irritation; No skin puncturing to reduce the chance of infection.
- Effectiveness**
Patented rotary adjustable closure, balanced tension, scar hyperplasia reduction
Well performed adhesion and cohesion;
The closing force is greater than 15 newtons and meets the wound closure requirement.
- Convenience**
Can be cut or superimposed application for various shapes of wounds;
Visible design for convenient wound observation and handling;
No need for suture or staple removal.

02 BIFIX® SPECIFICATION & MODEL



03 BIFIX® CLINICAL STUDIES


A clinical multicenter study of 144 cases showed that Bi-fix is superior in scar aesthetics, wound closure time, and postoperative removal pain to suture.



BIFIX® Surgical Incision Closure Set VS Stapler, traditional suture and absorbable suture


- Adjustable Wound closure, good blood per-fusion, good wound healing VS Large closing wound tension, insufficient blood per-fusion, easy to cause wound complications
- No needles, no foreign body stimulation, aesthetic wound healing VS secondary damage, foreign body stimulation, easy to cause tissue fibrosis and scar
- Simple and fast operation, operating time saving VS Long suturing time, complicated operating process
- No need for return and staple removal, it can be removed by the patients VS Need to remove the suture or staples, increase the burden of patients

04 BIFIX® CLINICAL APPLICATION \ OBSTETRICS




- Increase scar aesthetics
- Postoperative tension reduction
- Well-performed adhesion
- Low sensitization
- Precise adjustment
- Applicable for different shapes of wounds
- Cuttable or superimposed


05 BIFIX® CLINICAL APPLICATION \ ORTHOPEDICS



06 BIFIX® CLINICAL APPLICATION \ GENERAL SURGERY



07 BIFIX® CLINICAL APPLICATION \ OTHERS



Be Tolerant
To Diversity
Innovation
Technology

BIFIX®
Surgical Incision Closure Set



Zhejiang Haichuang Medical Device Co., Ltd.

Product name: Surgical Incision Closure Set
Trade mark: Bi-Fix®
Registration number: 20162650103
License number: 20160009
Packing size: Please see catalogue
Scope of application: Obstetrics and gynecology, general surgery, orthopedics, paediatric surgery, thoracic surgery, orthopedics and minimally invasive surgery
Storage: Store at room temperature, sealed, cool and dry place, avoid direct sunlight
Validity period: 30 months
Address: 1st Floor, No. 146, East Chaofeng Road, Yuhang Economic And Technological Development Zone, Yuhang District, Hangzhou, China
Contact: 01186-18616652496
Fax: 01186-571-87178683
Website: www.zjhycl.com

