



### **Shoulder Operation Back Plate of Operating Table**

- Surgical Fix Instrument ( Shoulder Operation Back Plate ), Beach Chair Position surgery, suitable for shoulder surgery and shoulder traction.
- The left and right backboards adopt a split design, which can be disassembled and positioned. Meet the requirements of medical staff for surgical posture.
- The Surgical Fix Instrument adopts plug-in design, which can be used and fixed with any brand of operating table.
- The width of the side rail of the operating table can be adjusted at will within the range of 500mm-600mm.

Daphne +86-18345168420 / daphne930120@gmail.com

- The angle of the Shoulder Operation Back Plate adopts the imported gas spring, and the back angle can be adjusted from  $-15^{\circ}$  /  $+80^{\circ}$ . It can be positioned at any position.
- The overall width of the back plate is 520mm and the height is 580mm. The cushion is made of high-density sponge with a thickness of 50mm.
- The head and neck support can be adjusted  $360^{\circ}$  universal and positioned arbitrarily.
- The height adjustment range of head and neck bracket is 0-320mm.
- It is equipped with an auxiliary trolley, which can be placed on the trolley and fixed to prevent falling off after the operation.