

Pain Therapy Patch DJTT-IC



Function and Feature

Pain Therapy Patch using electrode transfer exogenous weak electric current to the injured area, to increase the endogenous current, the injured tissue impedance will be reduced, the bioelectricity to re-enter the injured area, internal and external balance to restore injured cells, Pain Therapy Patch can stimulate and maintain many chemical and bioelectricity reaction of the cells in the healing process. The product does not contain chemical Medicine.

Portable, anytime and anywhere can be treated.

It works on the injured area continuously for 24 hours to promote the self-repair of the injured tissue and effectively relieve the pain.

Long - acting paste can be used repeatedly, the treatment time up to 500 hours.

Indications: low back pain, muscle cramps, rheumatism, shoulder pain, stiff neck, neck pain, knee and joint pain, tennis elbow, sciatic nerve pain, etc.



Battery	Built-in disposable 3V battery, wireless external power supply
Gel paste	It can be reused more than 30 times
Electrode cleaning	Rinse with water
Treatment current	21uA
Smart turn off	Remove the device and separate the electrode pieces
Treatment hours	Recommended to use 4 hours a day
Safety type	Class I BF type
Operating temperature	0°C~40°C
Working life	≥500 hours