



DENTAL CHAIR

unparalleled quality

美国 FDA 欧盟 CE MDR ISO13485 认证

AD-S6 睿智 26 款





互换操作功能

左右手的“互换操作功能”的实现，满足了特别医生对产品的使用要求，使得产品在使用上更加得心应手。

Left-right hand Operation

The dental chair supports left and right-hand interchangeable operation, fully considering the operating habits of different dentists, making the treatment process smoother.

进口亲肤型

采用“进口亲肤型”皮革+2CM记忆海绵，配备腰枕和头枕，皮质手感细腻乘坐舒适、同时兼备耐用、防酒精、抗老化、按摩加热功能

Comfortable and skin-friendly seat cushion

Imported leather + 2CM memory foam + double cushion design: delicate touch × comfortable for long periods of sitting × durable and anti-aging (alcohol-resistant/anti-aging)
Built-in massage and heating functions, providing warmth in the cold winter.



挑空式液晶操作台

主控台采用挑空式不锈钢托盘（可旋转 180°）设计及可高温消毒硅胶垫和不锈钢拉手控制；主控制板采用“医用级别”7寸液晶屏带 9 个记忆位操作界面、防触碰、防水、防尘，戴手套操作同样顺畅

Instrument Tray

180° rotatable stainless steel tray with high-temperature sterilizable silicone mat.
7-inch "medical-grade" LCD screen, touch-proof, waterproof, and dustproof, operable with medical gloves.
Stainless steel handle control;



高端医生座椅

靠背多角度调节

Doctor's Stool

The backrest can be adjusted at multiple angles



豪华护士椅

硅胶靠背

Nurse's Stool

Silicone backrest

※ 本公司保留不预先通知进行产品设计改良之权利，本资料图片仅供参考。
※ Our company reserves the right to make any improvement on product design without advance notice. the data pictures just for reference



广东雅登医疗器械有限公司

地址：佛山市南海区丹灶镇南沙社区工业
开发区逢拱路 12 号车间之二

销售电话：+86-0757-86603813

邮编：528216

传真：+86-0757-86603813

网址：www.dentalyd.com

邮箱：893666888@qq.com

Guangdong Yadeng Medical Apparatus Co., Ltd

Add: 2nd Factory, No.12 Fenggong Road, Nansha Industry
Zone, Foshan, 528216, China

Tel: +86-0757-86603813

Postcode: 528216

Fax: +86-0757-86603813

Web: www.dentalyd.com www.yadengdental.com

E-mail: A1@dentalyd.com A4@dentalyd.com



公众号



yadengdental