

Bring Technology to Life!



Viral DNA and RNA Extraction Kit

Tianlong's Viral DNA and RNA Extraction Kit is intended for rapidly extracting viral RNA/DNA from whole blood, serum, plasma, urine, swabs and other liquid samples. Realizes the high-throughput processing of parallel samples. The whole operation process is safe, convenient and no need for organic reagent. The extracted nucleic acids are highly purified and can be widely used in the fields of diagnostics, genomics research, disease detection, food safety and forensic identification, etc.



Features



Flexible options

more specifications for different scenarios



Various specimen

whole blood, serum,plasma, Urine, swabs and other liquid samples



Easy operation


Pre-filled design



High sensitivity

excellent purity and accurate detection

Application areas



Clinical diagnostics



Genomics research



Food safety



Forensic identification

Ordering information

Magnetic Beads method

Cat.No	Specification	Sample Type	Remark
T014H	64T/Kit (pre-filled), 16T/Plate x 4 Plates	Swab and other liquid samples	For COVID
T016H	20T/Box (Pre-filled) 1T/Strip		
T050H	64T/Kit (pre-filled), 16T/Plate x 4 Plates	Swab, Saliva and other liquid samples	For all viruses(5.0)
T051H	40T/Kit (pre-filled), 10T/Plate x 4 Plates		
T052H	20T/Kit (pre-filled), 5T/Plate x 4 Plates		
T325H	96T/Kit (pre-filled) 16T/Plate x 6 Plates	Serum, Plasma, Urine, Swab, etc.	*For all viruses(4.0) *For PANA9600S
T041H	48T/Kit (pre-filled), 8T/Plate x 6 Plates	Whole blood, Serum, Plasma, Tissue fluid, Urine, Swab medium, etc.	*For all viruses(4.0) *For GeneRotex48
T322H	40T/Kit (pre-filled), 10T/Plate x 4 Plates	Serum, Plasma, Urine, Swab, etc.	For all viruses(4.0)
T323H	20T/Kit (pre-filled), 5T/Plate x 4 Plates		
T324H	64T/Kit (pre-filled), 16T/Plate x 4 Plates		
T038H	20T/Box (Pre-filled) 1T/Strip	Serum; Plasma; Urine; Swab etc.	For all viruses
T104H	20T/Box (Pre-filled) 1T/Strip		
T111H	20T/Kit (pre-filled), 5T/Plate x 4 Plates	Serum, Plasma, Whole blood, etc.	For HCMV,EB
T112H	64T/Kit (pre-filled), 16T/Plate x 4 Plates		
T123H	20T/Box (Pre-filled) 1T/Strip		
T113H	20T/Kit (pre-filled), 5T/Plate x 4 Plates	Cytobrush collection medium, etc.	For HPV
T114H	64T/Kit (pre-filled), 16T/Plate x 4 Plates		
T115H	20T/Box (Pre-filled) 1T/Strip		

Spin column method

Cat.No	Specification	Sample Type	Remark
T108H	32T/Kit (Non-pre-filled)	Serum, Plasma, Urine,Swab etc.	For all viruses

version 2.0

All rights reserved by Tianlong. August 6, 2022



Tianlong Science and Technology

Mail: inquiry@medtl.com
Phone: 86 029 82682132
website: www.tlgenetech.cn
Address: No. 389 Zhuhong Road,Xi'an, China