

Rapid Test Kit(Colloidal Gold)



Introduction:

Hospitals, customs, health and epidemic prevention departments, pharmacies, personal and family use, etc.

Parameters:

No.	Name	Packing Parameters
1	Chikungunya IgG/IgM Combo Rapid Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
2	Combined Detection Kit For Influenza A/B (Colloidal Gold)	Packaging information: 20 pieces/box, 16 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 7 kg
3	Dengue Rapid Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
4	Helicobacter Pylori Antibody Rapid Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
5	Malaria P.Falciparum/P.Vivax Rapid Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
6	Typhoid Rapid Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
7	6-Monoacetylmorphine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
8	Amphetamine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
9	Benzodiazepine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
10	Buprenorphine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
11	MDMA Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
12	Cocaine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
13	Ketamine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
14	Morphine Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
15	Tetrahydrocannabinolic Acid Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
16	Luteinizing Hormone (LH) Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
17	Follicle-Stimulating Hormone (FSH) Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg
18	Beta-Human Chorionic Gonadotropin (β-HCG) Test Kit (Colloidal Gold)	Packaging information: 20 pieces/box, 30 boxes/carton Packaging size: 310*330*360 mm; Gross weight: 8 kg

Disposable Virus Sampling Tube Kit




Introduction:

In the current epidemic situation, virus sample collection is an important part of virus detection. The single-use virus sampling tube can collect, transport, Inactivate and store virus samples from specific parts of the human body.


Advantage:

- ①. Easy to operate and use.
- ②. Adding virus stabilizing ingredients can maintain the activity of the virus in a wide temperature range and reduce the decomposition rate of the virus (Active type).
- ③. Contains virus cleavage and virus nucleic acid preservation solution, which can quickly cleave the virus to release nucleic acid and store the nucleic acid stably (Inactive type).

Test Procedures:




1. According to different sampling requirements, use a sampling swab to sample at the corresponding part.




2. Place the swab after collecting the sample in a sterile sampling tube.



3. Put the swab into the preservation solution and break it from the breaking point.



4. Tighten the sterile sampling tube cover.



5. Mark the relevant sample information on the label of the sampling tube.

Product Components:

Item	Packing Information
10ml tube with 2ml medium; 50pcs/board	Total of 1200pcs, Packing size: 605*485*255mm, gross weight: 11.5kg
10ml tube with 3ml medium; 50pcs/board	Total of 1200pcs, Packing size: 605*485*255mm gross weight:12.5 kg
Remark	The swab can be optional