

Floor Type Low Speed Centrifuge



Features:

- * Whole metal structure and multi-layer explosion-proof, safe, durable, to ensure the safety of operators.
- * Embedded microprocessor control, DC brushless motor driving, maintenance-free, stable operation.
- * The 5-inch LCD touch screen, supporting switching between Chinese and English, minor languages can be customized.
- * All parameters are displayed real-time. The operation interface is intuitive and simple.
- * Electromechanical integration door lock, safe and reliable.
- * Integral food-grade silicone sealing ring is used to avoid aerosol overflow and ensure the safety of operator.
- * With self-test, blocking, overspeed, overtemperature and other alarm functions, excellent shock-absorbing performance, the motor runs smoothly.
- * Multiple groups of user-defined programs, easy to use common programs, boot for the last program.
- * Equipped with a variety of angle rotors and adapters, suitable for 0.2ml-750ml centrifugal tube, according to actual needs, to achieve multi-purpose.
- * Optional automatic rotor recognition.

Technical Parameters:

Model	BKC-DL6C
Max. Speed	6000rpm
Speed Accuracy	±10rpm
Max. RCF	5111×g
Max. Volume	4*750ml
Power Supply	AC220V, 50/60Hz
Timing Range	1s~99h59min, Continuous/Short-span
Noise	≤65dB
External Size(L*W*H)	500*635*820mm
Package Size(L*W*H)	700*580*1100mm
Net Weight	91kg
Gross Weight	105kg(220V), 110kg(110V)

Rotor:

NO.	Rotors Type	Volume	Max. Speed(rpm)	Max. RCF(×g)	Adapters		
5M00024		48*2/3/4/5/7ml vacuum blood collection tube	4000	2952	/		
5M00025		48*5/7ml vacuum blood collection tube		2952	2/3/4ml vacuum blood collection tube		
5M00026		72*2/3/4/5/7ml vacuum blood collection tube		2808	/		
5M00052		80*2/3/4/5/7ml vacuum blood collection tube		2808	/		
5M00027		96*2/3/4/5/7ml vacuum blood collection tube		3104	/		
5M00029		16*10/15ml		2957	2/3/4ml blood collection tube, 5/7ml blood collection tube, 15ml(▲)		
5M00030		24*10/15ml		2952			
5M00031		32*10/15ml		2952			
5M00080		48*15ml		3282	/		
5M00056		36*10ml		2862	/		
5M00057		64*10ml/15ml		3220	/		
5M00032		4*50ml		4500	3962	50ml(▲), 20ml, 15ml, 10ml, 5ml, 1.5/2ml	
5M00033		8*50ml		4000	2952		
5M00034		Swing Rotor		4*100ml	4500	3962	16*2/3/4ml blood collection tube, 16*5/7ml blood collection tube
5M00036				8*100ml	4000	3077	32*2/3/4ml blood collection tube, 32*5/7ml blood collection tube
5M00037	4*250ml(▲)		3578	/			
5M00038	4*250ml(flat bottom)		2780	48*1.5/2ml, 40*2/3/4/5/7ml blood collection tube, 48*2/3/4/5/7ml blood collection tube,			
5M00053	4*250ml(biosafety)		3041	24*10ml(▲), 28*10ml, 20*15ml, 4*50ml(▲), 8*50ml(●), 4*100ml			
5M00058	4*300ml		3220	/			
5M00039	4*500ml(flat bottom)		3130	48*2/3/4/5/7ml blood collection tube, 64*2/3/4/5/7ml blood collection tube, 36*10ml(●), 36*10ml(▲), 36*15ml(●), 36*15ml(▲), 12*50ml(●), 12*50ml(▲), 16*50ml(●), 16*50ml(▲), 48*15ml(▲), 4*250ml(flat bottom), 32*15ml(▲), 48*10ml(●)			
5M00054	4*500ml(biosafety)		3399				
5M00040	4*500ml(▲)		3578				
5M00041	4*750ml(flat bottom)		3399	108*2/3/4/5/7ml blood collection tube, 60*15ml(●), 60*15ml(▲), 24*50ml(▲), 24*50ml(●), 28*50ml(●), 28*50ml(▲)			
5M00055	4*750ml(biosafety)		3399	3399			
5M00047	Auto Uncap Swing Rotor		80*2/3/4/5/7ml blood collection tube(Auto uncap)	3220		/	
5M00048			108*2/3/4/5/7ml blood collection tube(Auto uncap)	3399		/	
5M00049			96*10ml vacuum blood collection tube(Auto Uncap)	3211		/	
5M00045	Swing		2*2*96 well	5000		2361	/
5M00043	Microplate Rotor		4*2*96 well			2987	/
5M00072		12*20ml	5000		3516	15ml, 10ml, 5ml, 1.5/2ml	
5M00073		12*15ml			3354	1.5/2ml, 2/3/4ml blood collection tube	
5M00074		8*15/20ml			2655		
5M00075		24*10ml			3354		
5M00076		Angle Rotor			18*10ml	3354	2/3/4ml blood collection tube, 5/7ml blood collection tube, 15ml(▲)
5M00077		12*10ml			3354		
5M00083		12*10/15ml			5111		
5M00069		24*10/15ml			6000	5031	50ml(▲), 20ml, 15ml, 10ml, 5ml, 1.5/2ml
5M00084		6*50ml			3904		

With mark(●) or no mark is round-bottom tube, with mark(▲) is sharp-bottom tube.