



MFZC1S



MFZC2S

Features

- The bed frame is formed by high-quality rigid welding, and the surface is electrostatically sprayed.
- The foot cover adopts wear-resistant and non-slip outer foot cover, which can better protect the foot of the bed and the ground.
- The bed surface is made of artificial leather, and the inner lining is made of 3cm thick sponge, which rebounds quickly.
- A storage rack is added at the bottom to place items, which is more convenient.(MFZC2S)

Technical Parameter

| Model | MFZC1S | MFZC2S |
|---------------------|-----------------------|--------|
| Product Size(W*D*H) | 1900*600*600mm(±10mm) | |
| Loading | ≤175kg | |
| Packing Size(W*D*H) | 1960*660*200mm | |
| Gross Weight | 30kg | 35kg |



MFZC3S



MFZC4S

Introduction

The examination bed is generally placed in the consultation room to check the specific condition of the patient. From small clinics to large hospitals, it is suitable for use in various departments of outpatient clinics, and can be used as lying or lying on patients during pre-diagnosis.

Features

- The bed frame is formed by high-quality rigid welding, and the surface is electrostatically sprayed.
- The foot cover adopts wear-resistant and non-slip outer foot cover, which can better protect the foot of the bed and the ground.
- The bed surface is made of artificial leather, and the inner lining is made of 3cm thick sponge, which rebounds quickly.
- A storage rack is added at the bottom to place items, which is more convenient.(MFZC4S)

Technical Parameter

| Model | MFZC3S | MFZC4S |
|---------------------|-----------------------|--------|
| Product Size(W*D*H) | 1800*600*600mm(±10mm) | |
| Loading | ≤175kg | |
| Packing Size(W*D*H) | 1860*660*200mm | |
| Gross Weight | 25kg | 30kg |