BIOBASE®



Features:

- * Excellent chemical resistance.
- * Fully autoclavable at 121℃.
- * Adjustable control knob makes setting the volume easy.
- * Made of PTFE, FEP, BSG,PP.
- * Vapor pressure Max. 500mbar, viscositymax. 500mm²/s, temperature max.40°C, density max.2.2g/cm².

Technical Parameters:

Cat.No.	Volume Range	Subdivision	Accuracy	Precision	Package	Gross Weight(kg)	Package Size(mm)
30013010	1ml	1	±2.00%, 20µl	1.00%, 10µl	20pcs/carton	14	595*430*275mm
30013020	0.5-5ml	0.1ml	±0.50%, 25µl	0.10%, 5µl	20pcs/carton	14	595*430*275mm
30013030	1-10ml	0.2ml	±0.50%, 50µl	0.10%, 10µl	20pcs/carton	14	595*430*275mm
30013040	2.5-25ml	0.5ml	±0.50%, 125µl	0.10%, 25µl	20pcs/carton	14	595*430*275mm
30013050	5-50ml	1.0ml	±0.50%, 250µl	0.10%, 50µl	20pcs/carton	14	595*430*275mm

Accesories:

Accessory Code	Description	Accessory Code	Description	
30013110	Adapt 45/40mm	30013150	Adapt 32/28mm	
30013120	Adapt 45/32mm	30013160	Bottle(Brown, 1L)	
30013130	Adapt 45/38mm	30013170	Filling tube(m)	
30013140	Adapt 32/25mm	30013180	Discharge tube(m)	

Automatic Micro Liquid Dispenser



Introduction:

The Automatic Micro Liquid Dispenser is equipped with a high-precision peristaltic pump and a robotic arm, which can quickly, accurately and consistently complete the automatic liquid dispensing operation, liberate the hands of the experimenter, and avoid mistakes while improving efficiency. The liquid dispensing accuracy of the product meets the minimum liquid dispensing volume as low as 5µl, and the filling time of 100µl in a 96-well plate is 10s.

Features:

- * Fast: 96-well plate 100µl adding time 10s
- * Accurate: 96-well plate 100µl liquid addition error ± 1.5%
- * Efficient: 10 times the efficiency of manual operation
- * Intuitive: Graphical operation interface, one-key liquid addition

Technical Parameters: 100 as a communication of the communication of the

Model	AMS-01		
Dispensing Volume	5-9999µl(1µl increment)		
Suitable Plate	8~96 well plate		
Dispensing Way	Whole plate or any line, the liquid volume of the specified line can be set		
Dispensing Speed	High-middle-low, can be adjustable		
Dispensing Pipe	8-channel		
Waste Liquid Handling	Movable water liquid slot		
Auxiliary Function	Shake, rinse,suck, pipe wash		
Display	7 inch touch display		
Dispensing Precision	20μl≤5%; 100μl≤1.5%(96-well plate)		
Dispensing Accuracy	20μl±1.5%; 100μl±1%(96-well plate)		
Power Supply	AC220/110V, 50/60Hz		
External Size	465*275*280mm		
Net Weight	10kg		
Packing Size	550*370*470mm		
Gross Weight	27kg		

225