

Optional Accessories:

									
Cat. No.	18900027	18900040	18900155	18900028	18900025	18900038	18900185	18900026	18900200
Model	SK330.1	SK330.3	SK330.5	SK330.2	SK180.1	SK180.3	SK180.4	SK180.2	SK180.6
Descriptions	Universal attachment with 4 bars	Lengthways roller attachment	Petri dish platform with non-slip mat and fixed string (8 pcs)	Fixing clip attachment	Universal attachment with 3 bars	Lengthways roller attachment	Petri dish platform with non-slip mat	Fixing clip attachment	Dish attachment with non-slip mat
	Size(W*D): 320*320mm	Size(W*D): 320*320mm	Size(W*D): 300*285mm		Size(W*D): 250*245mm	Size(W*D): 250*245mm	Size(W*D): 240*215mm		Size(W*D): 280x280mm
Used with Shakers									
SK-R330-Pro	√	√	√	√	/	/	/	/	/
SK-D3309-Pro	√	√	√	√	/	/	/	/	/
SK-O330-Pro/ SK-L330-Pro	√	√	√	√	/	/	/	/	/
SK-O180-Pro/ SK-L180-Pro	/	/	/	/	√	√	√	√	√

Email: export@biobase.com

								
Cat. No.	18900029	18900030	18900031	18900032	18900033	18900036	18900037	18900280
Descriptions	Fixing clip for flask volume 25ml, used with SK330.2 (56pcs)&SK 180.2 (36pcs)	Fixing clip for flask volume 50ml, used with SK330.2 (28pcs)&SK 180.2 (18pcs)	Fixing clip for flask volume 100ml, used with SK330.2 (28pcs) & SK180.2 (18pcs)	Fixing clip for flask volume 200ml/250ml, used with SK330.2 (13pcs) & SK180.2 (8pcs)	Fixing clip for flask volume 500ml, used with SK330.2 (9pcs) & SK180.2 (8pcs)	Bar for SK330.1 universal attachment	Bar for SK180.1 universal attachment	Connect sticks for double dishes shaker
Used with Shakers								
SK-R330-Pro	√	√	√	√	√	√	/	√
SK-D3309-Pro	√	√	√	√	√	√	/	√
SK-O330-Pro/ SK-L330-Pro	√	√	√	√	√	√	/	√
SK-O180-Pro/ SK-L180-Pro	√	√	√	√	√	/	√	√

Wrist Shaker



Features:

- * PID microprocessor control, LCD display.
- * Maintenance-free spherical bearings.
- * 8 customized programs, with different speed and time setting.
- * Automatic operation, auto-stop, timing, time display, parameters memory and recovery function.
- * Simulating manual oscillation from mild to strong shaking.
- * Suitable for flask, reagent bottles, tubes, jars, narrow mouth bottles, retort, separating funnel, straight mouth bottles.
- * Leakage current and over current safety protection system.
- * Can be used in incubators, warm rooms, environmental chambers and refrigerators.

Technical Parameters:

Model	SK-W810
Angle Shaking	12°
Rotary Speed	50~600rpm
Rotary Accuracy	±1rpm
Operation Mode	Continuous or timing
Max. Opening of the Clips	50mm
Max. loading	Each side is 6.8 Kg
Ambient Temp.	4°C~40°C
Ambient Humidity	20%~80%RH
Motor	Variable frequency motor
Material	Main body and arm 304 stainless steel
Standard Accessory	Main body, arm stand, flask clip*12 pcs
Consumption	60W
Power Supply	AC110/220V±10%, 50/60Hz
Main Body Size(W*D*H)	272*452*320mm
Complete Machine Size(W*D*H)	810*452*320mm
Package Size(W*D*H)	590*500*600mm
Net/Gross Weight	30kg/45kg