



MFCYA101

Introduction

The toilet chair is mainly suitable for the elderly or those who can not squat for a long time due to inconvenience of legs and feet. It can provide physical support for defecation and reduce the occurrence of dizziness and fainting accidents caused by standing after a long squat.

Technical Parameter

Model	MFCYA101
Seat height	380-510mm
Loading Capacity	100kg
External Size	570*465*(670-800)mm
Package Size	500*225*710mm
Net Weight	3.3kg
Gross Weight	4kg



Electric Shift Machine  
MFYW201(Upgrade model)



Manual Shift Machine  
MFYW102(Upgrade model)

Features

- It can be used with a toilet, and it can be used in many places, including bedrooms and toilets.
- Extended all-inclusive arc-shaped armrests, safer anti-collision and more convenient to hold.
- Thickened sponge cushion, more comfortable to use.
- Compared with similar products, it is more cost-effective.

Technical Parameter

Model	MFYW201(Upgrade model)	MFYW102 (Upgrade model)
External Size	630*530*(870-1100) mm	750*500*(750-900)mm
Brake	Foot Switch	
Driving Method	Electric drive pusher	Hydraulic pressure drive
Passing Width	560mm	500mm
Lifting Height	230mm	150mm
Load-bearing	120kg	
Power Supply	100-240V 50/60Hz	
Packing Size	750*540*410mm	760*550*420mm
Gross Weight	28kg	37kg