

## **Karl Fischer Titrator**

#### Features:

- \* LCD display screen.
- \* Support KF Titration Mode and Titer Detection Mode.
- \* Features of auto-filling, auto-purging and auto-mixing of the reagents ensure safe handling of Karl Fischer chemicals.
- \* Selectable units including mg, mg/L,%,ppm, etc.
- \* Settable parameters, including measurement unit, polarization current, stirring rate, titration rate stop volume, endpoint potential, stop criterion and etc.
- \* Store up to 200 titration data sets (GLP-compliant).
- \* Data can be easily transferred to printer via RS-232 commnication
- \* Reset feature automatically resumes all settings back to factory default



#### **Technical Parameters:**

Model		KFT-40C
Volumetric Titration	Range	0.1mg~250mg
	Resolution	0.1mg
	Polarization Current Accuracy	1μΑ±0.2μΑ; 50μΑ±10μΑ
	Repeatability	±0.5%
General Parameters	Power Supply	110~220V, 50/60Hz
	External Size(mm)	340*400*400
	Packing Size (mm)	540*490*480
	Net/Gross Weight (kg)	10/12

#### Features:

- \* High quality motor with excellent precision and accuracy
- \* Electronic control provides a fatigue-free operation
- \* Remote control panel prevents manual disturbance during operation
- \* Built-in magnetic stirrer for a complete and easy operation



#### **Technical Parameters:**

Model	dTrite	
Cat No.	7017100001	
Volume Range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10μL	
Volume Accuracy	R=0.2% CV=0.07%	
Velocity	16 stages	
Operating Temperature Range	10-30℃	
Quality Control	DIN EN ISO 8655	
Control Type	External control, including the stirrer and the burette program	
Power Supply	100~240VAC, 50/60 Hz	
Package Size(W*D*H)	300*240*220mm	
Gross Weight	1.62kg	

# Water Activity Meter

### Application:

Water activity meter, indicates the free water content when food is stable, measures the possibility of microbial proliferation. Food formulators use water activity to predict how much moisture migration will affect their product, so that to choose suitable package and storage method to reduce the usage of antiseptics.





#### Features:

- \* The measured object is unlimited, liquid solids and powder etc
- \* Printer output as standard accessory
- \* Chemical resistant. (only for BWA-6)

#### **Technical Parameters:**

Model		BWA-3A	BWA-6	
Working Temp.		0~50℃	-10°C~50°C	
Working Humidity		0~95%RH		
Display		LCD display		
Measuring Range	Temperature	0~50℃	-10℃~50℃	
	Water Activity	0~0.980	0~1.000	
Measurement Accuracy	Temperature	±0.5℃	±0.5℃	
	Water Activity	±0.015 (@25°C)	±0.012 (@23℃±5℃)	
Measurement Time		5~40min		
Measurement Points		1	1~4	
Temperature Coefficient		≤0.0005/°C		
Activity Repeatability		≤0.015	≤0.008	
Print Function		a. Print at the end of the measurement b. Real-time print		
Correction Mode		Single-point calibration	Multi-point joint calibration, single-point calibration	
Interface		1	RS232, USB	
Standard Accessories		Micro printer, one sensor	Data cable, micro printer, one sensor	
Consumption		15W	10W	
Power Supply		AC110/220V±10%, 50/60Hz		
Package Size(W*D*H)mm		550*430*240	740*430*240	
Gross Weight(kg)		11.5	15.5	

435 436