



MFCYA101

Introduction

The toilet chair is mainly suitable for the elderly or those who can not squat for a long time due to inconvenience of legs and feet. It can provide physical support for defecation and reduce the occurrence of dizziness and fainting accidents caused by standing after a long squat.

Technical Parameter

Model	MFCYA101
Seat Height	380~510mm
Loading Capacity	100kg
External Size(W*D*H)	570*465*(670~800)mm
Package Size(W*D*H)	510*230*720mm
Net Weight	3.3kg
Gross Weight	4kg



Electric Transfer Lift
MFYW201(Upgrade model)



Manual Transfer Lift
MFYW102(Upgrade model)

Features

- It can be used with a toilet, and it can be used in many places, including bedrooms and toilets.
- Extended all-inclusive arc-shaped armrests, safer anti-collision and more convenient to hold.
- Thickened sponge cushion, more comfortable to use.
- Compared with similar products, it is more cost-effective.

Technical Parameter

Model	MFYW201(Upgrade model)	MFYW102 (Upgrade model)
External Size(W*D*H)	630*530*(870~1100) mm	630*530*(750~900)mm
Brake	Foot switch	
Driving Method	Electric drive pusher	Hydraulic pressure drive
Passing Width	560mm	560mm
Lifting Height	230mm	150mm
Load-bearing	120kg	
Power Supply	AC100~240V 50/60Hz	
Packing Size(W*D*H)	760*550*420mm	760*550*420mm
Gross Weight	33kg	33.5kg