

Hotplate Magnetic Stirrer







Features:

- * High-quality cold-rolled steel plate outer shell, the surface is sprayed with environmental protection spray processing, good corrosion resistance.
- * Good performance, no noise, no vibration, fast heating
- * Artistic and reasonable, easy to maintain

 * LED display, real-time temperature display

Technical Parameters: Diobase.cc / www.biobase.com

Model	HMS-101S
Display	LED
Power	600W
Heating Temp. Range	R.T.~100℃
Stirring Power	2.5W
Stirring Capacity	50~2000ml
Work Plate Size	Ф200*100
Speed Range	0~2600rpm
Power Supply	220V, 50Hz
External Size	340*340*290mm
Package Size	390*390*360mm
Net Weight	6.4kg
Gross Weight	10kg

Note:

- 1. Please use a three-phase safety socket and properly grounded.
- 2. The instrument should be kept clean, it should be cut off the power supply when long-term not using.
- 3. Do not open the temperature control switch when the stainless steel container is not added with thermal oil and the temperature sensor is not connected, so as to avoid damage to the heating element and constant temperature meter.