

Amylose Tester



Introduction:

AT-II can be used for fast and online detection of the amylose quality of rice, corn, wheat and other cereals, as well as for testing and evaluating the quality of imported rice, improving the detection level and efficiency of grain departments and food production departments, and controlling grain quality and cost.

Features:

- * High sensitivity and good stability.
- * Intuitive interface, simple operation and convenient use.
- * Control software based on Windows operating system.
- * High degree of automation, fast detection speed and good repeatability.

Technical Parameters:

| Model | AT-II |
|-------------------|---------------------------------------|
| Testing Object | rice, corn, wheat, all kinds of crops |
| Testing Parameter | the content of the amylase |
| Testing Time | about 1 min |
| Wavelength | 650nm |
| Error | 0.1% |
| Temp Range | 20~30℃ |
| Power Supply | 220V 50/60Hz |
| Package Size | 470*600*350mm |
| Gross Weight | 14kg |