



MFSJY-01



MFSJY-02

The toilet is only for scene display.

### Introduction

Lift and tilt the seat plate through the electric push rod to assist the user to stand. Solve the elderly, pregnant women, the injured and other knee can not force the patient group, and solve the problem of down, up and so on.

### Technical Parameter

Model	MFSJY-01	MFSJY-02
Seat Height (Highest Point)	470-560mm	450-600mm
Handrail Spacing	450mm	
Internal Width	440mm	
Lifting Angle	0-14°	0-30°
Loading Capacity	150kg	
Power Supply	100-240V, 50/60Hz	
External Size	650*450*820mm	650*450*800mm
Package Size	650*550*820mm	
Net Weight	18kg	20kg
Gross Weight	21kg	23kg



MFMY101



MFMY102

### Introduction

Bath chair is mainly used for the elderly or people who cannot adapt to standing for a long time to reduce the occurrence of dizziness and fainting accidents caused by standing for a long time.

### Technical Parameter

Model	MFMY101	MFMY102
Adjustable Height	680-800mm	
Loading Capacity	150kg	
Material	Aluminium alloy, ABS	
External Size	510*490*690mm	470*500*730mm
Package Size	400*120*650mm	
Net Weight	2.6kg	3.1kg
Gross Weight	3.4kg	3.9kg