

## Laminar Flow Pathology Workstation



360° rotatable and adjustable spotlight, convenient for taking materials



Foldable hot and cold faucet



Formaldehyde perfusion universal faucet, The telescopic faucet is convenient for cleaning the countertop



Real-time display of intelligent control parameters, Human body induction switch is easy to operate

### Features:

- \* The open table design is adopted, and the front and rear of the operation area are unobstructed to ensure the maximum experimental space.
- \* Airflow pattern design with top air supply and bottom exhaust.
- \* Intelligent control system adopts LCD operation panel.
- \* The whole machine is made of 304 stainless steel, easy to clean.
- \* Standard formaldehyde filling device, built-in liquid tank, with liquid filling accessories and low liquid level alarm system, convenient for filling formaldehyde.

### Technical Parameters:

| Model                             | QCT-1800L  |
|-----------------------------------|--|
| External Size                     | 1800*800*2000mm  |
| Operating Table Size              | 1455*490mm   |
| Operating Table Adjustable Height | 800-950mm  |
| Airflow Pattern                   | Top air supply and bottom exhaust  |
| Standard Function                 | Table lifting, formaldehyde filling device   |
| Optional Function                 | Pulverizer, water heater, plasma system, mobile display stand, sampling plate exhaust air fan, new air fan |
| Power Supply                      | Standard: 220V 50/60Hz, optional: 110V 50/60Hz   |
| Net Weight                        | 350kg  |
| Packing Size                      | 1980*960*2250mm  |
| Gross Weight                      | 380kg  |