



SYIV100-MFLG813

Introduction

Wheelchair is an important mobile tool for home rehabilitation, turnover transportation, medical treatment, and outing activities for the wounded, sick, and disabled. It not only satisfies the mobility of the physically disabled and people with mobility impairments, but more importantly, it is convenient for family members to move and take care of patients. Allow patients to use wheelchairs to exercise and participate in social activities.

Features

- The stepless adjustment lever allows you to adjust the angle of back, full or half-lying as you wish.
- The pull-out type enlarges the toilet bowl. The toilet bowl is deepened and enlarged to prevent side leakage. The pull-out type design is clean and hygienic and does not release odor.
- The aluminum alloy wheel frame is equipped with air-free high-quality PU tires, which are easy to care for and maintenance-free.
- It folds compactly and can be placed in the trunk, making it easy to carry when traveling to meet various travel needs.

Technical Parameter

Model	SYIV100-MFLG813
Rear Wheel Configuration	PU wheel
Instrument Size (W*D*H)	1190*670*1200mm
Folded Size (W*D*H)	880*370*1000mm
Seat Size (W*D*H)	460*400*390mm
Seat Height	570mm
Seat Width	470mm
Seat Depth	470mm
Front Wheel Diameter	8 inches
Rear Wheel Diameter	24 inches
Loading Capacity	100kg
Net Weight	18.3kg
Gross Weight	21kg
Package Size(W*D*H)	750*330*720mm