

## Non-invasive Water Instrument




### Introduction:

It is a multi-functional injection penetration beauty device that employs TDA non-penetrating pulse air wave technology to deliver nutrient solutions to the deepest layer of the skin, thereby achieving the effects of hydrating, whitening, moisturizing, reducing wrinkles and antioxidation.

### TDA non-penetrating pulse air wave technology

TDA technology delivers hyaluronic acid, oxygen molecules, minerals, trace elements required by the skin, and nutrients which rich in moisture molecules (with the world's smallest molecular weight; one spoonful of the stock solution can cover an area equal to two football pitches) through the epidermis, penetrating layer by layer into every skin cell down to the periosteum layer. This enables the active ingredients to spread rapidly and evenly across the entire skin, delivering effective and long-lasting skincare results.



### Features:

- \* The operation time is short and the operation is simple.
- \* Highly secure, with no recovery period and without affecting daily life.
- \* The process is comfortable, no anesthesia is required, no needles are used, and no injections are given.
- \* The effect is highly and lastingly significant.

### Technical Parameters:

Model	BK-HDN01
Mode	Face/body
Handle Pressure	5kg/cm <sup>2</sup> ~10kg/cm <sup>2</sup>
Handle Depth	0.2mm~5mm
Power Supply	AC110V~240V, 50/60Hz
Consumption	500W
Standard Accessory	Needle*7
Optional Accessory	Trolley
External Size(W*D*H)	470*370*260mm
Net Weight	17kg
Packing Size(W*D*H)	600*435*480mm
Gross Weight	20kg