
Vertical Electrophoresis Tank



Application:


Vertical electrophoresis is used for separation, purification and preparation of nucleic acid and protein samples in life science research.

Features:

* High-transparency polycarbonate injection molding, durable and easy to observe.
* The integrated design of the glass cushion strip ensures that there is no leakage of glue.
* The special glue making rack makes the glue making operation simple and reliable.
* Open the cover and cut off the power to ensure the safety of the experiment.
* The upper cover limit function ensures that the polarity of the electrode will not be connected incorrectly.
* The power-off function of opening the lid ensures the safety of the experiment.
* It can make double-board glue at the same time, with clear strips.
* 99.99% high-purity platinum electrode, the best conductivity, strong corrosion resistance.
* Recommended power supply: BEP-300I, BEP-600I, BEP-600D.

Technical Parameters:

Model	BK-VET01
Glass Tray (W*L)	Thick type: 101*83mm, thin type: 101*73mm
Gel Trays(W*L)	Hand glue: 83*73mm
Comb Thickness	Standard 1.0mm, optional 0.75/1.5mm
Comb Throughput	Standard 10, 15 samples, optional 12 samples
Amount of Gel	1~2 pcs
NO. of Samples	10~30pcs
Buffer Volume	Internal tank: 130ml; external tank: 750ml
Standard Accessories	Body tank *1, Cap*1, Electrode core*1 set, Gel casting base*2, Gel frame*2, Single glue baffle*1, Gel shovel*1, Glass plates*5 pairs, 1.0mm thickness 10 wells comb*2, 1.0mm thickness 15 wells comb*2, glass drying rack*1
External Size (L*W*H)	160*115*150mm
Package Size(L*W*H)	390*250*270mm
Gross Weight	2.75kg



Application:

Vertical electrophoresis is used for separation, purification and preparation of nucleic acid and protein samples in life science research.

Features:

* High-transparency polycarbonate injection molding, durable and easy to observe.
* The integrated design of the glass cushion strip ensures that there is no leakage of glue.
* The special glue making rack makes the glue making operation simple and reliable.
* Open the cover and cut off the power to ensure the safety of the experiment.
* The upper cover limit function ensures that the polarity of the electrode will not be connected incorrectly.
* The power-off function of opening the lid ensures the safety of the experiment.
* It can make double-board glue at the same time, with clear strips.
* 99.99% high-purity platinum electrode, the best conductivity, strong corrosion resistance.
* Recommended power supply: BEP-300I, BEP-600I, BEP-600D

Technical Parameters:

Model	BK-VET02
Glass Tray (W*L)	Thick type: 101*83mm, thin type: 101*73mm
Gel Trays(W*L)	Hand glue: 83*73mm, precast glue: 86*68mm
Comb Thickness	Standard 1.0mm, optional 0.75/1.5mm
Comb Throughput	Standard 10 samples, optional 12/15samples
Amount of Gel	1~4 pcs
NO. of Samples	10~40pcs
Buffer Volume	700ml(1~2 pcs), 1000mll(3~4 pcs),
Standard Accessories	Body tank *1, Cap*1, Electrode core*1 set, Gel casting base*2, Gel frame*4,Single glue baffle*1, Gel shovel*1, Glass plates*5 pairs, 1.0mm thickness 10 wells comb*5, glass drying rack*1
External Size (L*W*H)	160*120*180mm
Package Size(L*W*H)	400*360*300mm
Gross Weight	3.75kg